# Homocysteine HPLC Kit

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## 17. LITERATUR


1. Role of Plasma Homocyst(e)ine in Arterial Occlusive Diseases. (1994) *Clinical Chemistry*, Vol. **40** No. 6, pp. 857-858.

2. Association between Plasma homocysteine concentrations and extracranial carotid-artery stenosis. Selhub, J. et al. (1995). *The New England Journal of Medicine*, Feb. 95, Vol. **332** No. 5, pp. 286-291.
3. Plasma Homocyst(e)ine as a Risk Factor for Early Familial Coronary Artery disease. Wu, L. et al. (1994). *Clinical Chemistry*, Vol. **40** No. 4, pp. 552-561.
4. Homocyst(e)ine and arterial occlusive diseases. Malinow, M. (1994). *Journal of Intern Medicine*, Vol. **236**, pp. 603-617.
5. Hormone replacement therapy may reduce high serum homocysteine in postmenopausal women. van der Mooren, M. (1994). *European Journal of Clinical Investigation*, Vol. **24**, pp. 733-736.
6. Homocysteine metabolism in pregnancies complicated by neural-tube defects. Mills, J. et al. (1995). *The Lancet*, Jan. 1995, Vol. **345**, pp. 149-151.

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

# Homocysteine HPLC Kit

*For the determination of homocysteine in plasma and serum*

Valid from 2015-12-23

**REF** **KC2801**

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>15</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>4. MATERIAL SUPPLIED</b>	<b>16</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>6. STORAGE AND PREPARATION OF REAGENTS</b>	<b>17</b>
<b>7. PRECAUTIONS</b>	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>18</b>
<b>9. ASSAY PROCEDURE</b>	<b>18</b>
<i>Procedural notes</i>	18
<i>Test procedure</i>	18
<i>Chromatographic conditions</i>	19
<b>10. TREATMENT OF THE COLUMN</b>	<b>19</b>
<b>11. RESULTS</b>	<b>19</b>
<i>Calculation</i>	19
<i>Typical chromatogram</i>	20
<b>12. LIMITATIONS</b>	<b>20</b>
<b>13. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<i>Controls</i>	21
<b>14. PERFORMANCE CHARACTERISTICS</b>	<b>21</b>
<i>Precision and reproducibility</i>	21
<i>Linearity</i>	21
<i>Detection limit</i>	21
<i>Recovery</i>	21
<b>15. DISPOSAL</b>	<b>21</b>
<b>16. TROUBLESHOOTING</b>	<b>22</b>
<b>17. REFERENCES</b>	<b>23</b>
<b>18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>

## 1. INTENDED USE

This HPLC application is intended for the quantitative determination of homocysteine in plasma and serum. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Homocysteine exists in plasma as protein-bound and, so-called, free forms. About 70 percent of the total plasma homocysteine is bound to albumin. The free forms include homocystine (homocysteine disulphide) and mixed disulphide. Just minute amounts of reduced homocysteine are found in plasma. Total homocysteine in plasma, including protein-bound and free forms, has been referred to as homocyst(e)ine [H(e)]. Studies for peripheral vascular, cerebrovascular and coronary artery disease (CAD), evaluating H(e) as a cardiovascular risk factor, have been consistently demonstrated that there are more persons with high H(e) among patients than among healthy persons. Increased concentrations of free forms of homocysteine have also generally been reported among patients with occlusive artery disease, with one notable exception of a negative finding in relation to CAD.

## 3. PRINCIPLE OF THE TEST

The first step in determining homocysteine includes sample preparation with additional derivatisation. Therefore, the internal standard is added to the calibrator, controls and samples. Reduction and derivatisation is carried out in one step (10 min at 60°C). During the reduction, homocysteine is cleaved from albumin and homocystine is reduced into two homocysteines. The derivatisation reagent transforms homocysteine into a fluorescent product. Higher molecular substances are removed by precipitation and centrifugation. 20 ml of the supernatant are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a reversed phase column. One run lasts 7 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak areas/heights by the internal standard method.

### Summary

The application of homocysteine for HPLC makes a fast and precise determination of the amino acid possible in an easy way. The kit includes all reagents for preparation and separation of the samples, except the column.

As for many other parameters, the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even

laboratories without experience in high performance liquid chromatography to use this technique for clinical chemical routines quickly and precisely. Mostly, a one-point calibration is sufficient for calibrating the test system – unlike immunoassays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. With short test series, the one-point calibration is much more economic than 6-point calibration for immunoassays.

#### 4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KC2801LM	MOPHA	Mobile phase	1000 ml
KC2801KA	CAL	Calibrator (lyoph. 1 ml; see label for concentration)	1 vial
KC2801IS	INT STD	Internal standard (lyoph. 12 ml)	1 vial
KC2801RE	RECSOL	Reconstitution solution	2x 15 ml
KC2801RL	REDSOL	Reduction solution (lyoph. 2.4 ml)	1 vial
KC2801DL	DER	Derivatisation solution	12 ml
KC2801FR	PREC	Precipitation reagent	12 ml
KC2801KO	CTRL 1 CTRL 2	Control 1 and 2 (lyoph. 250 µl.; see specification data sheet for concentration)	2x 3 vials

The HPLC column (KC2801RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- 1.5 ml reaction tubes (e.g. Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with fluorescence detector
- Reversed phase C<sub>18</sub> column
- Water bath / thermo shaker
- Vortex

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 6. STORAGE AND PREPARATION OF REAGENTS

- **The lyophilised calibrator (CAL)** is stable at **-20 °C** until the expiry date stated on the label. Before use, the CAL has to be reconstituted with **1 ml reconstitution solution (RECSOL)**. The concentration of homocysteine slightly changes from lot to lot, the exact concentration is stated on the label. **Calibrator (reconstituted CAL) has to be aliquotted and stored at -20 °C until use.**
- **The lyophilised internal standard (INT STD)** is stable at **-20 °C** until the expiry date stated on the label. Before use, the INT STD has to be reconstituted with **12 ml reconstitution solution (RECSOL)**. **Internal standard (reconstituted INT STD) has to be aliquotted and stored at -20 °C until use.**
- **The lyophilised controls 1 and 2 (CTRL 1 and CTRL 2)** are stable at **-20 °C** until the expiry date stated on the label. Before use, they have to be reconstituted with each **250 µl reconstitution solution (RECSOL)**. The concentration of homocysteine slightly changes from lot to lot, the exact concentration is stated on the specification data sheet.
- **The lyophilised reduction solution (REDSOL)** is stable at **-20 °C** until the expiry date stated on the label. Before use, the REDSOL has to be reconstituted with **2.4 ml reconstitution solution (RECSOL)**. **Reduction solution (reconstituted REDSOL) can be stored for 3 months at 2–8 °C.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## 7. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- The precipitation reagent consists of an acid. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.
- Reagents should not be used beyond the expiration date stated on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

Since the concentration in serum samples is depending on clotting time, EDTA-plasma should be preferred as sample. Please store the sample directly after blood collection at 2–8 °C. After blood collection, the sample should be centrifuged within the next 30 minutes (conditions: 2000 g, 10 min, 2–8 °C). Afterwards the sample should be stored at 2–8 °C.

For long time storage (> 1 week), keep the samples at -20 °C to -80 °C.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Pipet each <b>50 µl sample, calibrator or control 1 or 2</b> into an 1.5 ml reaction tube.
Add each <b>50 µl internal standard</b> .
Add each <b>20 µl reduction solution</b> .
Add each <b>100 µl derivatisation solution (DER)</b> and mix.
Incubate for <b>10 min at 60 °C</b> , then let the tubes cool down.
Add each <b>100 µl precipitation reagent (PREC)</b> and mix well.
<b>Incubate for 5 min at 2–8 °C</b> .
Centrifuge for <b>5 min at 10 000 g</b> .
The samples prepared like this are stable for at least <b>6 days at 2–8 °C</b> .
Inject <b>20 µl</b> supernatant into the HPLC system.

### Chromatographic conditions

Column material:	MZ Inertsil ODS-2; 5 µm
	MZ PerfectBond ODS-2; 5 µm
	Bischoff Prontosil Eurobond; 5 µm
Column dimension:	125 mm × 4 mm
Flow rate:	0.7–1.0 ml/min
	Please refer to the quality certificate of the column
Fluorescence detection:	Excitation: 385 nm
	Emission: 515 nm
Temperature:	30 °C
Injection volume:	20 µl
Running time:	7 min

It is recommended to use a guard column to extend column life.

## 10. TREATMENT OF THE COLUMN

After analysis, the column should be flushed with 30 ml ultra pure water (1 ml/min) and stored in 50% methanol in water (~ 30 ml, flow 0.6 ml/min). Before use, the system should be equilibrated with ~ 30 ml mobile phase (MOPHA).

## 11. RESULTS

### Calculation


$$F = \frac{\text{Peak height sample} \times \text{calibrator concentration}^*}{\text{Peak height internal standard of the sample}} \times F = \text{sample concentration}$$

$$F = \frac{\text{Peak height internal standard of the calibrator}}{\text{Peak height calibrator}}$$

\* see label

**Tip:** Alternatively, the peak area instead of the peak height can be used for quantification.

## Typical chromatogram

## 12. LIMITATIONS

Since the concentration in serum samples is depending on clotting time, EDTA-plasma should be preferred as sample.

Do not use whole blood.

## 13. QUALITY CONTROL

### Reference range

Homocysteine values: (Resch, "Homocystein", ISBN 3-920328-17-5)

< 15 $\mu\text{mol/l}$	Normal value
15-30 $\mu\text{mol/l}$	Lack of vitamins
31-100 $\mu\text{mol/l}$	Indication of heterocygote homocysteinemia
> 100 $\mu\text{mol/l}$	Indication of homocygote homocysteinemia

After a methionine loading test an increase in homocysteine concentration is normal. Concentrations more than 30  $\mu\text{mol/l}$  are indicating a heterocygote homocysteinemia.

We recommend each laboratory to establish its own reference range. The above mentioned values are only for orientation and may vary from other published data.

### *Controls*

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## **14. PERFORMANCE CHARACTERISTICS**

### *Precision and reproducibility*

#### **Intra-Assay CV**

2.5 % (8.9  $\mu\text{mol/l}$ ) [n = 12]

1.6 % (16.5  $\mu\text{mol/l}$ ) [n = 12]

#### **Inter-Assay CV**

5.2 % (8.7  $\mu\text{mol/l}$ ) [n = 12]

2.9 % (17.4  $\mu\text{mol/l}$ ) [n = 12]

### *Linearity*

up to 500  $\mu\text{mol/l}$

### *Detection limit*

0.6  $\mu\text{mol/l}$

### *Recovery*

98 %

## **15. DISPOSAL**

The mobile phase (MOPHA), reduction solution (REDSOL), internal standard (INT STD), and derivatisation solution (DER) must be disposed as non-halogenated solvents. The precipitation reagent (PREC) can be neutralized with NaOH to pH 7.0 and disposed as salt solution.

**Important:** Reaction will produce heat, be careful!

Please refer to the appropriate national guidelines.

## 16. TROUBLESHOOTING

Problem	Possible cause	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system

Problem	Possible cause	Solution
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

## 17. REFERENCES

1. Role of Plasma Homocyst(e)ine in Arterial Occlusive Diseases. (1994) *Clinical Chemistry*, Vol. **40** No. 6, pp. 857-858.
2. Association between Plasma homocysteine concentrations and extracranial carotid-artery stenosis. Selhub, J. et al. (1995). *The New England Journal of Medicine*, Feb. 95, Vol. **332** No. 5, pp. 286-291.
3. Plasma Homocyst(e)ine as a Risk Factor for Early Familial Coronary Artery disease. Wu, L. et al. (1994). *Clinical Chemistry*, Vol. **40** No. 4, pp. 552-561.
4. Homocyst(e)ine and arterial occlusive diseases. Malinow, M. (1994). *Journal of Intern Medicine*, Vol. **236**, pp. 603-617.
5. Hormone replacement therapy may reduce high serum homocysteine in postmenopausal women. van der Mooren, M. (1994). *European Journal of Clinical Investigation*, Vol. **24**, pp. 733-736.
6. Homocysteine metabolism in pregnancies complicated by neural-tube defects. Mills, J. et al. (1995). *The Lancet*, Jan. 1995, Vol. **345**, pp. 149-151.

## 18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date stated on kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

**Used symbols:**

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		