# IDKmonitor<sup>®</sup> Infiximab total ADA ELISA

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# **IDKmonitor<sup>®</sup> Infliximab total ADA ELISA**

***For the in vitro determination of total  
human antibodies against infliximab (e. g. REMICADE<sup>®</sup>)  
in EDTA plasma and serum***

Valid from 2019-08-21

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>6. PREPARATION OF THE ASSAY</b>	<b>17</b>
<i>Sample storage</i>	17
<i>Preparation of samples and controls</i>	17
<b>7. ASSAY PROCEDURE</b>	<b>18</b>
<i>Principle of the test</i>	18
<i>Test procedure</i>	18
<b>8. RESULTS</b>	<b>19</b>
<b>9. LIMITATIONS</b>	<b>20</b>
<i>Measurement range</i>	20
<i>Biotin interference</i>	20
<b>10. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Accuracy – Precision</i>	20
<i>Accuracy – Trueness</i>	21
<i>Linearity</i>	21
<i>Analytical sensitivity</i>	22
<i>Analytical specificity – cross-reactivity</i>	22
<i>Analytical specificity – Interferences</i>	22
<i>High-Dose-Hook effect</i>	22
<b>12. PRECAUTIONS</b>	<b>23</b>
<b>13. TECHNICAL HINTS</b>	<b>23</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>15. LITERATURE</b>	<b>24</b>

## 1. INTENDED USE

This Immundiagnostik AG assay is an enzyme immunoassay intended for the determination of human antibodies against TNF $\alpha$  blocker infliximab (e.g. REMICADE®) in the presence of infliximab in EDTA plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis or rheumatoid arthritis are often treated with anti-TNF $\alpha$  antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF $\alpha$  therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF $\alpha$  antibody [1–3]. Several factors influence the trough level, among them dosage and frequency of anti-TNF $\alpha$  blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation of anti-drug antibodies, ADA). It is thought that ADA functionally neutralise the therapeutic antibodies or induce their rapid elimination [4]. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF $\alpha$  antibody application [1, 5].

The IDKmonitor® Infliximab total ADA ELISA for the detection of total antibodies against infliximab (e.g. REMICADE®) measures free and bound antibodies against infliximab. This assay allows a reliable determination of ADA even in the presence of infliximab; therefore it is ideal for therapy monitoring when a measurable infliximab concentration is expected, for example shortly after the last infusion. In combination with the drug level determination, the IDKmonitor® Infliximab total ADA ELISA is an opportunity for the treating physician to monitor and optimise the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9654	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 9654	CONJ	Conjugate concentrate, peroxidase-labelled	1 x 600 $\mu$ l
K 9654	TRACER	Tracer concentrate, biotinylated	1 x 600 $\mu$ l
K 9654	CTRL NEG	Negative control, lyophilised (see specification for range)	4 x 1 vial

Cat. No.	Label	Kit components	Quantity
K 9654	CTRL POS	Positive control, lyophilised (see specification for range)	4 x 1 vial
K 9654	CTRL CUT-OFF	Cut-off control, lyophilised	4 x 1 vial
K 9654	ASYBUF	Assay buffer, ready-to-use	2 x 15 ml
K 9654	ABBUF	Antibody dilution buffer, ready-to-use	1 x 10 ml
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. PREPARATION AND STORAGE OF REAGENTS

- **Bring all reagents to room temperature (15–30 °C) prior to use.**
- To run the assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** should be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt con-

centration in the concentrate. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

- **Preparation of controls** see chapter 6.
- **Preparation of the conjugate and tracer: A few minutes before use, tracer concentrate (TRACER) and conjugate concentrate (CONJ) have to be diluted 1:12 in antibody dilution buffer (ABBUF)** (e.g. 3000 µl antibody dilution buffer + 300 µl tracer + 300 µl conjugate). TRACER and CONJ are stable at **2–8 °C** until expiry date stated on the label. **Tracer (1:12 diluted TRACER) and conjugate (1:12 diluted CONJ) are not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

## 6. PREPARATION OF THE ASSAY

### *Sample storage*

**Undiluted samples** can be stored for **2 months at -20 °C** and for **7 days at 2–8 °C or room temperature**. **Avoid more than 3 freeze-thaw-cycles.**

**Diluted samples** are **not stable** and **cannot be stored**.

### *Preparation of samples and controls*

1.	<p>Dilute samples <b>1:10</b> in <b>assay buffer (ASYBUF)</b> by pipetting <b>25 µl sample</b> into a reaction tube and adding <b>225 µl assay buffer</b>. Mix well. Addition of assay-buffer to all the samples should be performed without pause since this step dissociates the antibody–therapeutic antibody complexes.</p> <p>Reconstitute the <b>controls</b> with <b>250 µl assay buffer</b> and vortex. Carry out this step simultaneously with sample dilution in order to ensure equal treatment of controls and samples.</p>
2.	<p>Incubate controls and diluted samples in reaction tubes for <b>20 min</b> with <b>shaking</b> on a horizontal shaker at room temperature.</p> <p><b>CAUTION:</b> incubation time begins <b>upon addition of assay buffer</b>.</p>

3.	Add <b>60 µl tracer/conjugate/antibody dilution buffer solution</b> (see preparation of reagents) to <b>250 µl</b> control/diluted sample. Vortex and incubate for <b>1 hour</b> with <b>shaking</b> at room temperature (15–30 °C).
----	--

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA serves for the determination of antibodies against TNF $\alpha$ -Blocker infliximab (e.g. REMICADE®). During sample preparation, the anti-drug antibodies (ADA) are separated from the therapeutic antibody in order to acquire free ADA. By adding the peroxidase conjugate (peroxidase labelled therapeutic antibody) and the tracer (biotinylated therapeutic antibody), the unmarked therapeutic antibodies are replaced and the marked antibodies can form a complex with the ADA. This complex binds via biotin to the streptavidin coated microtiter plate. It is detected via the peroxidase conjugate with the peroxidase converting the substrate TMB to a blue product. The enzymatic reaction is stopped by adding an acidic solution. The samples convert from blue to yellow. The colour change should be measured in a photometer at 450 nm. The interpretation is made using the cutoff control.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of controls/samples on a protocol sheet.

Take as many microtiter strips (PLATE) as needed from kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until the expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	<b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Pipet <b>100 µl of preincubated controls/samples</b> into the wells of the microtiter plate.

3.	Cover the strips and incubate for <b>1.5 hours shaking*</b> on a horizontal shaker at room temperature (15–30 °C).
4.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5.	Add <b>100 µl substrate</b> (SUB) into each well.
6.	Incubate for <b>10–20 minutes</b> at room temperature (15–30 °C) until a sufficiently large difference in colour occurs**.
7.	Add <b>100 µl stop solution</b> (STOP) into each well and mix shortly by using the <b>shake function*</b> of the microplate reader.
8.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference.

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

Cut-off = 10 AU/ml = OD cut-off control

Samples which have a higher average optical density (OD) than the cut-off control are positive.

Immundiagnostik AG recommends linear regression using a linear ordinate and abscissa to calculate the results.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a manual control of the paired values should be made.

### Sample calculation for a positive sample

Average OD of patient's sample	0.735
Average OD of cut-off control	0.085 = 10 AU/ml
Concentration patient's sample	$\frac{0.735 \times 10 \text{ AU/ml}}{0.085} = 86.47 \text{ AU/ml}$

## 9. LIMITATIONS

### *Measurement range*

The lower limit of the measurement range is the LoB.

LoB see chapter "Performance Characteristics".

Samples with concentrations lower than the measurement range cannot be clearly quantified.

### *Biotin interference*

Samples containing a biotin concentration of < 100 ng/ml show a change of the results of  $\leq 25\%$ . Higher concentrations of biotin can lead to falsely low results. Patients taking > 5 mg biotin per day should wait at least 24 hours after taking biotin to have their samples collected. Results of patients taking biotin supplements or receiving a high-dose biotin therapy should generally be interpreted along with the total clinical picture.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n=41**

The repeatability was assessed with 2 serum samples under **constant** parameters (same operator, measurement system, day and kit lot).

Sample	Mean value [AU/ml]	CV [%]
1	54.08	3.2
2	219.54	4.7

### Reproducibility (Inter-Assay)

The reproducibility was assessed with 3 serum samples under **varying** parameters (different operators, measurement systems, days and kit lots).

Sample	Mean value [AU/ml]	CV [%]
1 (n=32)	56.02	8.8
2 (n=32)	62.15	9.2
3 (n=20)	9.49	6.1

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, known amounts of anti-Infliximab antibodies were added to 4 serum samples and measured. The recovery was found between 80.01 and 118.18%.

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A by two serial dilutions of 2 serum samples.

For anti-infliximab antibodies in serum and plasma, the method has been demonstrated to be linear from 14.23 to 497.00 AU/ml, showing a non-linear behaviour of less than  $\pm 20\%$  in this interval. The following table shows one representative example for each sample:

Sample	Dilution	Expected [AU/ml]	Obtained [AU/ml]	Recovery [%]
1	1:2	497.00	497.00	100.00
	1:4	248.50	247.07	99.42
	1:8	124.25	129.74	104.42
	1:16	62.13	63.00	101.41
	1:32	31.06	32.71	105.31
	1:64	15.53	16.35	105.25

Sample	Dilution	Expected [AU/ml]	Obtained [AU/ml]	Recovery [%]
2	1:2	370.48	370.48	100.00
	1:4	185.24	199.17	107.52
	1:8	92.62	98.74	106.61
	1:16	46.31	55.37	119.57
	1:32	23.16	28.84	124.56

### *Analytical sensitivity*

Limit of blank, LoB	3.12 AU/ml
Limit of detection, LoD	6.93 AU/ml
Limit of quantitation, LoQ	10 AU/ml

The evaluation was performed according to the CLSI guideline EP17-A2. The specified accuracy goal for the LoQ was 10% CV.

### *Analytical specificity – cross-reactivity*

The specificity of the antibody was tested by measuring the cross-reactivity against anti-Adalimumab antibodies in patient samples. Therefore, 3 patient samples, which have previously been tested positive on anti-adalimumab antibodies, were tested. There was no cross-reactivity observed.

### *Analytical specificity – Interferences*

Different substances that might interfere when using the IDKmonitor® Infliximab total ADA ELISA K 9654 were tested. Therefore, positive as well as negative serum samples were supplemented with either drugs (maximum daily dose) or serum components (doses according to CLSI guideline EP7-A2) and then measured.

No interferences with the following drugs were found: Azathioprin, Pantoprazol, Me-salazin, Clarithromycin, Levofloxacin, iron supplements, acetylsalicylic acid, vitamin D, Gabapentin, multivitamin preparation or Ibuprofen.

No interferences with the following serum components were found: hemoglobin, bilirubin or triglycerides.

### *High-Dose-Hook effect*

No High-Dose-Hook effect was observed for antibody concentrations up to 981 AU/ml.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDKmonitor® is a trademark of Immundiagnostik AG.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

## 15. LITERATURE

1. Afif W, Loftus EV, Jr., Faubion WA, Kane SV, Bruining DH, Hanson KA, Sandborn WJ (2010). Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *Am J Gastroenterol* **105**(5): 1133-1139.
2. Kopylov U, Mazor Y, Yavzori M, Fudim E, Katz L, Coscas D, Picard O, Chowers Y, Eliakim R, Ben-Horin S (2012). Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-infliximab antibodies. *Inflamm Bowel Dis* **18**(9): 1628-1633.
3. Tak PP (2012). A personalized medicine approach to biological treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology* **51**(4): 600-609.
4. Ordas I, Mould DR, Feagan BG, Sandborn WJ (2012) Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clin Pharmacol Ther* **91**(4): 635-646.
5. Bender NK, Heilig CE, Droll B, Wohlgemuth J, Armbruster FP, Heilig B (2007). Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatol Int* **27**(3): 269-274.

**Used symbols:**

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

**Immundiagnostik AG**

Stubenwald-Allee 8a  
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)