# **IDKmonitor<sup>®</sup> Vedolizumab drug level LC-MS/MS Kit**

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

# **IDKmonitor<sup>®</sup> vedolizumab drug level LC-MS/MS Kit**

*For the in vitro determination of vedolizumab  
concentration (e. g. Entyvio<sup>®</sup>) in EDTA plasma and serum*

Valid from 2017-05-10

**REF** **KM9600**

**RUO**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

1. INTENDED USE _____	15
2. INTRODUCTION _____	15
3. MATERIAL SUPPLIED _____	15
4. MATERIAL SUPPLIED WITH THE VEDOLIZUMAB <i>ImmuTube</i> ® EXTRACTION KIT	16
5. MATERIAL SUPPLIED WITH THE VEDOLIZUMAB <i>ImmuTube</i> ® DIGEST KIT _____	16
6. MATERIAL REQUIRED BUT NOT SUPPLIED _____	17
7. STORAGE AND PREPARATION OF REAGENTS _____	17
8. SAMPLE PREPARATION WITH THE VEDOLIZUMAB <i>ImmuTube</i> ® EXTRACTION KIT FOR 80 SAMPLES _____	18
9. PEPTIDE PREPARATION WITH THE VEDOLIZUMAB <i>ImmuTube</i> ® DIGEST KIT FOR 80 SAMPLES _____	20
10. CHROMATOGRAPHIC CONDITIONS _____	21
11. MS/MS METHOD _____	22
12. CALCULATION _____	23
13. QUALITY CONTROL _____	23
14. PRECAUTIONS _____	23
15. TECHNICAL HINTS _____	23
16. DISPOSAL _____	23
17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE _____	24

## 1. INTENDED USE

This LC-MS/MS application is for the quantitative determination of the human  $\alpha_4\beta_7$  integrin therapeutic antibody vedolizumab (e.g. Entyvio®) in EDTA plasma and serum. For research use only. Not for use in diagnostic procedures.

## 2. INTRODUCTION

Patients with moderately to severely active Colitis ulcerosa and Crohn's disease, who have had an inadequate response to conventional therapy or to anti-tumor necrosis factor alpha (TNF $\alpha$ ) agents, can also be treated with vedolizumab. The humanised monoclonal therapeutic antibody vedolizumab binds  $\alpha_4\beta_7$  integrin on activated lymphocytes and stops them from migrating into the intestinal mucosa. Thus, vedolizumab suppresses the inflammatory response via a different mechanism than anti-TNF $\alpha$  agents and specifically targets inflammation in the gastrointestinal tract.

The clinical efficacy of vedolizumab therapy correlates with the trough level, that is the drug level just before the next application of vedolizumab. Several factors influence the trough level, among them dosage and frequency of vedolizumab infusion, disease activity, individual pharmacokinetics and immune reaction (formation of anti-drug antibodies, ADA). The IDKmonitor®LC-MS/MS application for the determination of vedolizumab (e.g. Entyvio®) measures reliably the effective drug concentration, giving the treating physician an opportunity to monitor and optimise the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KM9600LA	MOPHA A	Mobile phase A	1 x 250 ml
KM9600LB	MOPHA B	Mobile phase B	1 x 250 ml
KM9600KA	CAL 1 CAL 2 CAL 3	Calibrator 1 to 3, lyophilised (Concentrations and reconstitution details are given in the specification data sheet.)	5 x 3 vials each
KM9600KO	CTRL 1 CTRL 2	Control 1 and 2, lyophilised (Ranges and reconstitution details are given in the specification data sheet.)	5 x 2 vials each
KM9600RE	RECSOL	Reconstitution solution	3 x 15 ml
KM9600AC	ACTSOL	Activation reagent	1 x 2 ml

Cat. No.	Label	Kit components	Quantity
KM9600IS	INT STD	Internal standard, lyophilised (Concentrations and reconstitution details are given in the specification data sheet.)	5 x 1 vial
KM9600SO	SOL A	Solution A	2 x 25 ml
KM9610		Vedolizumab <i>ImmuTube</i> ® Extraction Kit	1 Kit
KM9620		Vedolizumab <i>ImmuTube</i> ® Digest Kit	1 Kit

As a first step for the application of the IDKmonitor® LC-MS/MS application for the determination of the therapeutic  $\alpha_4\beta_7$  integrin antibody vedolizumab, a tuning is necessary to estimate the optimal LC-MS/MS settings as well as to assess the sufficiency of the sensitivity. The UPLC separation column (KM9600RP) can be ordered separately from Immundiagnostik AG. Please ask for our single components price list.

#### 4. MATERIAL SUPPLIED WITH THE VEDOLIZUMAB *ImmuTube*® EXTRACTION KIT

Cat. No.	Label	Kit components	Quantity
KM9610SI	COLUMNS	<i>ImmuTube</i> ® columns for affinity purification of vedolizumab from the samples	80 columns
KM9610WL	WASHSOL	Wash buffer	1 x 300 ml
KM9610EL	ELUREAG	Elution reagent	1 x 20 ml

#### 5. MATERIAL SUPPLIED WITH THE VEDOLIZUMAB *ImmuTube*® DIGEST KIT

Cat. No.	Label	Kit components	Quantity
KM9620B	BUF ABC	Buffer ABC concentrate	5 x 1 vial à 2 ml
KM9620D	BUF D	Buffer D, lyophilised	5 x 1 vial
KM9620IA	BUF IA	Buffer IA, lyophilised	5 x 1 vial
KM9620ER	ENZSOL	Enzyme reagent	1 x 1,5 ml
KM9620E	ENZ	Enzyme, lyophilised	1 vial
KM9620SL	STOP	Stop solution	1 x 1,5 ml

## 6. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Glass tubes, suitable for LC-MS/MS
- 500 ml graduated cylinder, suitable for LC-MS/MS
- Calibrated precision pipettors and 10–1000 µl tips
- Multi-channel pipets or repeater pipets
- Centrifuge, 1 000 g
- Vortex
- Vacuum centrifuge or nitrogen distributor
- Standard single-use laboratory plastic vials, cups, etc.
- Thermal block, heatable up to at least 60 °C
- Incubator
- Upside-down shaker
- LC-MS/MS equipment
- RP-C18 column for peptide analysis, (e.g. Acquity BEH130 C18, 1,7 µm (2,1 × 50 mm))

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 7. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 5 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of mobile phases:** The **mobile phases (MOPHA A, MOPHA B)** have to be spiked with **0.2 % activation reagent (ACTSOL)** before use (e.g. 250 ml MOPHA A/MOPHA B + 0,5 ml ACTSOL), mix well. **Hint:** If possible, use glass pipets to pipet ACTSOL. The **mobile phases** are stable at **2–8 °C** until the expiry date stated on the label. **Mobile phases spiked with activation reagent** can be stored for four weeks in a closed flask at **2–8 °C**.  
**Attention:** The **activation reagent** must be added under the **fume hood**. All vials to be used must be absolutely clean, detergent-free and preferably made of glass suitable for LC-MS/MS.
- The **lyophilised controls (CTRL)** and **calibrators (CAL)** are stable at **-20 °C** until the expiry date stated on the label. Before use, the standards and con-

trols have to be reconstituted with each **100 µl of reconstitution solution (RECSOL)**. Mix the vial content well and allow it to dissolve for 10 minutes to ensure complete reconstitution. **Controls and calibrators** (reconstituted CTRL and CAL) **are not stable and cannot be stored.**

- The **lyophilised internal standard** (INT STD) is stable at **-20 °C** until the expiry date stated on the label. Before use, the INT STD has to be reconstituted with **solution A** (SOL A; see specification data sheet for volume). Mix the vial content well and allow it to dissolve for 10 minutes to ensure complete reconstitution. **Internal standard** (reconstituted INT STD) **is not stable and cannot be stored.**

**Attention: Do not prepare internal standard until directly before use!**

- **Preparation of buffers:**

- The **buffer concentrate ABC** (BUF ABC) is stable at **4 °C** until the expiry date stated on the label. Before use, the BUF ABC has to be diluted **1:20 with reconstitution solution** (RECSOL), mix well.
- The **lyophilised buffers D** (BUF D) and **IA** (BUF IA) are stable at **-20 °C** until the expiry date stated on the label. Before use, they have to be reconstituted with each **500 µl of buffer ABC**, mix well. **Buffers D and IA** (reconstituted BUF D and BUF IA) **are not stable and cannot be stored.**
- The **lyophilised enzyme (ENZ)** is stable at **-20 °C** until the expiry date stated on the label. **Directly before use**, it is reconstituted with **1 ml of enzyme reagent** (ENZSOL). **Enzyme** (reconstituted ENZ) **can be stored at -20 °C for 6 months.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

## **8. SAMPLE PREPARATION WITH THE VEDOLIZUMAB ImmuTube® EXTRACTION KIT FOR 80 SAMPLES**

Serum and plasma samples can be used with this assay. All samples must be vortexed before use in the test.

Quality control samples should be analysed with each run.

1.	Bring <b>all reagents and samples</b> to <b>room temperature</b> (18–26 °C) and mix well.
2.	Put the <i>ImmuTubes</i> ® into an appropriate rack. Assure that the included suspension does not stick to the lid. Label the <i>ImmuTubes</i> ®.
3.	Unscrew the lids of the <i>ImmuTubes</i> ® and <b>open</b> the outlet by breaking the endpiece at the predetermined breaking point. Put the <i>ImmuTubes</i> ® into appropriate tubes. <b>Wash</b> the matrix (add 500 µl wash buffer (WASHSOL) and centrifuge for 30 s at 1000 <i>g</i> ). <b>Repeat</b> the washing step <b>twice</b> .
4.	Put <b>500 µl wash buffer</b> (WASHSOL) into the <i>ImmuTubes</i> ®. <b>Close</b> the <i>ImmuTubes</i> ® firmly by recapping the outlet with the upside-down endpiece. Speedily add <b>25 µl calibrator/sample/control</b> (CAL/SAMPLE/CTRL). <b>Close</b> the <i>ImmuTubes</i> ® by screwing the lid on and <b>mix</b> the content carefully by gentle inversion.
5.	Mix for <b>1 hour</b> by slow <b>end-over-end rotation</b> at room temperature. After the incubation, shortly centrifuge the columns to collect matrix remains in the lid.
6.	Unscrew the lid of the <i>ImmuTubes</i> ® and open the outlet. Put the <i>ImmuTubes</i> ® into plastic test tubes and centrifuge for 30 s at 1000 <i>g</i> . Add <b>500 µl wash buffer</b> (WASHSOL) and centrifuge for 30 s at 1000 <i>g</i> . <b>Repeat</b> this step <b>once with wash buffer</b> and <b>three times with ultra pure water</b> . Discard flow-through if necessary.
7.	Place <i>ImmuTubes</i> ® in already labelled new 1.5 ml test tubes (Immundiagnostik AG recommends Protein LowBind Tubes from Eppendorf). Add <b>100 µl ELUREAG</b> (elution reagent for <i>ImmuTubes</i> ®), centrifuge for 30–60 s at 1000 <i>g</i> and collect the eluate containing the vedolizumab in the test tube. <b>Repeat</b> this step once with <b>50 µl ELUREAG</b> , collect the eluate in the same test tube.
8.	Dry the <b>eluate</b> under a nitrogen stream at 37 °C or in a vacuum centrifuge.

Dried eluate is stable for two weeks at -20 °C.

## 9. PEPTIDE PREPARATION WITH THE VEDOLIZUMAB *ImmuTube*® DIGEST KIT FOR 80 SAMPLES

Before each incubation, spin the tubes shortly to pool all reagents at the bottom of the tube.

1.	All solutions for the enzymatic reaction should always be <b>prepared directly before use</b> .
2.	Reconstitute the dried eluate in <b>20 µl buffer D</b> (BUF D) by intense vortexing. If necessary, incubate for 5 min in an ultrasonic bath.
3.	Incubate for <b>1 hour</b> in a thermal block at <b>60°C</b> .
4.	Add <b>20 µl buffer IA</b> (BUF IA), vortex shortly and <b>shake gently</b> for <b>30 min</b> at room temperature <b>in the dark</b> .
5.	Add <b>100 µl buffer ABC</b> (BUF B) to the samples.
6.	Now <b>reconstitute</b> the <b>lyophilised enzyme</b> in <b>1 ml enzyme reagent</b> (ENZ-SOL) at <b>2–8°C</b> .
7.	<b>Preparation of the enzyme solution:</b> For each sample, fill the reconstituted <b>enzyme</b> up to a <b>final volume of 50 µl</b> with <b>buffer ABC</b> (BUF ABC), e.g. per sample 10 µl enzyme solution + 40 µl buffer ABC (BUF ABC). Aliquot the remaining enzyme and store it immediately at -20°C until the next usage.
8.	Add <b>50 µl enzyme solution</b> to each sample and incubate in an incubator or thermal block at <b>37°C over night</b> (~ 16 h).
9.	Stop the reaction by adding <b>10 µl stop solution</b> (STOP).
10.	Reconstitute lyophilised <b>internal standard</b> (INT STD) in solution A and add <b>200 µl</b> to each sample (final volume: 400 µl). Vortex shortly and put the samples in vials suitable for LC-MS/MS. <b>Inject 25 µl</b> into the HPLC system.

## 10. CHROMATOGRAPHIC CONDITIONS

Example of a comparable UPLC method

**Column material:** e.g. Acquity UPLC® BEH130 C18 (Waters); 1.7 µm

**Column dimension:** 2.1 × 50 mm

**Flow rate:** 0.1 ml/min

**Column temperature:** 45°C

**Injection volume:** 25 µl

**Running time:** 25 min

**Gradient:**

Time	Mobile phase A	Mobile phase B
0 min	97 % A	3 % B
5.1 min	97 % A	3 % B
15.1 min	55 % A	45 % B
15.2 min	20 % A	80 % B
17.0 min	20 % A	80 % B
17.1 min	97 % A	3 % B
25.0 min	97 % A	3 % B

It is recommended to use a guard column/filter to extend the column's life.

After the analysis, the separation column can be washed with ~ 50 ml MOPHA A using the above mentioned gradient. The column can be stored in 80 % acetonitrile.

## 11. MS/MS METHOD

The MS/MS method shown below is an example for a TSQ Vantage Triple Quadrupol-mass spectrometer (Thermo Fisher Scientific)\*.

<b>Mode:</b>	SRM
<b>Polarity:</b>	ESI+
<b>Capillary temperature [°C]:</b>	350
<b>Vaporizer temperature [°C]:</b>	400
<b>Sheath gas pressure [Arb]:</b>	60
<b>Aux gas pressure [Arb]:</b>	15
<b>Spray voltage [V]:</b>	4000
<b>Resolution in 1st and 3rd Quadrupol [FWHM]:</b>	0,7
<b>Collision gas [mTorr]:</b>	1,5
<b>Cycle time [s]:</b>	0.5–1

\* These settings are sample settings and can be different depending on the manufacturer and the machine used. Before analysis, tuning is necessary to estimate the optimal LC-MS/MS settings as well as to assess the sufficiency of the sensitivity.

*SRM transitions (m/z)*

**Vedolizumab, peptide for quantification: MH1+: 2237.0357 (MH2+: 1119.022)**

Transitions			
Precursor ion	Fragment ion	S-LENS voltage	Collision energy
1119.022	429.213	231	37
1119.022	956.472	231	28
1119,022	728.340	231	28
1119.022	1281.571	231	33

**Internal Standard (isotopic labelled peptide: MH1+: 2245.28 (MH2+: 1123.029))**

Transitions			
Precursor ion	Fragment ion	S-LENS voltage	Collision energy
1123.029	429.213	231	37
1123.029	956.472	231	28
1123.029	1001.965	231	26
1123.029	1289.585	231	33

## 12. CALCULATION

Linear regression can be used as model for evaluation of the results. The three calibrator concentration points are connected by a straight line. The samples can be calculated using this obtained line.

## 13. QUALITY CONTROL

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## 14. PRECAUTIONS

- For research use only
- Quality control samples should be analysed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

## 15. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date stated on kit label.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according the enclosed manual.

## 16. DISPOSAL

Mobile phase (MOPHA A, MOPHA B), solution A (SOL A), activation reagent (ACTSOL), reconstitution solution (RECSOL) and stop solution (STOP) have to be disposed as non-halogenated solvents.

## 17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the kit package are for research use only.
- The guidelines for laboratories should be followed.
- *IDKmonitor*® and *ImmuTube*® are trademarks of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

### Used symbols:

	Temperature limitation		Catalogue Number
	For research use only		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use