Malondialdehyd HPLC Kit

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
Malondialdehyde HPLC Kit

*For the determination of malondialdehyde
in plasma, serum and urine*

Valid from 2019-02-12

KC1900

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

| | |
|---|-----------|
| 1. INTENDED USE | 17 |
| 2. SUMMARY AND EXPLANATION OF THE TEST | 17 |
| 3. PRINCIPLE OF THE TEST | 18 |
| <i>Summary</i> | 18 |
| 4. MATERIAL SUPPLIED | 18 |
| 5. MATERIAL REQUIRED BUT NOT SUPPLIED | 19 |
| 6. PREPARATION AND STORAGE OF REAGENTS | 19 |
| 7. SPECIMEN COLLECTION AND PREPARATION | 19 |
| 8. ASSAY PROCEDURE | 20 |
| <i>Procedural notes</i> | 20 |
| <i>Sample preparation</i> | 20 |
| <i>Chromatographic conditions</i> | 21 |
| 9. TREATMENT OF THE COLUMN | 21 |
| 10. RESULTS | 21 |
| <i>Calculation</i> | 21 |
| <i>Typical chromatogram</i> | 22 |
| 11. LIMITATIONS | 22 |
| 12. QUALITY CONTROL | 22 |
| <i>Expected values</i> | 22 |
| <i>Controls</i> | 22 |
| 13. PERFORMANCE CHARACTERISTICS | 23 |
| <i>Accuracy – Precision</i> | 23 |
| <i>Accuracy – Trueness</i> | 23 |
| <i>Analytical specificity</i> | 24 |
| 14. DISPOSAL | 25 |
| 15. TROUBLESHOOTING | 25 |
| 16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE | 26 |
| 17. REFERENCES | 27 |

1. INTENDED USE

The Immundiagnostik AG assay is intended for the quantitative determination of malondialdehyde in plasma, serum and urine. This assay is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

In the last years the damage effects of lipid peroxidation products were studied intensively. These are formed when free radicals will overcome the radical-trapping-mechanisms of the body and reacting with unsaturated fatty acids.

The reaction of polyunsaturated fatty acids (PUFA's) with activated oxygen species results in lipidhydroperoxides (primary lipid peroxidation products) which were degraded to secondary lipid peroxidation products like alkanes (e.g. ethane and pentane), aliphatic aldehydes (e.g. malondialdehyde [MDA] and 4-hydroxynonenale [4-HNE]).

Primary and secondary lipid peroxidation products have an influence on a lot of molecules responsible for correct cell function.

So lipidhydroperoxides easily pass the nuclear membrane and can react with nucleic acids. Also proteins can be attacked on their thiol groups changing their stereochemistry and function.

Moreover lipidhydroperoxides interfere with chemical and physical properties of the cell membrane. The fluidity decreases and rigidity increases. The so influenced cell membrane can't maintain their barrier function and intracellular potassium ions leak out as well as intracellular enzymes are lost. If erythrocytes are afflicted haemolysis takes place. In this case haemoglobin can initiate or propagate the lipid peroxidation.

Secondary lipid peroxidation products like MDA or 4-HNE can react with DNA as well, in particular with the bases guanin and adenin. These DNA aberrations lead to erroneous transcriptions and therewith to altered gene products. Peptide bonds are broken through the impact of MDA in proteins. Aldehydes react with amino groups of proteins building Schiff-bases, elementary disturbing correct function of proteins.

All these toxic features of oxidised fatty acids are discussed in the pathogeneses of many diseases and dysfunction of organs. In particular, arteriosclerosis, tumor diseases, rheumatic diseases and reperfusion damage of organs after ischämische processes.

3. PRINCIPLE OF THE TEST

The first step in determining malondialdehyde is a sample preparation with derivatisation. The derivatisation reagent transforms malondialdehyde into a fluorescent product. Afterwards the pH is optimised through to addition of a reaction solution. 20 µl of the supernatant are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 4 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak heights by the external standard method.

Summary

Besides many other parameters the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is mostly sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

4. MATERIAL SUPPLIED

| Cat. No. | Label | Kit components | Quantity |
|----------|----------------|--|-------------|
| KC 1900 | MOPHA | Mobile phase | 1 000 ml |
| KC 1900 | CAL | Calibrator (lyophilised 0.25 ml, concentration see label) | 5 vials |
| KC 1900 | DER | Derivatisation solution | 100 ml |
| KC 1900 | REABUF | Reaction solution | 50 ml |
| KC1900 | CTRL1 CTRL2 | Control 1 and 2 (lyophilised 0,25 ml; concentration see specification) | 2 x 3 vials |

For reorders of single components, use the catalogue number followed by the label as product number.

HPLC column (Cat.No.: KC1900, label: Column) as well as individual components can be ordered separately from Immundiagnostik AG. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Vortex
- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C18-column - Bischoff Prontosil Eurobond C₁₈ 5 µm, 125 x 4 mm
- Water bath or heating block for heating at 95 °C

* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** (CAL) in **250 µl** aqua bidist. Take aliquots and store at -20 °C. Reconstituted calibrator is stable for at least 2 weeks at -20 °C. Avoid repeated freeze-thaw circles. The concentration of malondialdehyde might have minor changes from lot to lot.
- Reconstitute **controls** (CTRL1, CTRL2) in **250 µl** aqua bidist. **Storage** details are given in the **specification data sheet**.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

7. SPECIMEN COLLECTION AND PREPARATION

Plasma, serum and urine are suited for this test system. After venipuncture the sample should be stored at 2–8 °C immediately. Samples are stable for at least 24 h at 2–8 °C and 2 weeks at -20 °C.

8. ASSAY PROCEDURE

Procedural notes

- Control samples should be analysed with each run.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according to the enclosed manual.
- Calibrator (CAL) and controls (CTRL1, CTRL2) contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The derivatisation solution (DER) contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.

Sample preparation

Important: work with a water blank as a fluorescent compound is produced from the derivatisation solution (DER) that has the same retention time as the MDA-derivatisation product. The obtained blank value must be subtracted from all preparations.

Pipet each **20 µl water** (blank), **sample, calibrator** (CAL) **or control 1 or 2** (CTRL1, CTRL2) into a 1.5 ml reaction tube.

Add each **1 ml derivatisation solution** (DER) and vortex for **15 seconds**.

Incubate for **60 min** at 95 °C.

Keep incubation time and temperature constant as only at these conditions the given MDA-concentrations for calibrator (CAL) and controls (CTRL1, CTRL2) are valid.

It is possible that the lids of the reaction tube open during incubation. This can be avoided by putting a weight on the reaction tubes.

Cool down the solution (**min. 15 min at 2–8 °C**) and centrifuge for 5 min.

Take **500 µl** supernatant + **500 µl** reaction solution (REABUF) and mix thoroughly on a vortex mixer.

Inject **20 µl** of the mixture into the HPLC system.

The derivatised sample is stable at 2–8 °C for at least 4 days and at room temperature for at least 12 hours.

Chromatographic conditions

| | | | |
|--------------------------|--|------------|--------|
| Important: | No recirculation of the eluent is allowed for this test-system | | |
| Column material: | Bischoff Prontosil Eurobond, 5 µm | | |
| Column dimension: | 125 mm x 4 mm | | |
| Flow rate: | 0.8–1.2 ml/min | | |
| Detection: | Fluorescence: | Excitation | 515 nm |
| | | Emission | 553 nm |
| Temperature: | 30 °C | | |
| Injection volume: | 20 µl | | |
| Running time: | 4 min | | |

It is recommended that a guard column is used to extend column life.

9. TREATMENT OF THE COLUMN

After the analysis the column should be flushed with 30 ml ultrapure water (1 ml/min) and stored in 50 % methanol in aqua bidest (ca. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. **30 ml MOPHA** (mobile phase).

Important: Do not recirculate the MOPHA (mobile phase) in this test system.


10. RESULTS

Calculation

Take into account that the mentioned heights of calibrator and patient are the heights from which the blank value has been subtracted.

$$\text{Concentration sample} = \frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height calibrator}}$$

Typical chromatogram

11. LIMITATIONS

Do not use whole blood.

Strong haemolytic and lipaemic samples often show pathological concentrations. Do not measure such samples.

12. QUALITY CONTROL

Expected values

Serum, heparin plasma: $1.97 \pm 0.41 \mu\text{mol/l}$

EDTA plasma: $0.61 \pm 0.24 \mu\text{mol/l}$

It is recommended that each laboratory should establish its own normal range. Above mentioned values are only for orientation and may vary from other published data.

Controls

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the

same assay one or more values of the quality control sample are outside the acceptable limits.

13. PERFORMANCE CHARACTERISTICS

Accuracy – Precision

Repeatability (Intra-Assay); n = 12

The repeatability was assessed with one control sample under **constant** parameters (same operator, measurement system, day and kit lot).

| Mean value [$\mu\text{mol/l}$] | CV [%] |
|----------------------------------|--------|
| 5.98 | 4.0 |

Reproducibility (Inter-Assay); n = 21

The reproducibility was assessed with 2 control samples under **varying** parameters (different operators, measurement systems, days and kit lots).

| Sample | Mean value [$\mu\text{mol/l}$] | CV [%] |
|--------|----------------------------------|--------|
| 1 | 2.58 | 9.3 |
| 2 | 5.51 | 9.7 |

Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, three serum samples and three control samples were spiked with known amounts of Malondialdehyde three times. The mean values are shown in the table below.

The recovery for Malondialdehyde was found between 96.38 bis 102.95 %.

| Sample | Spike [$\mu\text{mol/l}$] | Expected [$\mu\text{mol/l}$] | Obtained [$\mu\text{mol/l}$] | Recovery [%] |
|--------|-----------------------------|--------------------------------|--------------------------------|--------------|
| 3.62 | 5 | 8.62 | 8.79 | 101.92 |
| 4.66 | 5 | 9.66 | 9.79 | 101.38 |
| 3.16 | 5 | 8.16 | 8.39 | 102.76 |
| 2.50 | 5 | 7.50 | 7.23 | 96.38 |
| 4.09 | 5 | 9.09 | 9.17 | 100.92 |
| 1.95 | 5 | 6.95 | 7.16 | 102.95 |

Lower detection limit

The limit of detection (LoD) is defined as 3 times the background noise. It is calculated by the formula below:

$$\frac{(3 \times \text{peak height background noise}) \times \text{concentration calibrator } [\mu\text{mol/l}]}{\text{Peak height calibrator}} = \text{LoD } [\mu\text{mol/l}]$$

LoD Malondialdehyde: 0.064 $\mu\text{mol/l}$

The limit of quantitation (LoQ) is defined as 10 times the background noise. It is calculated by the formula below:

$$\frac{(10 \times \text{peak height background noise}) \times \text{concentration calibrator } [\mu\text{mol/l}]}{\text{Peak height calibrator}} = \text{LoQ } [\mu\text{mol/l}]$$

LoQ Malondialdehyde: 0.213 $\mu\text{mol/l}$

Upper detection limit & linearity

The upper limit of detection states up to which concentration a method results in a linear signal. Therefore, 6 samples with known concentrations were measured two times.

Between the range from 2.5 to 105 $\mu\text{mol/l}$ non-linearity was found below 20%.

| Sample [$\mu\text{mol/l}$] | Signal height Experiment 1 | Signal height Experiment 2 |
|---------------------------------|-------------------------------|-------------------------------|
| 2.5 | 8805 | 10862 |
| 10 | 34561 | 33645 |
| 15 | 56909 | 54170 |
| 25 | 102935 | 97879 |
| 55 | 249409 | 236687 |
| 105 | 489665 | 474167 |

It should be noted that the detection limits depend on the instrument as well as on the application.

Analytical specificity

There were found no interferences to other blood components.

14. DISPOSAL

The mobile phase (MOPHA), derivatisation solution (DER) and reaction solution (RE-ABUF) must be disposed as non-halogenated solvent.

Please refer to the appropriate national guidelines.

15. TROUBLESHOOTING

| Problem | Possible reason | Solution |
|--------------------------|--|---|
| No signal | No or defect connection to evaluation system | Check signal cord and connection |
| | Detector lamp is altered | Change lamp |
| No peaks | Injector is congested | Check Injector |
| Doublepeaks | Dead volume in fittings and / or column | Renew fittings and / or column |
| Contaminating peaks | Injector dirty | Clean injector |
| | Contamination at the head of the column | Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase |
| | Air in the system | Degas pump |
| | Autosampler vials contaminated | Use new vials or clean them with methanol |
| Broad peaks, tailing | Precolumn / column exhausted | Use new precolumn / column |
| Variable retention times | Drift in temperature | Use a column oven |
| | Pump delivers imprecise | Check pump, degas the system |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |

| Problem | Possible reason | Solution |
|------------------------|---|--------------------------------------|
| Baseline is drifting | Detector lamp did not yet reach working temperature | Wait |
| | Detector lamp is too old | Renew lamp |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |
| | Pump delivers imprecise | Check pump, degas the system |
| Baseline is not smooth | Pump delivers imprecise | Check pump, degas the system |
| | Detector flow cell is dirty | Clean flow cell |

16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Do not mix plugs and caps from different reagents.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

17. REFERENCES

1. Draper H.H., Hadley M. (1990). A review of recent studies on the metabolism of exogenous and endogenous malondialdehyde. *xenobiotika* **20**; 9; 901-907.
2. Griesmacher A. et al. (1995). Enhanced serum levels of thiobarbituric-acid-reactive substances in diabetes mellitus. *Am J Med* **98**; 469-475.
3. Valenzuela A (1990). The biological significance of malondialdehyde determination in the assessment of tissue oxidative stress. *Life sciences* **48**; 301-309

Used symbols:

| | | | |
|---|------------------------------------|--|-----------------------------------|
|  | Temperature limitation |  | Catalogue Number |
|  | In Vitro Diagnostic Medical Device |  | To be used with |
|  | Manufacturer |  | Contains sufficient for <n> tests |
|  | Lot number |  | Use by |
|  | Attention |  | Consult instructions for use |
|  | Consult specification data sheet | | |

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com