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# Tyrosin ELISA Kit

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---

|   |           |
|---|-----------|
| <b>1. INTENDED USE</b>                                  | <b>16</b> |
| <b>2. PRINCIPLE OF THE TEST</b>                         | <b>16</b> |
| <b>3. MATERIAL SUPPLIED</b>                             | <b>17</b> |
| <b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>17</b> |
| <b>5. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>18</b> |
| <b>6. PRECAUTIONS</b>                                   | <b>19</b> |
| <b>7. SPECIMEN COLLECTION AND PREPARATION</b>           | <b>19</b> |
| <b>8. ASSAY PROCEDURE</b>                               | <b>20</b> |
| <i>SAMPLE PREPARATION PROCEDURE</i>                     | 20        |
| <i>TEST PROCEDURE</i>                                   | 21        |
| <b>9. EVALUATION OF RESULTS</b>                         | <b>22</b> |
| <i>EXPECTED RESULTS</i>                                 | 24        |
| <b>10. PERFORMANCE CHARACTERISTICS</b>                  | <b>24</b> |
| <i>CROSS REACTIVITY</i>                                 | 24        |
| <i>PRECISION AND REPRODUCIBILITY</i>                    | 24        |
| <i>SENSITIVITY</i>                                      | 25        |
| <i>RECOVERY</i>   | 25        |
| <i>LINEARITY</i>  | 25        |
| <b>11. LIMITATIONS</b>                                  | <b>26</b> |
| <b>12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>26</b> |

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Hersteller

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# Tyrosine ELISA Kit



*For the determination of tyrosine in human EDTA plasma and serum*

**For research use only**

Valid from 25.06.2012

K 7015

## 1. INTENDED USE

This ELISA Kit is intended for the quantitative determination of tyrosine in human EDTA plasma and serum. It is for research use only.

## 2. PRINCIPLE OF THE TEST

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for tyrosine derivatization.

Afterwards, the treated samples and a polyclonal tyrosine antiserum are incubated in wells of a microtiter plate coated with tyrosine-derivative (tracer). During the incubation period, the target tyrosine in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The tyrosine in the sample displaces the antibodies out of the binding to the tracer. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the tyrosine concentration in the sample.

During the second incubation step, a peroxidase conjugate is added to each microtiter well to detect the anti-tyrosine antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the tyrosine concentration in the sample; this means, high tyrosine concentration in the sample reduces the concentration of tracer-bound antibody and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Tyrosine present in the patient samples is determined directly from this curve.

### 3. MATERIAL SUPPLIED

| Catalog No           | Content          | Kit Components                             | Quantity        |
|----------------------|------------------|--|-----------------|
| K7015MTP             | PLATE            | One holder with precoated strips           | 12 x 8 wells    |
| K7015ST              | STD              | Standards (0, 3, 10, 30, 100, 300 $\mu$ M) | 6 x 200 $\mu$ l |
| K7015KO1<br>K7015KO2 | CTRL 1<br>CTRL 2 | Controls                                   | 2 x 200 $\mu$ l |
| K7015WP              | WASHBUF          | Wash buffer concentrate (10-fold)          | 2 x 100 ml      |
| K7015AK              | AB               | Anti-tyrosine antibody (lyo.)              | 2 x 1 vial      |
| K7015K               | 2.AB             | POD antibody (concentrate 200-fold)        | 65 $\mu$ l      |
| K7015CSP             | 2.ABDIL          | Conjugate stabilizing buffer               | 12 ml           |
| K7015RP              | REABUF           | Reaction buffer                            | 2 x 100 ml      |
| K7015DR              | DER              | Derivatization reagent                     | 2 x 13.3 mg     |
| K7015LM              | DMSO             | Dimethylsulfoxide (DMSO)                   | 2 ml            |
| K7015AP              | ASYBUF           | Assay buffer concentrate (10-fold)         | 2 x 3 ml        |
| K7015TMB             | SUB              | TMB substrate                              | 15 ml           |
| K7015AC              | STOP             | Stop solution                              | 15 ml           |

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Precision pipettors and disposable tips to deliver 10-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm (reference wave length 620 or 690 nm)

- \* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25°C (≤18.2 MΩ cm).

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The wash buffer concentrate is stable at **2-8°C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8°C for one month**.
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20°C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use. They can be re-frozen up to 3 times.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (13.3 mg) in 0.8 ml DMSO**. Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use**. Discard any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.
- Dilute the **assay buffer concentrate (ASYBUF)** with reaction buffer (REABUF) **1:10** before use: Add **27 ml REABUF** to the **3 ml ASYBUF concentrate** in the bottle, mix well. The ELISA kit can be separated into two performances by providing 2 x 3 ml assay buffer concentrate. Diluted assay buffer can be stored at **2-8°C for 4 weeks**.

- Dissolve the content of one vial of **anti-tyrosine antibody (AB)** in **5.5 ml of diluted wash buffer**. The ELISA kit can be separated into two performances by providing two AB vials. Diluted anti-tyrosine antibody can be stored at **2-8°C for 4 weeks**.
- Dilute the **POD antibody (2.AB) 1:200** with conjugate stabilizing buffer (2.ABDIL) (e.g. **60 µl 2.AB + 12 ml 2.ABDIL, prepare only the required amount**). The undiluted POD antibody (2.AB) is stable at 2-8°C until the expiry date stated on the label. Diluted POD antibody can be stored at **2-8°C for 5 days**.
- All other test reagents are ready for use. Test reagents are stable until the expiry date (see label) when stored at 2-8°C.

## 6. PRECAUTIONS

- It is for research use only.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

## 7. SPECIMEN COLLECTION AND PREPARATION

- Serum and EDTA-Plasma is suited for this test system. The samples are stable for up to 6 hours at 2-8°C. For longer storage keep samples frozen at -20°C.
- Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.
- Samples with visible amounts of precipitates should be centrifuged.
- The EDTA plasma and serum samples are diluted for derivatization (see sample preparation procedure).
- For sample preparation, a derivatization reagent (DER) for derivatization of tyrosine is added (details are given in the sample preparation procedure).

## 8. ASSAY PROCEDURE

### Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed according to the enclosed manual.

### *Sample preparation procedure*

Dilute standards (STD), controls (CTRL) and samples (SAMPLE) by factor **1:41** as follows:

#### **25 µl STD/ CTRL/ SAMPLE + 1000 µl REABUF (reaction buffer)**

Please note: Samples from patients with tyrosine supplementation (e.g. in depletion studies) may have to be more diluted.

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials) as follows:

|   |
|---|
| 1. Bring all reagents and samples to <b>room temperature (18-26°C)</b> .  |
| 2. Add <b>100 µl of diluted standards (STD)/ controls (CTRL)/ samples (SAMPLE)</b> in the corresponding vial.   |
| 3. Add <b>25 µl</b> of freshly prepared <b>derivatization reagent (DER)</b> into each vial (STD, CTRL, SAMPLE), <b>mix well</b> and incubate for <b>45 min at room temperature (18-26°C)</b> on a shaker (180-240 rpm). |
| 4. Afterwards add <b>1000 µl of assay buffer (ASYBUF)</b> into each vial, mix well and incubate for <b>45 min at room temperature (18-26°C)</b> on a shaker (180-240 rpm).  |

**2 x 25 µl of each treated sample (STD, CTRL, SAMPLE)** are used in the ELISA as duplicates.

### *Test procedure*

|     |  |
|-----|--|
| 5.  | Mark the positions of standards (STD)/ controls (CTRL)/ samples (SAMPLE) in duplicate on a protocol sheet.   |
| 6.  | Take as many microtiter plate strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.  |
| 7.  | <b>Wash the microtiter plate strips 5 times by dispensing 250 µl of diluted wash buffer (WASHBUF) into each well.</b> After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.                |
| 8.  | For the analysis in duplicate, take <b>2 x 25 µl</b> of derivatized <b>standards/ controls/ samples (STD/ CTRL)/ SAMPLE</b> out of the vials and add into the respective wells of the microtiter plate.  |
| 9.  | Add <b>100 µl</b> of diluted <b>anti-tyrosine antibody (AB)</b> into each well. Cover the plate tightly.   |
| 10. | Incubate overnight ( <b>15-20 hours</b> ) at <b>2-8°C</b> .  |
| 11. | Aspirate the contents of each well. Wash each well <b>5 times</b> by dispensing <b>250 µl of diluted wash buffer</b> into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution. |
| 12. | Add <b>100 µl</b> of diluted <b>POD antibody (2. AB)</b> into each well.   |
| 13. | Cover the plate and incubate for <b>1 hour at room temperature</b> (18-26°C) on a horizontal shaker (180-240 rpm).   |
| 14. | Aspirate the contents of each well. Wash each well <b>5 times</b> by dispensing <b>250 µl of diluted wash buffer</b> into each well. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution. |

|     |   |
|-----|---|
| 15. | Add <b>100 µl</b> of <b>TMB substrate (SUB)</b> into each well.   |
| 16. | Incubate for <b>8-12 min at room temperature</b> (18-26°C) in the dark*.  |
| 17. | Add <b>100 µl of stop solution (STOP)</b> into each well, mix thoroughly.   |
| 18. | Determine <b>absorption</b> immediately with an ELISA reader at <b>450 nm</b> . If the highest extinction of the standards ( <b>STD</b> ) is above the range of the photometer, absorption must be measured immediately at <b>405 nm</b> and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used |

\* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

## 9. EVALUATION OF RESULTS

If the test is performed in strict compliance with the manufacturer's instructions (i.e. with the exact volumes for standards, controls, samples, and with correct sample treatment), standards, controls, and blood samples are equally diluted. Therefore, **no dilution factor is required for the calculation of results.**

Exception: For samples with higher dilutions the deviating factor must be considered.

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

### 1. 4-parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.01).


Plausibility of the measured pairs of values should be examined before automatically evaluating the results. If this option is not available in the used program, the pairs of values should be controlled manually.

#### Controls

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside the acceptable limits.

The concentration of controls and patient samples can be determined directly from the calibration curve in  $\mu\text{mol/l}$ . In the following, an example of a calibration curve is given, do not use it for the calculation of your results.

#### Example of calibration curve

### Expected results

Based on internal studies with plasma samples of evidently healthy persons (n=146) a mean value of 58  $\mu\text{mol/l}$  was calculated. The standard deviation was 14.4  $\mu\text{mol/l}$ .

**Plasma mean value  $\pm$  2 x standard deviation:** **58  $\pm$  28.8  $\mu\text{mol/l}$**   
**Normal range:** **29 – 87  $\mu\text{mol/l}$**

We recommend each laboratory to develop its own normal range. The values mentioned above are indicative only and can deviate from other published data.

## 10. PERFORMANCE CHARACTERISTICS

### Cross reactivity

L-Phenylalanine < 2 %  
 L-Tryptophan < 0.5 %

### Precision and reproducibility

| Intra-assay (n=19) |                                |                                   |
|--------------------|--------------------------------|-----------------------------------|
| Sample             | Tyrosine [ $\mu\text{mol/l}$ ] | coefficient of variation (CV) [%] |
| 1                  | 31.4                           | 7,6                               |
| 2                  | 125.3                          | 7.4                               |

| Inter-assay (n=8) |                                |                                   |
|-------------------|--------------------------------|-----------------------------------|
| Sample            | Tyrosine [ $\mu\text{mol/l}$ ] | coefficient of variation (CV) [%] |
| 1                 | 69.7                           | 5.8                               |
| 2                 | 104.6                          | 4.6                               |

### Sensitivity

The sensitivity was set as  $B_0 + 2SD$ . The zero-standard was measured 82 times.

| Sample        | Tyrosine mean value [OD] | 2 x standard deviation (SD) | Detection limit [ $\mu\text{mol/l}$ ] |
|---------------|--------------------------|-----------------------------|---------------------------------------|
| Zero-standard | 21.24                    | 0.15                        | 3.6                                   |

### Recovery

Two plasma samples were spiked with different tyrosine concentrations and measured in this assay. The analytical recovery rate was determined by the expected and measured tyrosine levels. The mean recovery rate for all concentrations was 93.2 % (n=5).

| Sample [ $\mu\text{mol/l}$ ] | Spike [ $\mu\text{mol/l}$ ] | Tyrosine expected [ $\mu\text{mol/l}$ ] | Tyrosine measured [ $\mu\text{mol/l}$ ] | Recovery [%] |
|------------------------------|-----------------------------|---|---|--------------|
| 76.8                         | 25                          | 101.8                                   | 94.5                                    | 92.8         |
| 76.8                         | 50                          | 126.8                                   | 123.3                                   | 97.2         |
| 108.8                        | 25                          | 133.8                                   | 116.1                                   | 86.8         |
| 108.8                        | 50                          | 158.8                                   | 152.4                                   | 96.0         |

### Linearity

The linearity of the ELISA was determined by the dilution of two spiked plasma samples. The mean linearity was 97.1 % (n=5).

| Sample [ $\mu\text{mol/l}$ ] | Dilution | Expected [ $\mu\text{mol/l}$ ] | Measured [ $\mu\text{mol/l}$ ] | Recovery [%] |
|------------------------------|----------|--------------------------------|--------------------------------|--------------|
| 76.8                         | 1:2      | 38.4                           | 34.0                           | 88.4         |
| 76.8                         | 1:4      | 19.2                           | 20.0                           | 103.9        |
| 108.8                        | 1:2      | 54.4                           | 50.6                           | 93.0         |
| 108.8                        | 1:4      | 27.2                           | 28.0                           | 102.9        |

## 11. LIMITATIONS

Samples with tyrosine concentrations above 300 µmol/l should be further diluted and re-assayed.

Hemolytic and lipemic samples may give erroneous results. Do not measure hemolytic and lipemic samples.

## 12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the test package are for research use only.
- Reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature, and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

### Used symbols:

|   |                        |   |                                   |
|---|------------------------|---|-----------------------------------|
|  | Temperature limitation |  | Catalogue Number                  |
|  | For research use only  |  | Contains sufficient for <n> tests |
|  | Manufacturer           |  | Use by                            |
|  | Lot number             |   |                                   |