# Zinc-Protoporphyrin/ Protoporphyrin HPLC Kit

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

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>15</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>4. MATERIAL SUPPLIED</b>	<b>16</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>6. STORAGE AND PREPARATION OF REAGENTS</b>	<b>17</b>
<b>7. PRECAUTIONS</b>	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>17</b>
<b>9. ASSAY PROCEDURE</b>	<b>18</b>
<i>Procedural notes</i>	18
<i>Test procedure</i>	18
<i>Chromatographic conditions</i>	18
<b>10. TREATMENT OF THE COLUMN</b>	<b>19</b>
<b>11. RESULTS</b>	<b>19</b>
<i>Calculation</i>	19
<i>Typical chromatogram</i>	19
<b>12. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<i>Controls</i>	20
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Accuracy – Precision</i>	20
<i>Accuracy – Trueness</i>	21
<i>Analytical sensitivity</i>	21
<i>Analytical specificity</i>	22
<b>14. DISPOSAL</b>	<b>22</b>
<b>15. TROUBLESHOOTING</b>	<b>22</b>
<b>16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>17. REFERENCES</b>	<b>24</b>

## 1. INTENDED USE

This HPLC application is intended for the quantitative determination of zinc-protoporphyrin and protoporphyrin in EDTA whole blood. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Zinc protoporphyrin (ZPP) is a metabolite formed in erythrocytes during hemoglobin synthesis. By iron deficiency during erythropoiesis, instead of incorporating a ferrous ion to form the hem precursor, zinc becomes an alternative metal ion in the protoporphyrin complex. As a result, zinc-protoporphyrin is produced instead of hemoglobin. As a consequence, in cases of iron deficiency anemia, the zinc protoporphyrin concentration in the erythrocytes is elevated.

### Indications

- Lead poisoning
- Iron deficiency
- Sickle cell anemia
- Sideroblastic anemia
- Anemia of chronic disease
- Vanadium exposure

## 3. PRINCIPLE OF THE TEST

The first step in the determination of zinc-protoporphyrin and protoporphyrin is precipitation of the higher molecular components. After their removal by centrifugation, the supernatant is injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 10 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas/heights by the external standard method.

Parallel measurements of the hemoglobin concentration are required, because the zinc-protoporphyrin concentration is related to the hemoglobin concentration.

### Summary

This HPLC technique provides an easy, fast and precise method for quantitative determination of zinc-protoporphyrin and protoporphyrin. Except the column, the kit contains all reagents necessary for sample preparation and separation in ready-to-use form.

As for many other parameters, the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical chemical routines quickly and precisely. Mostly, a one-point calibration is sufficient for calibrating the test system – unlike immunoassays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. With short test series, the one-point calibration is much more economic than 6-point calibration for immunoassays.

#### 4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KC2700LM	MOPHA	Mobile phase	1000 ml
KC2700KA	CAL	Calibrator (500 µl lyoph.; see specification data sheet for concentration)	5 vials
KC2700IS	INT STD	Internal Standard (2 ml lyoph.)	4 x 1 vials
KC2700RE	RECSOL	Reconstitution solution	15 ml
KC2700FR	PREC	Precipitation reagent	150 ml
KC2700KO	CTRL 1 CTRL 2	Control 1 and 2 (250 µl lyoph.; see specification data sheet for concentration)	2 x 3 vials

The HPLC column (KC2700RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- 2 ml reaction tubes (e. g. Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with fluorescence detector
- Reversed phase C<sub>18</sub> column
- Vortex

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 6. STORAGE AND PREPARATION OF REAGENTS

- **The lyophilised calibrator (CAL)** is stable at **-20°C** until the expiry date stated on the label. Before use, the CAL has to be reconstituted with **500 µl reconstitution solution (RECSOL)**. The concentration of zinc-protoporphyrin slightly changes from lot to lot, the exact concentration is stated on the specification data sheet. **Calibrator (reconstituted CAL) is not stable and cannot be stored.**
- **The lyophilised controls 1 and 2 (CTRL 1 and CTRL 2)** are stable at **-20°C** until the expiry date stated on the label. Before use, they have to be reconstituted with each **250 µl reconstitution solution (RECSOL)**. The concentration of zinc-protoporphyrin slightly changes from lot to lot, the exact concentration is stated on the specification data sheet. **Controls (reconstituted CTRL 1 and 2) are not stable and cannot be stored.**
- **The precipitation reagent (PREC)** has to reach **room temperature (15–30°C)** before use.
- **The lyophilised internal standard (INT STD)** is stable at **-20°C** until the expiry date stated on the label. Before use, the INT STD has to be reconstituted with **2 ml precipitation reagent (PREC)**. **Internal standard (reconstituted INT STD) is not stable and cannot be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

## 7. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents should not be used beyond the expiration date stated on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

EDTA whole blood is suited for this test. Before analysis, the erythrocytes are lysed by freezing and thawing in order to release zinc-protoporphyrin and protoporphyrin.

Zinc-protoporphyrin as well as protoporphyrin are light sensitive. Samples should be transported lightproof.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- The assay should always be performed according to the enclosed manual.
- Use mobile phase (MOPHA) for needle wash in order to avoid sample carry-over.

### *Test procedure*

Pipet each <b>200 µl sample, calibrator or control 1 or 2</b> into a 2 ml reaction tube.
Add each <b>10 µl internal standard</b> .
Add each <b>1.2 ml precipitation reagent (PREC)</b> and mix well (vortex for at least 15 s).
Centrifuge for <b>10 min at 10000 g</b> and take the supernatant.
Inject <b>100 µl</b> supernatant into the HPLC system.

### *Chromatographic conditions*

Column material: Bischoff Prontoasil 120-5-C18 ace EPS; 5 µm

Column dimension: 125 mm × 4 mm

Flow rate: 1.0 ml/min

Please refer to the quality certificate of the column

Fluorescence detection:: Excitation: 417 nm

Emission: 635 nm

Temperature: 30 °C

Injection volume: 100 µl

Running time: 10 min

It is recommended to use a guard column to extend column life.

## 10. TREATMENT OF THE COLUMN

After analysis, the column should be flushed with 30 ml ultra pure water (1 ml/min) and stored in 50% methanol in water (~ 30 ml, flow 0.5 ml/min). Before use, the system should be equilibrated with ~ 30 ml mobile phase (MOPHA).

## 11. RESULTS


### Calculation

$$\frac{\text{Peak height sample} \times \text{calibrator concentration}^*}{\text{Peak height internal standard of the sample}} \times F = \text{sample concentration}$$
$$F = \frac{\text{Peak height internal standard of the calibrator}}{\text{Peak height calibrator}}$$

\* see specification data sheet

**Tip:** Alternatively, the peak area instead of the peak height can be used for quantification.

### Typical chromatogram

## 12. QUALITY CONTROL

### *Reference range*

Zinc-Protoporphyrin: < 40 µmol/mol haem<sup>[1]</sup>

We recommend each laboratory to establish its own reference range. The values mentioned above are only for orientation and may deviate from other published data.

### *Controls*

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## 13. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n=10**

The repeatability was assessed according to CLSI guideline EP5-A2 with 2 EDTA whole blood samples under constant parameters (same operator, measurement system, day and kit lot).

	<b>Sample</b>	<b>Mean value [nmol/l]</b>	<b>CV [%]</b>
1	Zinc-Protoporphyrin	348.6	3.0
	Protoporphyrin	197.8	5.5
2	Zinc-Protoporphyrin	148.8	8.3
	Protoporphyrin	59.0	5.7

#### **Reproducibility (Inter-Assay); n=10**

The reproducibility was assessed according to CLSI guideline EP5-A2 with 2 EDTA whole blood samples under varying parameters (different operators, measurement systems, days and kit lots).

	Sample	Mean value [nmol/l]	CV [%]
1	Zinc-Protoporphyrin	341.1	8.7
	Protoporphyrin	270.6	9.7
2	Zinc-Protoporphyrin	133.3	9.1
	Protoporphyrin	76.9	8.8

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, 9 EDTA whole blood samples with known concentrations were measured and the recovery was calculated.

The recovery for Zinc-Protoporphyrin was found between 94.2 and 107.9%.

The recovery for Protoporphyrin was found between 100.5 and 106.2%.

### Lower detection limit

The limit of detection (LoD) is defined as 3 times the background noise. It is calculated by the formula below:

$$\frac{(3 \times \text{peak height background noise}) \times \text{concentration calibrator [nmol/l]}}{\text{Peak height calibrator}} = \text{LoD [nmol/l]}$$

LoD Zinc-Protoporphyrin: 13.81 nmol/l

LoD Protoporphyrin: 9.43 nmol/l

The limit of quantitation (LoQ) is defined as 10 times the background noise. It is calculated by the formula below:

$$\frac{(10 \times \text{peak height background noise}) \times \text{concentration calibrator [nmol/l]}}{\text{Peak height calibrator}} = \text{LoQ [nmol/l]}$$

LoQ Zinc-Protoporphyrin: 46.02 nmol/l

LoQ Protoporphyrin: 31.42 nmol/l

### Analytical sensitivity

#### Upper detection limit & linearity

The upper limit of detection states up to which concentration a method results in a linear signal.

Zinc-Protoporphyrin: 1000 nmol/l

Protoporphyrin: 1000 nmol/l

It should be noted that the detection limits depend on the instrument as well as on the application.

### *Analytical specificity*

There were no known interferences to other blood components found.

## 14. DISPOSAL

The mobile phase (MOPHA) and precipitation reagent (PREC) must be disposed as non-halogenated solvents.

Please refer to the appropriate national guidelines.

## 15. TROUBLESHOOTING

Problem	Possible cause	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column

Problem	Possible cause	Solution
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

## 16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date stated on kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of

the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.

- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

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### Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use