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Manual*For professional use only*

Histamine elimination ratio (HERO) ELISA

*For the in vitro determination of histamine elimination
in serum*

Gültig ab / Valid from 2025-01-20

K 8215


192

K 8215.10

10 x 192

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: + 49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Safety information

The assay has to be performed exclusively according to the instructions for use enclosed with the kit.

Important safety information for this product can be found in the chapter **WARNINGS AND PRECAUTIONS**.

Table of Contents

1. INTENDED USE	21
2. INTRODUCTION	21
3. MATERIAL SUPPLIED	21
4. MATERIAL REQUIRED BUT NOT SUPPLIED	22
5. STORAGE AND PREPARATION OF REAGENTS	23
6. STORAGE AND PREPARATION OF SAMPLES	23
7. ASSAY PROCEDURE	24
<i>Principle of the test</i>	24
<i>Incubation and derivatisation procedure</i>	24
<i>Test procedure</i>	26
8. RESULTS	28
9. LIMITATIONS	29
<i>Biotin interference</i>	29
10. QUALITY CONTROL	29
<i>Reference range</i>	29
11. PERFORMANCE CHARACTERISTICS	29
<i>Precision and reproducibility</i>	29
<i>Accuracy – Trueness</i>	30
<i>Linearity</i>	30
<i>Analytical sensitivity</i>	31
<i>Specificity</i>	32
12. WARNINGS AND PRECAUTIONS	32
13. TECHNICAL HINTS	32
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	33
15. REFERENCES	33
16. SYMBOLS	34

1. INTENDED USE

The histamine assay K 8215 is an enzyme-linked immunosorbent assay (ELISA) intended for professional laboratory users for the quantitative measurement of serum histamine elimination from patients of all ages and gender.

The assay is an *in vitro* diagnostic medical device, that can be used manually or by an automated platform.

2. INTRODUCTION

Histamine is a biogenic amine derived from the decarboxylation of histidine. It is synthesised in mast cells, basophils, platelets, histaminergic neurons and enterochromaffine cells, where it is stored in vesicles. After stimulation and release, histamine acts by binding to its 4 receptors (H1R, H2R, H3R and H4R) on target cells in various tissues.

Its effects include smooth muscle cell contraction, vasodilation, increased vascular permeability, increased mucous secretion, tachycardia, blood pressure changes and arrhythmias. As a vasoactive mediator, it plays a dominant role in allergic diseases such as allergic rhinitis (hay fever), allergic bronchial asthma, and urticaria.

Histamine is mainly degraded by two pathways: By the extracellular diaminoxidase (DAO) and by the predominantly intracellular histamine N-methyltransferase (HNMT).

Measurement of serum histamine depletion is of clinical interest in patients showing one of the above-mentioned allergic or pseudoallergic reactions to histamine-rich foods.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity for cat.no.	
			K 8215	K 8215.10
K 8215	PLATE	Microtiter plate, pre-coated	2 x 12 x 8 wells	20 x 12 x 8 wells
K 0015	COPLATE	Plate for derivatisation	2 x 12 x 8 wells	20 x 12 x 8 wells
K 8215	INCPLATE	Plate for incubation	12 x 8 wells	10 x 12 x 8 wells
K 8215	STD	Standards, ready-to-use (0; 1; 3; 10; 30; 120 ng/ml)	6 x 2 ml	5 x 6 x 2 ml

K 8215	CTRL 1	Control, ready-to-use (see specification for range)	1 x 2 ml	5 x 2 ml
K 8215	CTRL 2	Control, ready-to-use (see specification for range)	1 x 2 ml	5 x 2 ml
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10 x	2 x 100 ml	15 x 100 ml
K 8215	SAMPLEBUF	Sample buffer	1 x 22 ml	10 x 22 ml
K 8215	AB	Histamine antibody, peroxidase-labelled, ready- to-use	2 x 6 ml	20 x 6 ml
K 8215	REABUF	Reaction buffer, ready-to- use	1 x 70 ml	8 x 70 ml
K 8215	DER	Derivatisation reagent, lyophilised	2 vials	20 vials
K 0008.10	DMSO	Dimethylsulfoxide (DMSO)	2 x 10 ml	20 x 10 ml
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	2 x 15 ml	20 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	2 x 15 ml	20 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Calibrated precision pipets and 10-1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker (1100 rpm and 550 rpm with orbit 2 mm, for example model ShakeID 2, available at Immundiagnostik AG)
- Multi-channel pipets or repeater pipets
- Vortex
- Incubator 37 °C
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** can be used until the expiry date stated on the label when stored at **2-8 °C**. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.
- **DMSO** crystallises at 2-8 °C. Before use, bring to room temperature to dissolve the crystals.
- The **lyophilised derivatisation reagent (DER)** can be used until the expiry date stated on the label when stored at **2-8 °C**. Bring to room temperature before opening and dissolve the content of the vial in **DMSO** as stated on the label. Allow to dissolve for **15 min** and mix thoroughly with a vortex-mixer. Ensure that it is completely dissolved. **The derivatisation reagent** (reconstituted DER) can be stored at **2-8 °C for 2 months**. Bring to room temperature before reuse. Please note: DMSO attacks plastics but not polypropylene products and laboratory glass.
- All other test reagents are ready-to-use. Test reagents can be used until the expiry date stated on the label when stored at **2-8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

Serum samples are stable for 3 days at room temperature or for 14 days at 2-8 °C. For longer storage keep samples frozen at -20 °C.

Serum samples are analysed **undiluted**.

For sample preparation, samples are treated with a histamine-containing sample buffer, stored, and a derivatisation reagent is added (see assay procedure).

7. ASSAY PROCEDURE

Principle of the test

This ELISA is designed for the quantitative determination of histamine elimination in serum. The test is based on the method of competitive enzyme linked immunoassays.

In a first step, the serum samples are mixed with a sample buffer containing histamine and aliquots are treated at 4 °C or 37 °C for 24 hours.

The histamine concentrations of both aliquots are then measured in the ELISA.

The preparation of standards, controls, and aliquots includes the addition of a derivatisation reagent for histamine derivatisation. The prepared samples and a peroxidase-conjugated polyclonal histamine antibody are then incubated in wells of a microtiter plate coated with histamine derivative (tracer). During the incubation period, the target histamine in the sample competes with the tracer, immobilised on the wall of the microtiter wells, for the binding of the polyclonal antibodies.

After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The colour changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow colour is inverse proportional to the histamine concentration in the sample; this means, high histamine concentration in the sample reduces the concentration of tracer-bound antibody and lowers the photometric signal. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standards. Histamine, present in the samples, is determined directly from this curve.

Incubation and derivatisation procedure

The incubation plate (INCPLATE) is used to treat the aliquots at 37 °C. Alternatively, reaction vials (e.g. 1.5 ml polypropylene vials) can be used.

The derivatisation of standards, controls, 37 °C and 4 °C aliquots is carried out in the derivatisation plate (COPLATE). Alternatively, micro reaction vials made of polypropylene can be used.

We recommend preparing one derivatisation per standard, control, and aliquot and transferring it in duplicate determinations into the wells of the microtiter plates (PLATE).

Day 1:

1.	Add 100 µl serum into the respective wells of the INCPLATE .
2.	Add 200 µl sample buffer (SAMPLEBUF) to each sample in the INCPLATE . Close the INCPLATE with the cap mat and mix thoroughly for 1 minute on a horizontal shaker.
3.	Add 25 µl standard (STD)/ control (CTRL) in the wells of the COPLATE .
4.	Pipette 25 µl prepared sample from the INCPLATE into the COPLATE .
5.	Cover the COPLATE tightly with foil and store at 2-8 °C for 24 h (4 °C aliquot).
6.	Close the INCPLATE tightly with the cap mat and incubate at 37 °C for 24 h (37 °C aliquot).

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Example of plate layout COPLATE (2-8 °C), optimised for multichannel pipette

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italic = 4 °C aliquots

Day 2:

Bring **all reagents** of the kit **to room temperature** (15-30 °C) and mix well.

7.	Pipette 25 µl of the 37 °C aliquots from the INCPLATE into the respective wells of the COPLATE .
8.	Add 250 µl reaction buffer (REABUF) into each well (STD, CTRL, aliquots) of the COPLATE.
9.	Add 50 µl derivatisation reagent into each well (STD, CTRL, aliquots) of the COPLATE. Mix thoroughly on a horizontal shaker . We recommend shaking at 1100 rpm with an orbit of 2 mm to achieve complete mixing.
10.	Incubate the COPLATE for 30 min on a horizontal shaker at room temperature (15-30 °C).

Example of plate layout COPLATE (derivatisation), optimised for multichannel pipette:

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italic = 4 °C aliquots

Test procedure

Take as many microtiter strips (PLATE) as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2-8 °C. Strips are stable until the expiry date stated on the label.

11.	For the analysis in duplicate take 2 x 50 µl of the derivatised standards/ controls/ aliquots out of the COPLATE and add into the respective wells of the microtiter plates (PLATE).
12.	Add 50 µl histamine antibody (AB) into each well of the microtiter plate.
13.	Cover the strips tightly with foil and incubate for 1 hour at 18-25 °C on a horizontal shaker . We recommend shaking at 500 rpm with an orbit of 2 mm.

14.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
15.	Add 100 µl substrate (SUB) into each well.
16.	Incubate for 12-18 min* at 18-25 °C in the dark .
17.	Add 100 µl stop solution (STOP) into each well and mix well.
18.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference.

* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

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Example of plate layout PLATE 2, optimised for multichannel pipette:

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD1	STD1	a21	a21	<i>a21</i>	<i>a21</i>	a29	a29	<i>a29</i>	<i>a29</i>	a37	a37
B	STD2	STD2	a22	a22	<i>a22</i>	<i>a22</i>	a30	a30	<i>a30</i>	<i>a30</i>	a38	a38
C	STD3	STD3	a23	a23	<i>a23</i>	<i>a23</i>	a31	a31	<i>a31</i>	<i>a31</i>	a39	a39
D	STD4	STD4	a24	a24	<i>a24</i>	<i>a24</i>	a32	a32	<i>a32</i>	<i>a32</i>	a40	a40
E	STD5	STD5	a25	a25	<i>a25</i>	<i>a25</i>	a33	a33	<i>a33</i>	<i>a33</i>	<i>a37</i>	<i>a37</i>
F	STD6	STD6	a26	a26	<i>a26</i>	<i>a26</i>	a34	a34	<i>a34</i>	<i>a34</i>	<i>a38</i>	<i>a38</i>
G	CTRL1	CTRL1	a27	a27	<i>a27</i>	<i>a27</i>	a35	a35	<i>a35</i>	<i>a35</i>	<i>a39</i>	<i>a39</i>
H	CTRL2	CTRL2	a28	a28	<i>a28</i>	<i>a28</i>	a36	a36	<i>a36</i>	<i>a36</i>	<i>a40</i>	<i>a40</i>

bold = 37 °C aliquots

italic = 4 °C aliquots

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. Attention should be paid to the following points:

Derivatisation of the samples is possible in 30 minutes on a horizontal shaker or, after thorough mixing, in 90 minutes without shaking in an automated processor. In automated processing, the volumes of samples and diluents may be scaled while maintaining the respective dilutions. Device-specific minimum and maximum volumes must be taken into account. The homogeneity of the resulting dilution must be ensured.

For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

1. 4 parameter algorithm

It is recommended using a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the program used, the duplicate values should be evaluated manually.

Calculation of the histamine elimination ratio

For the calculation of the histamine elimination ratio (HERO), the difference of the histamine concentrations of 4°C aliquot and 37°C aliquot is divided by the concentration of the 4°C aliquot.

$$\text{HERO [\%]} = \frac{c(\text{Histamin})_{4^{\circ}\text{C aliquot}} - c(\text{Histamin})_{37^{\circ}\text{C aliquot}}}{c(\text{Histamin})_{4^{\circ}\text{C aliquot}}} \times 100$$

9. LIMITATIONS

Biotin interference

Samples containing a biotin concentration of ≤ 343 ng/ml show a change of the results of < 25 %. Higher concentrations of biotin can lead to false results. Patients taking > 5 mg biotin per day should wait at least 24 hours after taking biotin to have their samples collected. Results of patients taking biotin supplements or receiving a high-dose biotin therapy should generally be interpreted along with the total clinical picture.

10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference range

Based on internal studies with samples of apparently healthy persons ($n = 42$) a mean histamine elimination ratio of 45 % (median) was calculated. Weak, moderate and good elimination was defined as follows:

Weak Elimination:	HERO < 25 %
Moderate Elimination:	HERO 25-40 %
Good Elimination	HERO > 40 %

We recommend that each laboratory establishes its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Repeatability (Intra-Assay); n = 6

The repeatability was assessed with 2 samples under **constant** parameters (same operator, instrument, day and kit lot) in single determinations.

Sample	Mean value 37°C aliquot [%]	Mean value 4°C aliquot [%]	Mean value HERO [%]	CV [%]
1	1.8	23.0	92.2	3.4
2	8.2	20.7	60.2	6.3

Reproducibility (Inter-Assay); n = 6

The reproducibility was assessed with 6 samples from 2 persons under **varying** parameters (different operators, instruments, days and kit lots) in duplicate determinations.

Samples from	Mean value 37°C aliquot [%]	Mean value 4°C aliquot [%]	Mean value HERO [%]	CV [%]
A	2.8	18.1	84.4	8.0
B	7.5	18.1	57.8	11.1

Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, histamine spikes with known concentrations were added to 3 different samples.

sample [ng/ml]	spike [ng/ml]	expected [ng/ml]	obtained [ng/ml]	recovery [%]
0.35	3.10	3.45	3.61	104.7
	6.43	6.78	6.35	93.7
	9.77	10.12	8.65	85.5
0.49	3.01	3.50	3.48	99.6
	6.34	6.83	6.27	91.8
	9.68	10.16	8.67	85.4
0.59	2.94	3.53	3.36	95.1
	6.28	6.86	5.94	86.6
	9.61	10.20	8.02	78.6

Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed with a serial dilution of 4 spiked serum samples with a low level serum sample containing 0.6 ng/l histamine.

For histamine in serum, the method has been demonstrated to be linear from 3.0 to 30.5 ng/ml with a recovery rate of 83.6 to 117.9 %.

sample [ng/ml]	dilution	expected [ng/ml]	obtained [ng/ml]	recovery [%]
19.8	1:1.5	13.4	13.7	102.6
	1:2	10.2	10.5	103.5
	1:3	7.0	7.1	101.4
	1:4	5.4	5.6	104.2
	1:5	4.4	5.2	117.9
	1:6	3.8	4.0	106.2
	1:8	3.0	2.8	93.8
23.5	1:1.5	15.8	14.9	94.3
	1:2	12.0	10.4	86.1
	1:3	8.2	7.5	91.6
	1:4	6.3	6.1	97.5
	1:5	5.2	4.4	84.3
	1:6	4.4	4.4	99.9
	1:8	3.4	3.3	95.8
30.5	1:1.5	20.5	20.4	99.4
	1:2	15.5	13.0	83.6
	1:3	10.6	9.8	93.2
	1:4	8.1	8.5	104.9
	1:5	6.6	6.9	105.4
	1:6	5.6	6.3	113.3
	1:8	4.3	5.1	118.6
22.8	1:1.5	15.4	16.1	104.4
	1:2	11.7	10.1	86.2
	1:3	8.0	8.9	111.7
	1:4	6.1	5.3	86.2
	1:5	5.0	4.4	86.5
	1:6	4.3	4.0	93.2
	1:8	3.4	3.4	99.6

Analytical sensitivity

The matrix of the sample buffer was measured 80 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be 0.9 ng/ml.

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to histamine. The specificity is calculated in percent, in relation to the histamine-binding activity:

3-methylhistamine	< 0.1 %
tyramine	< 0.001 %
L-phenylalanine	< 0.0002 %
L-histidine	< 0.0002 %
L-tyrosine	< 0.0002 %
tryptamine	< 0.0002 %
5-hydroxyindoleacetic acid	< 0.0002 %
serotonin (5-hydroxytryptamine)	< 0.0002 %

12. WARNINGS AND PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Kit reagents contain ProClin or thimerosal as bactericides. ProClin and thimerosal are harmful to health and the environment. Substrates for enzymatic colour reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

Warning: Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.

- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during ELISA incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. Küfner MA, Schwelberger HG, Ulrich P, Hahn EG, Raithel M. Total histamine degradation capacity (THDC) as an important biological marker of histamine metabolism in human colonic mucosa. *Inflamm Res*. 2002 Apr;51 Suppl 1:S87-8. doi: 10.1007/pl00022461. PMID: 12013425
2. Küfner MA, Ulrich P, Raithel M, Schwelberger HG. Determination of histamine degradation capacity in extremely small human colon samples. *Inflamm Res*. 2001 Apr; 50 Suppl 2:S96-7. doi: 10.1007/PL00022422. PMID: 11411621.
3. Raithel M, Küfner MA, Ulrich P, Hahn EG. The involvement of the histamine degradation pathway by diamine oxidase in manifest gastrointestinal allergies. *Inflamm Res*. 1999 Apr; 48 Suppl 1: S75-6. doi: 10.1007/s000110050414. PMID: 10350171.

4. Kellner KH. Automatic Immunoassay For Biogenic Amines. European Patent 2 612 147 B1. March 30, 2016

16. SYMBOLS

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin