# IDKmonitor<sup>®</sup> Vedolizumab free ADA ELISA

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## 15. LITERATUR

1. Afff, W. et al., 2010. Clinical utility of measuring Vedolizumab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *The American journal of gastroenterology*, **105**(5), pp.1133–9.
2. Kopylov, U. et al., 2012. Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-Vedolizumab antibodies. *Inflammatory bowel diseases*, **18**(9), pp.1628–33.
3. Tak, P.P., 2012. A personalized medicine approach to biologic treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology (Oxford, England)*, **51**(4), pp.600–9.
4. Ordás, I. et al., 2012. Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clinical pharmacology and therapeutics*, **91**(4), pp.635–46.
5. Bender, N.K. et al., 2007. Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatology international*, **27**(3), pp.269–74.

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
# **IDKmonitor<sup>®</sup> Vedolizumab free ADA ELISA**

*For the in vitro determination of  
free human antibodies against Vedolizumab  
(e. g. ENTYVIO<sup>®</sup>) in EDTA plasma and serum*

Valid from 2025-09-25

**REF** K 9648

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>16</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>6. SAMPLE PREPARATION</b>	<b>17</b>
<i>Sample preparation</i>	17
<b>7. ASSAY PROCEDURE</b>	<b>17</b>
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
<b>8. RESULTS</b>	<b>19</b>
<b>9. LIMITATIONS</b>	<b>20</b>
<b>10. QUALITY CONTROL</b>	<b>20</b>
<i>Reference range</i>	20
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Accuracy – Precision</i>	20
<i>Linearity</i>	21
<i>Analytical sensitivity</i>	21
<i>Analytical specificity</i>	22
<b>12. PRECAUTIONS</b>	<b>22</b>
<b>13. TECHNICAL HINTS</b>	<b>22</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>23</b>
<b>15. REFERENCES</b>	<b>23</b>

## 1. INTENDED USE

This Immundiagnostik assay is intended for the determination of free human antibodies against vedolizumab (e.g. ENTYVIO®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Patients with moderately to severely active Colitis ulcerosa and Crohn's disease, who have had an inadequate response to conventional therapy or to anti-tumor necrosis factor alpha (TNF $\alpha$ ) agents, can also be treated with vedolizumab. The humanised monoclonal therapeutic antibody vedolizumab binds  $\alpha$ 4 $\beta$ 7 integrin on activated lymphocytes and stops them from migrating into the intestinal mucosa. Thus, vedolizumab suppresses the inflammatory response via a different mechanism than anti-TNF $\alpha$  agents and specifically targets inflammation in the gastrointestinal tract.

The clinical efficacy of vedolizumab therapy correlates with the trough level, that is the drug level just before the next application of vedolizumab. Several factors influence the trough level, among them dosage and frequency of vedolizumab infusion, disease activity, individual pharmacokinetics and immune reaction (formation of anti-drug antibodies, ADA).

The IDKmonitor® vedolizumab free ADA ELISA for the detection of antibodies against vedolizumab (e.g. ENTYVIO®) measures free antibodies against vedolizumab. A co-determination of rheuma factors or irregular antibodies can be excluded. In combination with the drug level determination the IDKmonitor® vedolizumab free ADA ELISA is an opportunity for the treating physician to monitor and optimise the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9648	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 9648	CONJ	Conjugate concentrate (peroxidase-labelled streptavidin)	1 x 200 $\mu$ l
K 9648	CTRL POS	Positive control, lyophilised	4 x 1 vial
K 9648	CTRL NEG	Negative control, lyophilised	4 x 1 vial
K 9648	CTRL CUT OFF	Cut-off control, lyophilised	4 x 1 vial

Cat. No.	Label	Kit components	Quantity
K 9648	AB	Detection antibody concentrate (therapeutic antibody, biotinylated)	1 x 200 µl
K 9648	ASYBUF	Assay buffer, ready-to-use	1 x 50 ml
K 0002.15	SUB	Substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Absorbent paper
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (e.g. 100 ml WASH-BUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be re-dissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is

stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

- The **lyophilised controls** (CTRL NEG, CTRL POS, CTRL CUT-OFF) are stable at **2–8 °C** until the expiry date stated on the label. **Reconstitution** details are given in the **specification data sheet**.
- **Preparation of the conjugate:** The **conjugate concentrate (CONJ)** has to be diluted **1:101** in assay buffer (e.g. 100 µl CONJ + 10 ml ASYUF). The CONJ is stable at **2–8 °C** until expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored**.
- **Preparation of the detection antibody:** The **detection antibody concentrate (AB)** has to be diluted **1:101** in assay buffer (e.g. 100 µl AB + 10 ml ASY-BUF). The AB is stable at **2–8 °C** until expiry date stated on the label. **Detection antibody** (1:101 diluted AB) **is not stable and cannot be stored**.
- All other test reagents are ready-to-use. Test reagents can be used until the expiry date (see label) when stored at **2–8 °C**.

## 6. SAMPLE PREPARATION

### *Sample preparation*

The samples are diluted **1:10 in assay buffer** (ASYBUF). For example:

- **25 µl** sample + **225 µl** assay buffer, mix well = **1:10**

For analysis, pipet **100 µl** of the **dilution** per well. For the recommended analysis in duplicate, 2x 100 µl are required.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the determination of free antibodies against vedolizumab (e.g. ENTYVIO®). In a first incubation step, the free anti-vedolizumab antibodies from the sample are bound to the vedolizumab F(ab) fragments coated on the plate. After a washing step to remove all unbound substances, the detection antibody (biotinylated therapy antibody) is added. After another washing step, the streptavidine conjugate is added to detect the anti-vedolizumab antibodies. After another washing step, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic stop solution is then added. The colour converts to yellow.

The absorbance of the colour compound is determined photometrically at 450 nm. The intensity of the colour is directly proportional to the amount of bound anti-vedolizumab antibodies (e.g. ENTYVIO®) from the sample. The results are evaluated by a cut-off control.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash the pre-coated microtiter plate <b>5 x</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
2.	Add <b>100 µl</b> of <b>CTRL NEG, CTRL POS, CTRL CUT-OFF</b> (controls, negative, positive or cut-off) and prepared samples in the wells of the microtiter plate.
3.	Seal the strips with foil and incubate <b>1 hour</b> on a <b>horizontal mixer</b> at room temperature (15–30 °C).
4.	Aspirate the content of the plate and wash each well <b>5 x</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
5.	Add <b>100 µl detection antibody</b> (1:101 diluted AB) into each well.
6.	Seal the strips with foil and incubate <b>1 hour</b> on a <b>horizontal mixer</b> at room temperature (15–30 °C).
7.	Aspirate the content of the plate and wash each well <b>5 x</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add <b>100 µl conjugate</b> (1:101 diluted CONJ) into each well.
9.	Seal the strips with foil and incubate for <b>1 hour shaking</b> on a horizontal mixer at room temperature.

10.	Aspirate the content of the plate and wash each well <b>5 x</b> with <b>250 µl wash buffer</b> . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
11.	Add <b>100 µl of TMB substrate solution</b> into each well.
12.	Incubate for <b>10–20 minutes</b> at room temperature in the <b>dark*</b> .
13.	Add <b>100 µl of stop solution</b> into each well and mix shortly.
14.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* The intensity of the colour change is temperature sensitive. We recommend to observe the procedure of the colour change and to stop the reaction upon good differentiation.

## 8. RESULTS

The analysis of the results is done using the cut-off control. Samples with a higher optical density (OD) as the OD of the cut-off control are positive. Samples with an OD lower than the OD of the cut-off control are negative.

$$\text{Cut-off} = 10 \text{ AU/ml} = \text{OD}_{\text{cut-off control}}$$

For the calculation of the sample concentrations, linear regression using a linear ordinate and abscissa is recommended.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### Sample calculation for a positive sample

average OD of patient's sample	0.735
average OD of cut-off control	0.065 = 10 AU/ml
Concentration of the patient's sample	$\frac{0,735 \times 10 \text{ AU/ml}}{0,065} = 113 \text{ AU/ml}$

## 9. LIMITATIONS

The lower limit of the measurement range can be calculated as:

$$LoB \times \text{sample dilution factor to be used}$$

Samples with concentrations lower than the measurement range cannot be clearly quantified.

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n = 20**

The repeatability was assessed with 2 serum samples under **constant** parameters (same operator, instrument, day and kit lot).

Sample	Mean value [AU/ml]	CV [%]
1	26.53	7.8
2	57.26	7.7

#### **Reproducibility (Inter-Assay); n = 21**

The reproducibility was assessed with 2 serum samples under **varying** parameters (different operators, instruments, days and kit lots).

Sample	Mean value [AU/ml]	CV [%]
1	59.18	11.8
2	30.79	10.0

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP06-A by two serial dilutions of one serum sample.

For anti-vedolizumab antibodies in serum and EDTA plasma, the method has been demonstrated to be linear from 12.49 to 148.16 AU/ml, showing a non-linear behaviour of less than  $\pm 20\%$  in this interval.

Sample	Dilution	Expected [AU/ml]	Obtained [AU/ml]	Recovery [%]
A	Undiluted	148.16	148.16	100.0
	1:1.5	98.78	107.15	108.5
	1:2.25	65.85	69.62	105.7
	1:3.38	43.90	49.24	112.2
	1:5.06	29.27	33.29	113.7
	1:7.59	19.51	22.34	114.5
	1:11.39	13.01	15.25	117.3
B	Undiluted	142.22	142.22	100.0
	1:1.5	94.81	101.90	107.5
	1:2.25	63.21	65.25	103.2
	1:3.38	42.14	47.59	113.0
	1:5.06	28.09	31.71	112.9
	1:7.59	18.73	21.65	115.6
	1:11.39	12.49	14.62	117.1

### Analytical sensitivity

The following values have been estimated based on the concentrations of the standards without considering possibly used sample dilution factors.

Limit of blank, LoB	4.60 AU/ml
Limit of detection, LoD	8.06 AU/ml
Limit of quantitation, LoQ	10.00 AU/ml

The evaluation was performed according to the CLSI guideline EP17-A2. The specified accuracy goal for the LoQ was 10 % CV.

### *Analytical specificity*

The specificity of the antibody was tested by measuring the cross-reactivity against the following therapeutic antibodies: anti-adalimumab, anti-infliximab, anti-etanercept, anti-golimumab and anti-ustekinumab. No cross-reactivity was observed.

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide or ProClin are hazardous to health and the environment. Substrates for enzymatic colour reactions may also cause skin and/or respiratory irritation. Any contact with the substances must be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDKmonitor® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Afif, W. et al., 2010. Clinical utility of measuring Vedolizumab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease. *The American journal of gastroenterology*, **105**(5), pp.1133–9.
2. Kopylov, U. et al., 2012. Clinical utility of antihuman lambda chain-based enzyme-linked immunosorbent assay (ELISA) versus double antigen ELISA for the detection of anti-Vedolizumab antibodies. *Inflammatory bowel diseases*, **18**(9), pp.1628–33.
3. Tak, P.P., 2012. A personalized medicine approach to biologic treatment of rheumatoid arthritis: a preliminary treatment algorithm. *Rheumatology (Oxford, England)*, **51**(4), pp.600–9.

4. Ordás, I. et al., 2012. Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms. *Clinical pharmacology and therapeutics*, **91**(4), pp.635–46.
5. Bender, N.K. et al., 2007. Immunogenicity, efficacy and adverse events of adalimumab in RA patients. *Rheumatology international*, **27**(3), pp.269–74.

### Used symbols:

Lot number

Catalogue number

*In Vitro* Diagnostic Medical Device

To be used with

Manufacturer

Contains sufficient for  
<n> tests

Temperature limitation

Use by

Consult product specification  
data sheet

Consult  
instructions for use

European Conformity

Do not re-use

Contains plasma derivatives  
or human blood

Irritant

Unique Device  
Identification