IDKmonitor[®] Adalimumab free ADA ELISA

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15. LITERATUR

1. W. Aff, E. V Loftus, W. a Faubion, S. V Kane, D. H. Bruining, K. a Hanson, and W. J. Sandborn, "Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease.," *The American journal of gastroenterology*, vol. **105**, no. 5, pp. 1133–9, May 2010.
2. L. L. A. Lecluse, R. J. B. Driessen, P. I. Spuls, E. M. G. J. de Jong, S. O. Stapel, M. B. A. van Doorn, J. D. Bos, and G.-J. Wolbink, "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis.," *Archives of dermatology*, vol. **146**, no. 2, pp. 127–32, Feb. 2010.
3. P. A. van Schouwenburg, T. Rispens, and G. J. Wolbink, "Immunogenicity of anti-TNF biologic therapies for rheumatoid arthritis.," *Nature reviews. Rheumatology*, vol. **9**, no. 3, pp. 164–72, Mar. 2013.
4. I. Ordás, D. R. Mould, B. G. Feagan, and W. J. Sandborn, "Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms.," *Clinical pharmacology and therapeutics*, vol. **91**, no. 4, pp. 635–46, Apr. 2012.
5. N. K. Bender, C. E. Heilig, B. Dröll, J. Wohlgemuth, F. P. Armbruster, and B. Heilig, "Immunogenicity, efficacy and adverse events of adalimumab in RA patients.," *Rheumatology international*, vol. **27**, no. 3, pp. 269–74, Jan. 2007

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
IDKmonitor[®] Adalimumab free ADA ELISA

***For the in vitro determination of free human antibodies
against adalimumab (e. g. HUMIRA[®])
in EDTA plasma and serum***

Valid from 2024-09-18

REF K 9652

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	15
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. PREPARATION AND STORAGE OF REAGENTS	16
6. STORAGE AND PREPARATION OF SAMPLES	17
7. ASSAY PROCEDURE	17
<i>Principle of the test</i>	17
<i>Test procedure</i>	18
8. RESULTS	19
9. LIMITATIONS	19
10. QUALITY CONTROL	20
<i>Reference range</i>	20
11. PERFORMANCE CHARACTERISTICS	20
<i>Accuracy – Precision</i>	20
<i>Linearity</i>	21
<i>Analytical sensitivity</i>	21
<i>Analytical specificity</i>	21
12. PRECAUTIONS	22
13. TECHNICAL HINTS	22
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23
15. REFERENCES	23

1. INTENDED USE

This enzyme-linked immunosorbent assay (ELISA) kit is intended for the determination of free anti-drug antibodies (ADA) against the therapeutic TNF α antibody adalimumab (e.g. HUMIRA®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

2. INTRODUCTION

Chronic inflammatory diseases like Crohn's disease, ulcerative colitis, rheumatoid arthritis or psoriasis are often treated with anti-TNF α antibodies which target directly the underlying inflammatory processes.

The clinical efficacy of an anti-TNF α therapy correlates with the trough level of the therapeutic antibody, the drug level just before the next application of the anti-TNF α antibody [1-3]. Several factors influence the trough level, among them dosage and frequency of anti-TNF α blocker infusion, disease activity, individual pharmacokinetics and immune reaction (formation anti-drug antibodies, ADA). It is thought that ADA functionally neutralize the therapeutic antibodies or induce their rapid elimination [3]. Consequences of ADA formation can be therapy failure and allergic reactions during anti-TNF α antibody application [1, 5].

The IDKmonitor® Adalimumab free ADA ELISA for the detection of antibodies against Adalimumab (e.g. HUMIRA®) measures free anti-Adalimumab antibodies. In combination with the drug level determination of Adalimumab, the IDKmonitor® Adalimumab free ADA ELISA is an opportunity for the treating physician to monitor and optimise the therapy early on.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9652	PLATE	Microtiter plate, precoated with (F(ab) ₂)	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	1 x 100 ml
K 9652	CONJ	Antibody concentrate, (therapy antibody, peroxidase-labelled)	1 x 200 μ l
K 9652	CTRL NEG	Negative control, lyophilised	4 x 1 vial
K 9652	CTRL POS	Positive control, lyophilised	4 x 1 vial
K 9652	CTRL CUT-OFF	Cut-off control, lyophilised	4 x 1 vial
K 9652	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 30 ml

Cat. No.	Label	Kit components	Quantity
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Absorbent paper
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles >0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.

- The **lyophilised controls (CTRL NEG, CTRL POS and CTRL CUT-OFF)** are stable at **2–8 °C** until the expiry date stated on the label. Reconstitution details are given in the specification data sheet. **Controls** (reconstituted CTRL NEG, CTRL POS and CTRL CUT-OFF) **are not stable and cannot be stored.**
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in **wash buffer** (100 µl CONJ + 10 ml wash buffer). The **CONJ** is stable at **2–8 °C** until the expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C.**

6. STORAGE AND PREPARATION OF SAMPLES

EDTA plasma and serum

EDTA plasma or serum samples must be diluted **1:10** before performing the assay, e.g. **25 µl** sample + **225 µl** sample dilution buffer (SAMPLEBUF), mix well.

For testing in duplicates, pipet **2 x 100 µl** of each prepared sample.

Sample storage

Undiluted samples can be stored for 1 month at -20 °C and for 7 days at 2–8 °C or room temperature.

Diluted samples are not stable and cannot be stored.

7. ASSAY PROCEDURE

Principle of the test

This enzyme immunoassay is a sandwich assay for the determination of free antibodies against adalimumab (e.g. HUMIRA®). In a first incubation step, the free anti-therapeutic antibodies from the sample are bound to the adalimumab F(ab)₂ fragments coated on the plate. To remove all unbound substances, a washing step is carried out. In a further incubation step, peroxidase-labelled adalimumab is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added. The colour converts to yellow. The absorbance of the colour compound is determined photometrically. The intensity of the colour is directly proportional to the amount of bound ADAs (here: anti-adalimumab antibodies) from the sample. The results are evaluated by a cut-off control.

Test procedure

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until the expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1.	Wash the microtiter strips 5 x with 250 µl wash buffer before use . After the final washing step, the inverted microtiter strips should be tapped on absorbent paper.
2.	Add 100 µl of CTRL NEG, CTRL POS, CTRL CUT-OFF (controls), and diluted samples in the wells of the microtiter plate.
3.	Seal the strips with foil and incubate for 1 hour at room temperature (15–30 °C) on a horizontal shaker*.
4.	Aspirate the content of the plate and wash each well 5 x with 250 µl wash buffer . After the final washing step, the inverted microtiter should be firmly tapped on absorbent paper.
5.	Add 100 µl conjugate into each well.
6.	Seal the strips with foil and incubate for 1 hour at room temperature (15–30 °C) on a horizontal shaker*.
7.	Aspirate the content of the plate and wash each well 5 x with 250 µl wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8.	Add 100 µl of TMB substrate solution into each well.
9.	Incubate for 10–20 minutes at room temperature in the dark**.
10.	Add 100 µl stop solution into each well and mix shortly.
11.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference.

* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

** The intensity of the colour change is temperature sensitive. We recommend to observe the procedure of the colour change and to stop the reaction upon good differentiation.

8. RESULTS

The analysis of the results is done using the cut-off control. Samples with a higher optical density (OD) as the OD of the cut-off control are positive. Samples with an OD lower than the OD of the cut-off control are negative.

Cut-off = 10 AU/ml = OD of cut-off control

For the calculation of the sample concentrations, linear regression using a linear ordinate and abscissa is recommended.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Sample calculation for a positive sample

average OD of patient's sample	0.735
average OD of cut-off control	0.065 = 10 AU/ml
Concentration of the patient's sample	$\frac{0.735 \times 10 \text{ AU/ml}}{0.065} = 113 \text{ AU/ml}$

9. LIMITATIONS

The lower limit of the measurement range is the LoB.

LoB see chapter "Performance Characteristics".

Samples with concentrations lower than the measurement range cannot be clearly quantified.

10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Provided control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Accuracy – Precision

Repeatability (Intra-Assay); n = 30

The repeatability was assessed with 3 serum-samples under **constant** parameters (same operator, instrument, day and kit lot).

Sample	Mean value [AU/ml]	CV [%]
1	87.56	2.0
2	40.81	2.4
3	45.87	2.1

Reproducibility (Inter-Assay); n = 18

The reproducibility was assessed with 3 serum samples under **varying** parameters (different operators, instruments, days and kit lots).

Sample	Mean value [AU/ml]	CV [%]
1	87.29	5.2
2	47.83	8.3
3	79.05	7.9

Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A by serial dilution of 2 serum samples.

For anti-Adalimumab antibodies in serum and plasma, the method has been demonstrated to be linear from 18.56 to 232.05 AU/ml, showing a non-linear behaviour of less than $\pm 20\%$ in this interval.

Sample	Dilution	Expected [AU/ml]	Obtained [AU/ml]	Recovery [%]
1	1:640	232.05	232.05	100.00
	1:1000	148.51	151.52	102.02
	1:1280	116.03	129.31	111.45
	1:2000	74.26	87.20	117.43
	1:4000	37.13	43.10	116.08
	1:8000	18.56	21.22	114.32
2	1:160	157.03	157.03	100.00
	1:320	78.51	84.86	108.08
	1:640	39.26	42.86	109.17
	1:1280	19.63	21.26	108.30

Analytical sensitivity

The following values have been estimated without considering possibly used sample dilution factors:

Limit of blank, LoB	1.615 AU/ml
Limit of detection, LoD	3.361 AU/ml
Limit of quantitation, LoQ	10 AU/ml

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 20 % CV.

Analytical specificity

The specificity of the antibody was tested by measuring the cross-reactivity against anti-Infliximab antibodies in patient samples. There was no cross-reactivity observed.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only
- Control samples should be analysed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide or ProClin are hazardous to health and the environment. Substrates for enzymatic colour reactions may also cause skin and/or respiratory irritation. Any contact with the substances must be avoided. Further safety information can be found in the safety data sheet, which is available from Immundiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.
Warning: Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not to assemble wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions as sealed ones.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.

- The assay should always be performed according the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- IDKmonitor® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

15. REFERENCES

1. W. Afif, E. V Loftus, W. a Faubion, S. V Kane, D. H. Bruining, K. a Hanson, and W. J. Sandborn, "Clinical utility of measuring infliximab and human anti-chimeric antibody concentrations in patients with inflammatory bowel disease.," *The American journal of gastroenterology*, vol. **105**, no. 5, pp. 1133–9, May 2010.
2. L. L. A. Lecluse, R. J. B. Driessen, P. I. Spuls, E. M. G. J. de Jong, S. O. Stapel, M. B. A. van Doorn, J. D. Bos, and G.-J. Wolbink, "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis.," *Archives of dermatology*, vol. **146**, no. 2, pp. 127–32, Feb. 2010.
3. P. A. van Schouwenburg, T. Rispens, and G. J. Wolbink, "Immunogenicity of anti-TNF biologic therapies for rheumatoid arthritis.," *Nature reviews. Rheumatology*, vol. **9**, no. 3, pp. 164–72, Mar. 2013.
4. I. Ordás, D. R. Mould, B. G. Feagan, and W. J. Sandborn, "Anti-TNF monoclonal antibodies in inflammatory bowel disease: pharmacokinetics-based dosing paradigms.," *Clinical pharmacology and therapeutics*, vol. **91**, no. 4, pp. 635–46, Apr. 2012.
5. N. K. Bender, C. E. Heilig, B. Dröll, J. Wohlgemuth, F. P. Armbruster, and B. Heilig, "Immunogenicity, efficacy and adverse events of adalimumab in RA patients.," *Rheumatology international*, vol. **27**, no. 3, pp. 269–74, Jan. 2007

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com