# IDKmonitor<sup>®</sup> Rituximab drug level ELISA

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## 15. LITERATUR

1. Falchi L, Ferrajoli A, Jacobs I, Nava-Parada P. An Evidence-based Review of Anti-CD20 Antibody-containing Regimens for the Treatment of Patients With Relapsed or Refractory Chronic Lymphocytic Leukemia, Diffuse Large B-cell Lymphoma, or Follicular Lymphoma. *Clinical lymphoma, myeloma & leukemia*. 2018 Aug;**18**(8):508–518.e14.
2. Cohen SB, Emery P, Greenwald MW, Dougados M, Furie RA, Genovese MC, et al. Rituximab for rheumatoid arthritis refractory to anti-tumor necrosis factor therapy: Results of a multicenter, randomized, double-blind, placebo-controlled, phase III trial evaluating primary efficacy and safety at twenty-four weeks. *Arthritis and rheumatism*. 2006 Sep 1;**54**(9):2793–806.
3. Singer O, McCune WJ. Update on maintenance therapy for granulomatosis with polyangiitis and microscopic polyangiitis. *Current opinion in rheumatology*. 2017 May;**29**(3):248–53.

4. Bohn J-P, Steurer M. Current and evolving treatment strategies in adult immune thrombocytopenia. *Memo*. 2018 Sep 15;**11**(3):241–6.
5. Tandan R, Hehir MK, Waheed W, Howard DB. Rituximab treatment of myasthenia gravis: A systematic review. *Muscle & nerve*. 2017 Aug;**56**(2):185–96.
6. Kuijpers TW, Bende RJ, Baars PA, Grummels A, Derks IAM, Dolman KM, et al. CD20 deficiency in humans results in impaired T cell-independent antibody responses. *The Journal of clinical investigation*. 2010 Jan;**120**(1):214–22.
7. Rudnicka D, Oszmiana A, Finch DK, Strickland I, Schofield DJ, Lowe DC, et al. Rituximab causes a polarization of B cells that augments its therapeutic function in NK-cell-mediated antibody-dependent cellular cytotoxicity. *Blood*. 2013 Jun 6;**121**(23):4694–702.

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
# **IDKmonitor<sup>®</sup> Rituximab drug level ELISA**

***For the in vitro determination of free rituximab  
concentration (e.g. MabThera<sup>®</sup>, Rituxan<sup>®</sup>)  
in EDTA plasma and serum***

Valid from 2022-02-09

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com) [www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>17</b>
<b>2. INTRODUCTION</b>	<b>17</b>
<b>3. MATERIAL SUPPLIED</b>	<b>17</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>18</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>19</b>
<b>7. ASSAY PROCEDURE</b>	<b>19</b>
<i>Principle of the test</i>	19
<i>Test procedure</i>	20
<b>8. RESULTS</b>	<b>21</b>
<b>9. LIMITATIONS</b>	<b>21</b>
<b>10. QUALITY CONTROL</b>	<b>22</b>
<i>Reference range</i>	22
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>22</b>
<i>Accuracy – Precision</i>	22
<i>Accuracy – Trueness</i>	23
<i>Analytical sensitivity</i>	23
<i>Linearity</i>	23
<b>12. PRECAUTIONS</b>	<b>25</b>
<b>13. TECHNICAL HINTS</b>	<b>25</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>26</b>
<b>15. REFERENCES</b>	<b>26</b>

## 1. INTENDED USE

The Immundiagnostik AG assay is an enzyme immunoassay intended for the quantitative determination of free chimeric anti-CD20 therapeutic antibody rituximab (e.g. MabThera®, Rituxan®) in EDTA plasma and serum. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Rituximab is a chimeric monoclonal antibody against the B cell surface antigen CD20. It is used for the therapy of malign lymphomas [1] and different autoimmune diseases (e.g. rheumatoid arthritis [2], granulomatosis with polyangiitis [3], immune thrombocytopenia [4], myasthenia gravis [5]).

CD20 (also known as human B-lymphocyte-restricted differentiation antigen or Bp35) is a surface antigen on normal and malign pre-B-lymphocytes and mature B-lymphocytes, but it is not expressed on hematopoietic stem cells, pro-B-lymphocytes, plasma cells or cells of other tissues. CD20 supports the B-cell immune response, in particular against T-cell-independent antigens [6] and maybe works as a calcium channel.

By binding to CD20, rituximab improves i. a. the effect of natural killer (NK) cells which induce cell death in the B-lymphocytes marked with rituximab [7]. Hereby the number of alive B-lymphocytes is clearly reduced. This is one of the aims of lymphoma therapy. In autoimmune disease therapy, the reduction of B-lymphocytes results in a reduction of autoantibodies, leading to a reduction of the symptoms.

Using this IDKmonitor® rituximab drug level ELISA kit, you can determine the serum or plasma level of the therapeutic antibody rituximab. The results help the treating physician to monitor and optimise the therapy early on.

## 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 9661	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 9661	CONJ	Conjugate, peroxidase-labelled, ready-to-use	1 x 15 ml
K 9661	STD	Calibrators, ready-to-use (see specification for concentrations)	1 x 6 vials
K 9661	CTRL 1	Control 1, ready-to-use (see specification for range)	1 x 1 vial

Cat. No.	Label	Kit Components	Quantity
K 9661	CTRL 2	Control 2, ready-to-use (see specification for range)	1 x 1 vial
K 9661	SAMPLEBUF	Dilution buffer, ready-to-use	1 x 100 ml
K 0002.15	SUB	Substrate (tetramethylbenzidin), ready-to-use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1 000 µl single-use tips
- Absorbent paper
- Thermo shaker (25 °C)
- Multi-channel pipets or repeater pipets
- Centrifuge
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the stock solution. The crystals must be redissolved at room

temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for one month**.

- All other test reagents are ready-to-use, except microtiter strips (PLATE) - see “Test procedure” in chapter “Assay Procedure”. Test reagents are stable until the expiry date (see label) when stored at 2–8 °C.

## 6. STORAGE AND PREPARATION OF SAMPLES

### EDTA plasma and serum

EDTA plasma or serum samples must be diluted **1:2 000** with dilution buffer (SAMPLEBUF) before performing the assay, e.g.

- **10 µl** sample + **490 µl** SAMPLEBUF = **1:50 (dilution I)**
- **10 µl** dilution I + **390 µl** SAMPLEBUF = **1:40 (dilution II)**.

This results in a **total dilution of 1:2 000** (= diluted SAMPLE).

For testing in duplicates, pipet **2 x 100 µl** per well of each prepared sample.

### Sample storage

Freshly collected EDTA plasma or serum can be stored for 7 days at 2–8 °C. At -20 °C, the samples are stable for 6 months.

Diluted EDTA plasma or serum samples can be stored for 7 days at 2–8 °C.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed to determine the quantity of free rituximab (therapeutic antibody against CD20) in EDTA plasma or serum samples. In a first incubation step, the free rituximab from the sample is bound to the specific monoclonal anti-rituximab antibody coated on the plate. To remove all unbound substances, a washing step is carried out. In a further incubation step, peroxidase-labeled antibody is added. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow. The intensity of the yellow colour is directly proportional to the concentration of free rituximab in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the values obtained from standard. The concentrations of free rituximab in the samples are determined directly from this curve.

### Test procedure

Bring all **reagents and samples to room temperature** (20–30 °C) and mix well.

Take as many microtiter strips as needed from kit. Store unused strips covered with the provided foil (FOL) together with the desiccant bag in the re-closed aluminium packaging at 2–8 °C and use them in the next 4 months.

For automated ELISA processors the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Add <b>100 µl of standards, controls or samples</b> (diluted SAMPLE) into the respective wells.
2.	Cover the plate tightly and incubate for <b>1 hour at 25 °C</b> on a thermo shaker at <b>300 rpm..</b>
3.	Discard the contents of each well. Wash each well <b>5 x</b> by dispensing <b>250 µl of wash buffer</b> into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
4.	Add <b>100 µl conjugate</b> (CONJ) into each well.
5.	Cover the plate tightly and incubate for <b>1 hour at 25 °C</b> on a thermo shaker at <b>300 rpm.</b>
6.	Discard the contents of each well. Wash each well <b>5 x</b> by dispensing <b>250 µl of wash buffer</b> into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
7.	Add <b>100 µl substrate</b> (SUB) into each well.
8.	Incubate for <b>10–15 min*</b> at room temperature (20–30 °C) <b>in the dark.</b>
9.	Add <b>100 µl stop solution</b> (STOP) into each well, mix.
10.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference.

\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### EDTA-plasma and serum samples

The obtained rituximab levels of EDTA plasma and serum samples have to be multiplied by the dilution factor of 2 000. The resulting value has the unit ng/ml. To convert the value to µg/ml, divide it by 1 000.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## 9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be further diluted and re-assayed. Please consider this higher dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

$$\text{highest concentration of the standard curve} \times \text{sample dilution factor to be used} \\ 250 \text{ ng/ml} * 2\,000 = 500\,000 \text{ ng/ml} = 500 \mu\text{g/ml}$$

The lower limit of the measurement range can be calculated as:

$$\text{LoQ} \times \text{sample dilution factor to be used:} \\ 2,54 \text{ ng/ml} * 2000 = 5080 \text{ ng/ml} = 5,080 \mu\text{g/ml}$$

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n = 32**

The repeatability was assessed with 2 serum samples under **constant** parameters (same operator, instrument, day and kit lot).

Sample	Mean value [ $\mu\text{g/ml}$ ]	CV [%]
1	38.96	2.9
2	210.15	7.2

#### **Reproducibility (Inter-Assay); n = 24**

The reproducibility was assessed with 2 serum samples under **varying** parameters (different operators, instruments, days and kit lots).

Sample	Mean value [ $\mu\text{g/ml}$ ]	CV [%]
1	38.99	9.4
2	231.42	9.6

### Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, rituximab spikes with known concentrations were added to 4 different serum samples.

Sample [µg/ml]	Spike [µg/ml]	Expected [µg/ml]	Obtained [µg/ml]	Recovery [%]
0.8	10.0	10.8	13.3	123.1
	50.0	50.8	60.6	119.3
	100.0	100.8	119.2	118.3
	300.0	300.8	255.2	84.8
1.4	7.5	8.9	8.4	94.4
	15.0	16.4	16.4	100.0
	30.0	31.4	30.0	95.5
Negative	10.0	10.0	12.4	124.0
	50.0	50.0	50.0	100.0
	100.0	100.0	98.6	98.6
	300.0	300.0	228.4	76.1
Negative	5.0	5.0	5.6	112.0
	7.5	7.5	8.6	114.7
	10.0	10.0	11.2	112.0

### Analytical sensitivity

The following values have been estimated based on the concentrations of the standard curve without considering possibly used sample dilution factors.

Limit of blank, LoB 1.59 ng/ml

Limit of detection, LoD 2.54 ng/ml

Limit of quantitation, LoQ 2.54 ng/ml

The evaluation was performed according to the CLSI guideline EP17-A2. The specified accuracy goal for the LoQ was 20 % CV.

### Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A with a serial dilution of 3 different serum samples.

For rituximab in serum and plasma, the method has been demonstrated to be linear from 3.6 to 91.2 ng/ml, showing a non-linear behaviour of less than  $\pm 20\%$  in this interval. The following values have been estimated based on the concentrations of the standard curve without considering possibly used sample dilution factors.

Sample	Dilution	Expected [µg/ml]	Obtained [µg/ml]	Recovery [%]
1	1:500	91.2	91.2	100.00
	1:1 000	43.6	43.5	99.89
	1:2 000	21.8	21.1	96.90
	1:4 000	10.9	9.4	86.34
	1:6 000	7.3	7.1	97.82
	1:8 000	5.4	5.6	102.87
	1:10 000	4.4	4.5	103.33
	1:12 000	3.6	4.2	115.73
2	1:4 000	76.6	76.6	100.00
	1:8 000	36.4	33.3	91.48
	1:12 000	24.3	22.0	90.66
	1:16 000	18.2	16.3	89.56
	1:20 000	14.6	13.5	92.72
	1:24 000	12.1	11.8	97.25
	1:30 000	9.7	9.4	96.84
	1:32 000	9.1	9.0	98.90
3	1:1 000	62.8	62.8	100.00
	1:2 000	30.3	30.6	101.16
	1:4 000	15.1	15.9	105.12
	1:6 000	10.1	10.1	100.17
	1:8 000	7.6	7.7	101.82
	1:10 000	6.1	6.2	102.48
	1:16 000	3.8	4.3	113.72

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are harmful to health and the environment. Substrates for enzymatic color reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found in the safety data sheet, which is available from Immunodiagnostik AG on request.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

**Warning:** Causes serious eye irritation. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still should be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.

- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE


- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDKmonitor*® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

## 15. REFERENCES

1. Falchi L, Ferrajoli A, Jacobs I, Nava-Parada P. An Evidence-based Review of Anti-CD20 Antibody-containing Regimens for the Treatment of Patients With Relapsed or Refractory Chronic Lymphocytic Leukemia, Diffuse Large B-cell Lymphoma, or Follicular Lymphoma. *Clinical lymphoma, myeloma & leukemia*. 2018 Aug;**18**(8):508–518.e14.
2. Cohen SB, Emery P, Greenwald MW, Dougados M, Furie RA, Genovese MC, et al. Rituximab for rheumatoid arthritis refractory to anti-tumor necrosis factor therapy: Results of a multicenter, randomized, double-blind, placebo-controlled, phase III trial evaluating primary efficacy and safety at twenty-four weeks. *Arthritis and rheumatism*. 2006 Sep 1;**54**(9):2793–806.
3. Singer O, McCune WJ. Update on maintenance therapy for granulomatosis with polyangiitis and microscopic polyangiitis. *Current opinion in rheumatology*. 2017 May;**29**(3):248–53.
4. Bohn J-P, Steurer M. Current and evolving treatment strategies in adult immune thrombocytopenia. *Memo*. 2018 Sep 15;**11**(3):241–6.

5. Tandan R, Hehir MK, Waheed W, Howard DB. Rituximab treatment of myasthenia gravis: A systematic review. *Muscle & nerve*. 2017 Aug;**56**(2):185–96.
6. Kuijpers TW, Bende RJ, Baars PA, Grummels A, Derks IAM, Dolman KM, et al. CD20 deficiency in humans results in impaired T cell-independent antibody responses. *The Journal of clinical investigation*. 2010 Jan;**120**(1):214–22.
7. Rudnicka D, Oszmiana A, Finch DK, Strickland I, Schofield DJ, Lowe DC, et al. Rituximab causes a polarization of B cells that augments its therapeutic function in NK-cell-mediated antibody-dependent cellular cytotoxicity. *Blood*. 2013 Jun 6;**121**(23):4694–702.

**Used symbols:**

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		Irritant

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)