**Kurzanleitung / Short manual**

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

*For professional use only*

# Vitamin B<sub>6</sub> HPLC Kit

*For the determination of vitamin B<sub>6</sub> in serum, plasma,  
whole blood and dried blood spots*

Valid from 2023-12-05

**REF** KC2151

**IVD** **CE**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## **Safety information**

These accessories are to be used exclusively in accordance with the enclosed instructions for use. Important safety information for this product can be found in chapter 8.

## Table of Contents

1. INTENDED PURPOSE _____	18
2. SUMMARY AND EXPLANATION OF THE TEST _____	18
3. PRINCIPLE OF THE TEST _____	18
4. MATERIAL SUPPLIED _____	19
5. MATERIAL REQUIRED BUT NOT SUPPLIED _____	20
6. STORAGE AND PREPARATION OF REAGENTS _____	20
7. SPECIMEN COLLECTION AND PREPARATION _____	20
<i>Dried blood spots</i> _____	21
8. PRECAUTIONS _____	21
9. GENERAL NOTES ON THE TEST _____	21
10. REFERENCES _____	22
11. SYMBOL EXPLANATION _____	22

## 1. INTENDED PURPOSE

This HPLC application is intended for the quantitative determination of vitamin B<sub>6</sub> in plasma, serum, EDTA whole blood and dried blood spots (DBS). For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin B<sub>6</sub> is the term for pyridoxin, pyridoxal and pyridoxamin and the appropriate phosphate products. All forms can be transformed into the active form pyridoxal-5-phosphate (here referred to as vitamin B<sub>6</sub>).

Vitamin B<sub>6</sub> functions as a coenzyme and is essential for more than 50 reactions in the protein metabolism thereby synthesizing, transforming or degrading amino acids. Vitamin B<sub>6</sub> supports the resorption of amino acids and their transport into the cells. Furtheron vitamin B<sub>6</sub> contributes to the synthesis of neuro transmitters and amine products (histamin).

Due to the fact that vitamin B<sub>6</sub> contributes to a variety of different reactions lack of vitamin B<sub>6</sub> results in various clinical pictures.

### Applications:

- determination of vitamin B<sub>6</sub> status
- homocysteinaemia
- skin diseases
- movement disorders
- Anaemia, depression

## 3. PRINCIPLE OF THE TEST

The first step in the determination of vitamin B<sub>6</sub> includes the sample preparation with additional derivatisation. During the precipitation, higher molecular substances are removed. The supernatant is used for derivatisation (5 min at 60°C), thereby transforming the vitamin B<sub>6</sub> into a fluorescent product. The sample is then injected into the HPLC system.

The separation via HPLC follows an isocratic method at 22–35 °C using a reversed phase column; one run lasts about 2.5 minutes. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak area/heights.

## Summary

This HPLC application for the determination of vitamin B<sub>6</sub> allows to determine the vitamin in an easy, fast and precise method. The kit includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As for many other parameters, the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical chemical routines quickly and precisely. Mostly, a one-point calibration is sufficient for calibrating the test system, unlike immunoassays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. With short test series, the one-point calibration is much more economic than 6-point calibration for immunoassays.

## 4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 0005.15	RECSOL	Reconstitution solution	1 x 15 ml
KC0003.40	PREC	Precipitation reagent (serum, plasma, whole blood, DBS)	1 x 40 ml
KC2151	MOPHA	Mobile phase, ready to use (can be recirculated)	1 x 1 000 ml
	CAL	Calibrator, lyophilised	4x
	DER	Derivatisation solution (contains KCN)	1 x 26 ml
	CTRL1	Control 1, lyophilised	4x
	CTRL2	Control 2, lyophilised	4x

For reorders of single components, use the catalogue number followed by the label as product number.

The HPLC column (KC2151RP), can be ordered separately from Immundiagnostik. To extend the lifetime of your HPLC column, pre-columns (KC2151VS) are highly recommended. These and also the pre-column holders (KC2151VH) can also be ordered from Immundiagnostik. In addition to the complete kits, all components can also be ordered separately. Please ask for our single component price list.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (e.g. Eppendorf)
- Centrifuge
- Various pipettes
- Vortex
- HPLC with fluorescence detector
- Reversed phase C<sub>18</sub> PFP column
- Thermoblock
- Ultrasonic bath

## 6. STORAGE AND PREPARATION OF REAGENTS

- The lyophilised calibrator (**CAL**) is stable at **-20°C** until the expiry date stated on the label. Before use, the CAL has to be reconstituted in x µl (x = see the enclosed product specification for the volume needed) reconstitution solution (**RECSOL**). The concentration of vitamin B<sub>6</sub> slightly changes from lot to lot, please refer to the enclosed product specification sheet for the exact amount.
- The lyophilised controls (**CTRL1 and CTRL2**) are stable at **-20°C** until the expiry date stated on the label. Before use, they have to be reconstituted with each x µl (x = see the enclosed product specification for the volume needed) reconstitution solution (**RECSOL**). The concentration of vitamin B<sub>6</sub> slightly changes from lot to lot, the exact concentration is stated on the specification data sheet.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8°C**.

## 7. SPECIMEN COLLECTION AND PREPARATION

Plasma, serum, EDTA-whole blood, CPDA-whole blood and dried blood spots are suitable for use in the assay.

### Sample storage

Vitamin B<sub>6</sub> is light- and temperature-sensitive; therefore, samples have to be protected from light, cooled and centrifuged immediately.

The samples are stable in the dark at 2–8°C for one week. For longer storage, samples can be frozen and kept at -20°C for up to 6 months.

## Dried blood spots

### Collection and storage of dried blood spots

50 µl whole blood dripped on a dried sample carrier cleared by Immundiagnostik AG are suitable as sample material after complete drying. We recommend DrySpot-ID (catalogue no. **DZ9020ID** or **DZ9021ID**) as dried blood spot carrier. The moistened cards are stable for 7 days at room temperature. For longer storage, store at -20 °C in a dry place.

## 8. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- As the derivatisation solution (**DER**) contains KCN, it should be pipetted under a fume hood.
- The precipitation reagent consists of an acid. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.

## 9. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Do not use reagents beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Follow the guidelines for medical laboratories.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be lodged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Analyse controls with each run.
- Always perform the assay according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities. Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- The assay should always be performed according to the enclosed manual.

## 10. REFERENCES

1. Ambrosch, A. et al., 2001. Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic medicine : a journal of the British Diabetic Association*, **18**(3), pp.185–92.
2. Dierkes, J., Domröse, U., et al., 2001. Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *Journal of renal nutrition : the official journal of the Council on Renal Nutrition of the National Kidney Foundation*, **11**(2), pp.67–72.
3. Dierkes, J., Westphal, S., et al., 2001. Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis*, **158**(1), pp.161–4.

## 11. SYMBOL EXPLANATION

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)