Porphyrins HPLC Kit

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
Porphyrins HPLC Kit

For the determination of porphyrins (free acids) in urine

Valid from 2023-01-12

REF **KC2601**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

| | |
|---|-----------|
| 1. INTENDED USE | 15 |
| 2. SUMMARY AND EXPLANATION OF THE TEST | 15 |
| 3. PRINCIPLE OF THE TEST | 15 |
| 4. MATERIAL SUPPLIED | 16 |
| 5. MATERIAL REQUIRED BUT NOT SUPPLIED | 17 |
| 6. STORAGE AND PREPARATION OF REAGENTS | 17 |
| 7. SPECIMEN COLLECTION AND PREPARATION | 17 |
| 8. ASSAY PROCEDURE | 18 |
| <i>Test procedure</i> | 18 |
| <i>Chromatographic conditions</i> | 18 |
| 9. TREATMENT OF THE COLUMN | 19 |
| 10. RESULTS | 19 |
| <i>Calculation</i> | 19 |
| <i>Typical chromatogram</i> | 19 |
| 11. LIMITATIONS | 20 |
| 12. QUALITY CONTROL | 20 |
| <i>Reference ranges</i> | 20 |
| 13. PERFORMANCE CHARACTERISTICS | 21 |
| <i>Precision and reproducibility</i> | 21 |
| <i>Linearity</i> | 21 |
| <i>Detection limit</i> | 21 |
| 14. DISPOSAL | 22 |
| 15. TROUBLESHOOTING | 22 |
| 16. PRECAUTIONS | 23 |
| 17. GENERAL NOTES ON THE TEST | 24 |
| 18. REFERENCES | 24 |

1. INTENDED USE

This HPLC application is intended for the quantitative determination of porphyrins in urine. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

The porphyrins are precursors of heme, which is a cofactor in hemoglobin, myoglobin and the cytochromes. They play an important role in the oxygen metabolism. Heme biosynthetic activity is located quantitatively in bone marrow and liver. The series of reactions leading to heme synthesis begins with the condensation of succinyl coenzyme A (CoA) and glycine and ends with the insertion of an iron atom into a molecule of protoporphyrin IX.

There are several genetic defects of heme biosynthesis described. In all of them an increase of porphyrins is recognized.

The laboratory diagnosis and classification of porphyrin disorders are based primarily on the measurement of the excretion of porphyrins and porphyrin precursors, mainly in urine. Uro-, copro-, penta- and tricarboxyporphyrin are increased in autosomal dominant acute porphyrin disorders. In chronic porphyrin disorders, including *porphyria cutanea tarda*, uro- and hepta-porphyrin are elevated. In chronic lead intoxication the coproporphyrins are slightly increased (0.5–2 $\mu\text{mol}/24\text{ h}$ urine), acute lead poisoning is characterized by extremely high excretion of total porphyrins (up to 15 $\mu\text{mol}/24\text{ h}$ urine).

Indications

- Hereditary hepatic porphyria (e. g. *porphyria variegata*)
- *Porphyria cutanea tarda* in chronic liver disease and alcoholic liver syndrome
- Hepatic porphyria in prostate and liver tumors
- Acute and chronic lead poisoning

3. PRINCIPLE OF THE TEST

The first step in the determination of porphyrins includes an easy sample preparation. The pH value of the sample and the calibrator is adjusted below 2.5 by addition of 20 μl of hydrochloric acid. After centrifugation, 100 μl of the supernatant is injected into the HPLC system.

The separation via HPLC follows a gradient method, using a reversed phase column at 30 °C; one run lasts 25 minutes. The chromatograms are recorded with a fluorescence detector. The quantification is performed with the delivered calibrator. The concentration is calculated via integration of the peak areas or peak heights.

Summary

This HPLC application provides an easy, fast and precise method for quantitative determination of porphyrins. The kit contains all reagents necessary for sample preparation and separation except the column.

As for many other parameters, the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical chemical routines quickly and precisely. Mostly, a one-point calibration is sufficient for calibrating the test system – unlike immunoassays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. With short test series, the one-point calibration is much more economic than 6-point calibration for immunoassays.

4. MATERIAL SUPPLIED

| Cat. No. | Label | Kit components | Quantity |
|-----------|--------|--|--------------|
| K 0005.15 | RECSOL | Reconstitution solution | 1 x 15 ml |
| KC2601 | MOPHAA | Mobile phase A | 1 000 ml |
| | MOPHAB | Mobile phase B | 2 x 1 000 ml |
| | CAL | Calibrator; lyophilised (see specification data sheet for concentration) | 4 x |
| | HCL | Hydrochloric acid | 3 ml |
| | CTRL1 | Control1; lyophilised | 4 x |
| | CTRL2 | Control2; lyophilised | 4 x |

For reorders of single components, use the catalogue number followed by the label as product number.

The HPLC column (KC2601RP), can be ordered separately from Immundiagnostik. To extend the lifetime of your HPLC column, pre-columns (KC2601VS) are highly recommended. These and also the pre-column holders (KC2601VH) can also be ordered from Immundiagnostik.

In addition to the complete kits, all components can also be ordered separately. Please ask for our single component price list.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (e.g. Eppendorf)
- Centrifuge
- Various pipettes
- Vortex
- HPLC gradient pump with fluorescence detector
- Reversed phase C₁₈ column

6. STORAGE AND PREPARATION OF REAGENTS

- **The lyophilised calibrator (CAL)** is stable at **-20°C** until the expiry date stated on the label. Before use, the **CAL** has to be reconstituted with **x µl (x = see the enclosed product specification for the volume needed) reconstitution solution (RECSOL)**. The concentration of porphyrins slightly changes from lot to lot, the exact concentration is stated on the specification data sheet.
- **The lyophilised controls 1 and 2 (CTRL1 and 2)** are stable at **-20°C** until the expiry date stated on the label. Before use, they have to be reconstituted each with **x µl (x = see the enclosed product specification for the volume needed) reconstitution solution (RECSOL)**. The concentration of porphyrins slightly changes from lot to lot, the exact concentration is stated on the specification data sheet.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2–8°C**.

7. SPECIMEN COLLECTION AND PREPARATION

Urine is suited for this test system.

The porphyrins are light and temperature-sensitive. Therefore, protect samples from light and cool immediately after collection.

8. ASSAY PROCEDURE

Test procedure

| | |
|---|---|
| Add 20 µl hydrochloric acid (HCL) to the calibrator and control 1 or 2 in the original flask and mix well. Pipet 250 µl into a reaction tube. | Mix 250 µl of a cool and dark stored 24 h urine (without additives), with 20 µl hydrochloric acid (HCL) (pH should be below 2.5) and mix well. |
| Centrifuge for 5 min at 10000 g . | |
| Inject 100 µl of the supernatant into the HPLC | |

Chromatographic conditions

| | | |
|-------------------------|--|--------|
| Column material: | Bischoff Prontosil 120-5-C18 ace EPS; 5 µm | |
| Column dimension: | 125 mm × 4 mm | |
| Flow rate: | 0.75 ml/min | |
| Fluorescence detection: | Excitation: | 400 nm |
| | Emission: | 620 nm |
| Temperature: | 30 °C | |
| Injection volume: | 100 µl | |
| Run time: | 25 minutes | |
| Gradient: | | |

| Time | % B | % A |
|----------|-----|-----|
| 0 min | 38 | 62 |
| 1.5 min | 38 | 62 |
| 9.5 min | 80 | 20 |
| 11.5 min | 80 | 20 |
| 11.6 min | 95 | 5 |
| 15 min | 95 | 5 |
| 15.1 min | 38 | 62 |
| 24 min | 38 | 62 |

9. TREATMENT OF THE COLUMN

It is recommended to use a guard column to extend column life.

After analysis, the column should be flushed with 30 ml ultra pure water (1 ml/min) and stored in 50% methanol in water (~ 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ~ 30 ml mobile phase (MOPHA).

10. RESULTS


Calculation

$$\text{Sample concentration} = \frac{\text{Peak height sample} \times \text{calibrator concentration}^*}{\text{Peak height calibrator}}$$

* see specification data sheet

Tip: Alternatively, the peak area instead of the peak height can be used for quantification.

Typical chromatogram

Column: Prontosil 120-5-C18 ace-EPS

Order of peaks (from the left to the right)

- 1) 8-carboxyl Porphyrin (Uro-)
- 2) 7-carboxyl Porphyrin (Hepta-)
- 3) 6-carboxyl Porphyrin (Hexa-)
- 4) 5-carboxyl Porphyrin (Penta)
- 5) 4-carboxyl Porphyrin (Copro) I
- 6) 4-carboxyl Porphyrin (Copro) III
- 7) Mesoporphyrin IX

11. LIMITATIONS

Do not use blood samples for analysis.

12. QUALITY CONTROL

Reference ranges

The normal range is related to 24 h urine.

| | Mean ($\mu\text{g}/24\text{ h}$) | Area ($\mu\text{g}/24\text{ h}$) |
|--|------------------------------------|------------------------------------|
| 8-carboxyl porphyrin (Uro-) | 7.0 | 0–33 |
| 7-carboxyl porphyrin (Hepta-) | < detection limit | 0–10 |
| 6-carboxyl porphyrin (Hexa-) | < detection limit | 0–7 |
| 5-carboxyl porphyrin (Penta) | < detection limit | 0–5 |
| 4-carboxyl porphyrin (Copro) Coprop. I 17 – 31 % (Lit. 3) Coprop. III 69 – 83 % (Lit. 3) | 44.8 | 0–120 |

We recommend each laboratory to establish its own reference range. The above mentioned values are only for orientation and may deviate from other published data.

Controls

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV:

| | | |
|--------------------|-------------------|----------|
| Uroporphyrin | 5.6 % (35.8 µg/l) | [n = 10] |
| Heptaporphyrin | 2.5 % (25.3 µg/l) | [n = 10] |
| Hexaporphyrin | 2.6 % (27.0 µg/l) | [n = 10] |
| Pentaporphyrin | 2.6 % (31.7 µg/l) | [n = 10] |
| Coproporphyrin I | 2.5 % (24.9 µg/l) | [n = 10] |
| Coproporphyrin III | 2.6 % (57.1 µg/l) | [n = 10] |

Inter-Assay CV:

| | | |
|--------------------|-------------------|---------|
| Uroporphyrin | 6.6 % (52.5 µg/l) | [n = 5] |
| Heptaporphyrin | 6.0 % (30.1 µg/l) | [n = 5] |
| Hexaporphyrin | 5.0 % (30.8 µg/l) | [n = 5] |
| Pentaporphyrin | 5.1 % (35.9 µg/l) | [n = 5] |
| Coproporphyrin I | 7.5 % (24.0 µg/l) | [n = 5] |
| Coproporphyrin III | 4.7 % (65.9 µg/l) | [n = 5] |

Linearity

up to 1 mg/l

Detection limit

| | |
|--------------------|-----------|
| Uroporphyrin | 0.60 µg/l |
| Heptaporphyrin | 0.55 µg/l |
| Hexaporphyrin | 0.49 µg/l |
| Pentaporphyrin | 0.51 µg/l |
| Coproporphyrin I | 0.51 µg/l |
| Coproporphyrin III | 0.50 µg/l |

14. DISPOSAL

The mobile phase (MOPHA) must be disposed as non-halogenated solvent. The hydrochloric acid (HCL) solution can be neutralized with NaOH and if the pH value is neutral it can be disposed of as salt solution.

Important: Reaction will produce heat, be careful!

Please refer to the appropriate national guidelines.

15. TROUBLESHOOTING

| Problem | Possible cause | Solution |
|--------------------------|--|---|
| No signal | No or defect connection to evaluation system | Check signal cord and connection |
| | Detector lamp is too old | Change lamp |
| No peaks | Injector is congested | Check injector |
| Double peaks | Dead volume in fittings and / or column | Renew fittings and / or column |
| Contaminating peaks | Injector dirty | Clean injector |
| | Contamination at the head of the column | Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase |
| | Air in the system | Degas pump |
| | Auto sampler vials contaminated | Use new vials or clean them with methanol |
| Broad peaks, tailing | Precolumn / column exhausted | Use new precolumn / column |
| Variable retention times | Drift in temperature | Use a column oven |
| | Pump conveys inaccurately | Check pump, degas the system |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |

| Problem | Possible cause | Solution |
|------------------------|---|--------------------------------------|
| Baseline is drifting | Detector lamp did not reach working temperature yet | Wait |
| | Detector lamp is too old | Change lamp |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |
| | Pump conveys inaccurately | Check pump, degas the system |
| Baseline is not smooth | Pump conveys inaccurately | Check pump, degas the system |
| | Detector flow cell is dirty | Clean flow cell |

16. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- The hydrochloric acid (HCL) solution contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- Reagents should not be used beyond the expiration date shown on kit label.


17. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The assay should always be performed according to the enclosed manual.
- Plugs and caps of different reagents should not be swapped.
- Quality control guidelines should be observed.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

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Used symbols:

| | | | |
|---|--|---|------------------------------------|
|  | Temperature limitation |  | Catalogue number |
|  | In Vitro Diagnostic Medical Device |  | To be used with |
|  | Manufacturer |  | Contains sufficient for <n> tests |
|  | Lot number |  | Use by |
|  | Contains plasma derivatives or human blood |  | Consult instructions for use |
|  | Consult specification data sheet |  | Do not re-use |
|  | Unique Device Identification |  | Contains material of animal origin |
|  | Medicinal substance |  | Contains material of human origin |

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com