# MutaREAL<sup>®</sup> ApoE

## Real-Time-PCR-Kit

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# MutaREAL<sup>®</sup> ApoE

## Real-Time-PCR Kit

*For the analysis of the amino acid substitution  
Cys130Arg and Cys176Arg in the ApoE gene*

Valid from 2024-10-31

**KF290732**  
**KF290796**

32/96

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1</b>	<b>INTENDED USE</b>	<b>17</b>
<b>2</b>	<b>INTRODUCTION</b>	<b>17</b>
<b>3</b>	<b>PRINCIPLE OF THE TEST</b>	<b>17</b>
<b>4</b>	<b>PACKAGE CONTENTS</b>	<b>18</b>
<b>5</b>	<b>EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER</b>	<b>18</b>
<b>6</b>	<b>TRANSPORT, STORAGE AND STABILITY</b>	<b>18</b>
<b>7</b>	<b>WARNINGS AND PRECAUTIONS</b>	<b>19</b>
<b>8</b>	<b>SAMPLE MATERIAL</b>	<b>20</b>
<b>9</b>	<b>REAL-TIME-PCR</b>	<b>20</b>
	9.1 <i>Important points before starting</i>	20
	9.2 <i>Procedure</i>	20
	9.3 <i>Instrument settings</i>	22
<b>10</b>	<b>DATA ANALYSIS</b>	<b>22</b>
	10.1 <i>Evaluation and classification in ApoE isoforms</i>	24
<b>11</b>	<b>TROUBLESHOOTING</b>	<b>24</b>
<b>12</b>	<b>KIT PERFORMANCE</b>	<b>25</b>
	12.1 <i>Analytical Sensitivity</i>	25
	12.2 <i>Precision</i>	25
	12.3 <i>Diagnostic Sensitivity</i>	25
<b>13</b>	<b>LIMITATIONS OF THE METHOD</b>	<b>26</b>
<b>14</b>	<b>ABBREVIATIONS AND SYMBOLS</b>	<b>27</b>
<b>15</b>	<b>LITERATURE</b>	<b>27</b>

## 1 INTENDED USE

The MutaREAL® ApoE Real-Time PCR Kit is a FRET-based molecular biology test for the investigation of point mutations in the ApoE gene, which are responsible for the exchange of the amino acids Cys130Arg (formerly pos. 112) and Cys176Arg (formerly pos. 158).

## 2 INTRODUCTION

ApolipoproteinE (ApoE) belongs to the family of apolipoproteins, which play an important role in lipid metabolism and influences the formation of very-low-density lipoproteins (VLDL). The most common ApoE alleles are ApoE2 (E2), ApoE3 (E3) and ApoE4 (E4). These are defined by the variations Cys130Arg and Cys176Arg. E3 carries cysteine at position 130 and arginine at 176, E2 carries cysteine at both positions, whereas E4 carries arginine at both positions. The different ApoE isoforms have an influence on cholesterol levels. [1]

## 3 PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® ApoE real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay includes two specific primers flanking the target sequence and two hybridisation probes that bind adjacent to the target sequence. One of the hybridisation probes is labelled with a donor fluorophore and, after appropriate excitation, transfers its energy to the acceptor fluorophore with which the other hybridisation probe is labelled when they are in close proximity. After the energy transfer, the acceptor dye emits light with a longer wavelength. Energy transfer can only occur when both hybridisation probes have bound to the target sequence. The amount of hybridised probe pairs and thus the fluorescence signal increases with the amount of amplified PCR product. Here, the fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of amplification by melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. One of the hybridisation probes binds to a part of the target sequence that is present in the wild type and the mutation. The second hybridisation probe spans the mutation site. As the temperature rises, the mismatched and thus less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and thus the fluorescence signal decreases only at a higher temperature.

## 4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF290732) or 96 (KF290796) reactions.

Table 1: Components of the MutaREAL® ApoE Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 876 µl	2 x 1313 µl
Detection mix 1 (Cys130Arg)	yellow	1 x 368 µl	3 x 368 µl
Detection mix 2 (Cys176Arg)	brown	1 x 368 µl	3 x 368 µl
Positive control	red	1 x 30 µl	1 x 90 µl
Negative control	green	1 x 150 µl	1 x 150 µl

## 5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Roche LightCycler® 1.5, 2.0, 480 II or PRO real-time PCR system
  - \* CE conformity only exists if one of the mentioned devices is used.
- DNase/RNase-free Roche LightCycler® capillaries or 96-well plates/strips (white)
- Roche LightCycler® Cooling Block
- Sterile reaction tubes
- Calibrated precision pipettes with variable volumes from 2 to 1000 µl
- Sterile single-use pipette tips with aerosol barrier
- (Table) centrifuge (without specific requirements)
- Vortex mixer (without specific requirements)
- Water, PCR grade
- EDTA plasma tubes for sample collection

## 6 TRANSPORT, STORAGE AND STABILITY

The MutaREAL® ApoE real-time PCR kit is transported frozen on dry ice. All components are to be stored protected from light at a minimum of -20°C immediately after receipt. Avoid more than five freeze-thaw cycles of the detection mixes (make aliquots if necessary). The components can be stored at 4–8°C after initial thawing for up to three months. Do not use after the expiry date indicated on the package.

Be sure to protect the detection mixes from direct sunlight during the entire test period.

## 7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

## 8 SAMPLE MATERIAL

Starting material for the MutaREAL® ApoE real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

## 9 REAL-TIME-PCR

### 9.1 Important points before starting

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8°C) or on ice while working.

### 9.2 Procedure

For amplification, two reaction tubes are needed per sample plus two additional reaction tubes per master mix for the negative and the positive controls. The following table shows the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10% to compensate for inaccuracies. The master mix is pipetted as described in tables 2 and 3:

#### Master mix 1 (Cys130Arg)

Table 2: Preparation of master mix 1 (Cys130Arg)

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 1 (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary / each reaction tube.
- For the negative control, add **2 µl** of the supplied negative control (**green**).
- For the positive control, add **2 µl** of the supplied positive control (**red**).
- For each sample to be analyzed, add **2 µl** of the extracted genomic DNA to the corresponding capillary / reaction tube.

### Master mix 2 (Cys176Arg)

Table 3: Preparation of master mix 2 (Cys176Arg)

Reagenz	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 2 ( <b>brown</b> )	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix ( <b>blue</b> )	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary / each reaction tube.
- For the negative control, add **2 µl** of the supplied negative control (**green**).
- For the positive control, add **2 µl** of the supplied positive control (**red**).
- For each sample to be analyzed, add **2 µl** of the extracted genomic DNA to the corresponding capillary / reaction tube.

LightCycler® 1.5 / 2.0: Close the capillaries with the lids, transfer them into the LightCycler® carousel and centrifuge them in the LightCycler® centrifuge (if a tabletop centrifuge is used, centrifuge the capillaries in the inserts of the cooling block at 3000 rpm for 15 s). Then transfer the carousel to the LightCycler® device and start the PCR program described in 9.3.

LightCycler® 480 II /PRO: Seal the wells with a sealing foil and centrifuge the plate (2000 rpm for 15 s). Then transfer the plate into the LC 480 II / PRO and start the PCR program described in 9.3.

### 9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 4.

For the Roche LightCycler® 480 II the detection format *'Mono Color HybProbe'* and the filter *'Red 640 (498-640)'* must be selected.

Table 4: Real-Time-PCR thermal profile

Description	Time	Temperature	Heating rate	Cycles	Acquisition
Initial Denaturation	180 s	95 °C	max.	1	none
Denaturation	10 s	95 °C	max.	45	none
Primer annealing	25 s	55 °C	max.		<b>single</b>
Elongation	25 s	72 °C	max.		none
Melting curve	60 s	95 °C	max.	1	none
	30 s	60 °C	max.	1	none
	30 s	50 °C	max.	1	none
	30 s	40 °C	max.	1	none
	0 s	80 °C	0.1 °C/s	1	<b>constant</b>
Cooling	30 s	40 °C	max.	1	-

\* Roche LightCycler 480 II: 6 Acquisitions per °C

## 10 DATA ANALYSIS

Add an analysis of the type *'Tm calling'* for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is 640 nm.

### ApoE Cys130Arg

Temperature wild type allele (Cys)\*: 51.5 °C (+/-2 °C)

Temperature mutation allele (Arg)\*: 61.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **pink curve** - negative control, **red curve** - homozygous **Arg / Arg**, **blue curve** - homozygous **Cys / Cys**, **green curve** - heterozygous **Cys / Arg**.


Fig. 1: Evaluation of master mix 1 - Cys130Arg.

### ApoE Cys176Arg

Temperature wild type allele (Cys)\*: 51.0°C (+/-2°C)

Temperature mutation allele (Arg)\*: 61.0°C (+/-2°C)

The following graph shows the typical results for the possible genotypes: **yellow curve** - negative control, **blue curve** - homozygous **Cys / Cys**, **red curve** - homozygous **Arg / Arg**, **grey curve** - heterozygous **Cys / Arg**.


Fig. 2: Evaluation of master mix 2 - Cys176Arg.

\* Depending on the PCR device used, the  $T_m$  values determined may deviate by more than +/- 2°C. In this case use the positive control (red) as a reference and apply +/- 2°C to these  $T_m$  values.

The provided positive control contains a template which is **heterozygous** for both mutations (ApoE E2/E4).

## 10.1 Evaluation and classification in ApoE isoforms

Table 5: Evaluation and classification in ApoE isoforms

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## 11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on [info@immundiagnostik.com](mailto:info@immundiagnostik.com).

### No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than five freeze cycles or have been stored at 2-8°C for more than three month. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

## 12 KIT PERFORMANCE

### 12.1 Analytical Sensitivity

The limit of detection (LoD) of the MutaREAL® ApoE Real-Time PCR kit was determined using serial dilutions of DNA-fragments containing the mutation target sequences in a Roche LightCycler® 2.0 real time PCR instrument. The LoD of the MutaREAL® ApoE Real-Time PCR kit is 1000 genome copies per reaction each.

### 12.2 Precision

The precision of the MutaREAL® ApoE Real-Time PCR kit was determined as intra-assay variability and inter-assay variability.

Variability data are expressed by standard deviation and coefficient of variation. The data are based on  $T_m$  calculation analyses of defined concentrations of the ApoE mutation target sequences.

<b>ApoE Cys130Arg</b>	<b>copies/ µl</b>	<b>Standard deviation</b>	<b>Coefficient of variation [%]</b>
Intra-Assay Variability	3000	1.65	2.19
Inter-Assay Variability	3000	2.97	3.78

<b>ApoE Cys176Arg</b>	<b>copies/ µl</b>	<b>Standard deviation</b>	<b>Coefficient of variation [%]</b>
Intra-Assay Variability	3000	2.54	4.49
Inter-Assay Variability	3000	2.37	5.89

### 12.3 Diagnostic Sensitivity

The diagnostic sensitivity of Real-Time-PCR assays is mainly dependent on the DNA extraction method used to isolate DNA from various biological specimens. DNA extraction reagents are not part of the Immundiagnostik AG Real-Time-PCR kits. Immundiagnostik AG Real-Time-PCR kits include a positive and negative control. Therefore, Immundiagnostik AG guarantees the analytical sensitivities and specificities of the Real-Time PCR kits, performed with eluted DNA from reference materials and ring trial samples and with synthetic nucleic acid fragments. Immundiagnostik AG does not guarantee diagnostic sensitivities. If diagnostic sensitivities are mentioned in manuals of Immundiagnostik AG Real-Time-PCR kits, the data are strictly correlated to a specific nucleic acid extraction method that has been used during the validation of the respective kits and cannot be transferred to other extraction methods.

It is the responsibility of the user to qualify the extraction methods used for DNA isolation.


### **13 LIMITATIONS OF THE METHOD**

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

## 14 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix 1		Upper limit of temperature
	Detection mix 2		Manufacturer
	Positive control		Lot number
	Negative control		Consult instructions for use
	<i>In vitro</i> diagnostic medical device		Content
	Use by YYYY-MM-DD		

## 15 LITERATURE

1. Eichner et al., Am J Epidemiol, 2002, 155(6):487-95

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)