# MutaREAL<sup>®</sup> TPMT

## Real-Time-PCR-Kit

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# MutaREAL<sup>®</sup> TPMT

## Real-Time-PCR Kit

*For the analysis of the allele variants  
\*2, \*3A, \*3B and \*3C of the TPMT gene.*

Valid from 2023-02-22

**KF291332**  
**KF291396**

32/96

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1</b>	<b>INTENDED USE</b>	<b>15</b>
<b>2</b>	<b>INTRODUCTION</b>	<b>15</b>
<b>3</b>	<b>PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>4</b>	<b>PACKAGE CONTENTS</b>	<b>16</b>
<b>5</b>	<b>EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER</b>	<b>16</b>
<b>6</b>	<b>TRANSPORT, STORAGE AND STABILITY</b>	<b>16</b>
<b>7</b>	<b>WARNINGS AND PRECAUTIONS</b>	<b>17</b>
<b>8</b>	<b>SAMPLE MATERIAL</b>	<b>17</b>
<b>9</b>	<b>REAL-TIME-PCR</b>	<b>18</b>
	9.1 <i>Important points before starting</i>	18
	9.2 <i>Procedure</i>	18
	9.3 <i>Instrument settings</i>	20
<b>10</b>	<b>DATA ANALYSIS</b>	<b>20</b>
	10.1 <i>Special cases - TPMT G460A and A719G</i>	22
<b>11</b>	<b>TROUBLESHOOTING</b>	<b>22</b>
<b>12</b>	<b>LIMITATIONS OF THE METHOD</b>	<b>23</b>
<b>13</b>	<b>ABBREVIATIONS AND SYMBOLS</b>	<b>23</b>
<b>14</b>	<b>LITERATURE</b>	<b>24</b>

## 1 INTENDED USE

The MutaREAL® TPMT Real-Time PCR Kit is a FRET-based molecular biology test for the investigation of alleles \*2, \*3A, \*3B and \*3C (corresponding to point mutations G238C, G460A and A719G) of the TPMT gene.

## 2 INTRODUCTION

Thiopurine-based drugs, such as the immunosuppressant azathioprine and the cancer drugs mercaptopurine and thioguanine, are often used for the treatment of chronic inflammatory diseases. The enzyme thiopurine methyltransferase (TPMT) is responsible for the methylation of thiopurines. This step is part of the biotransformation.

The mutations G238C, G460A and A719G (\*2, \*3A, \*3B and \*3C) are associated with reduced TPMT enzyme activity. Low TPMT enzyme activity leads to an accumulation of toxic degradation products after ingestion of thiopurine-containing drugs, which can cause severe side effects such as myelosuppression. [1], [2], [3]

## 3 PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® TPMT real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay includes two specific primers flanking the target sequence and two hybridisation probes that bind adjacent to the target sequence. One of the hybridisation probes is labelled with a donor fluorophore and, after appropriate excitation, transfers its energy to the acceptor fluorophore with which the other hybridisation probe is labelled when they are in close proximity. After the energy transfer, the acceptor dye emits light with a longer wavelength. Energy transfer can only occur when both hybridisation probes have bound to the target sequence. The amount of hybridised probe pairs and thus the fluorescence signal increases with the amount of amplified PCR product. Here, the fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of amplification by melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. One of the hybridisation probes binds to a part of the target sequence that is present in the wild type and the mutation. The second hybridisation probe spans the mutation site. As the temperature rises, the mismatched and thus less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and thus the fluorescence signal decreases only at a higher temperature.

## 4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF291332) or 96 (KF291396) reactions.

Table 1: Components of the MutaREAL® TPMT Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Enzyme mix	blue	1 x 1313 µl	3 x 1313 µl
Detection mix 1 (TPMT*2, G238C)	yellow	1 x 368 µl	3 x 368 µl
Detection mix 2 (TPMT*3B, G460A)	brown	1 x 368 µl	3 x 368 µl
Detection mix 3 (TPMT*3C, A719G)	purple	1 x 368 µl	3 x 368 µl
Positive control 1 (TPMT*2 heterozygous)	red	1 x 25 µl	1 x 75 µl
Positive control 2 (TPMT*3B heterozygous)	orange	1 x 25 µl	1 x 75 µl
Positive control 3 (TPMT*3C heterozygous)	white	1 x 25 µl	1 x 75 µl
Negative control	green	1 x 150 µl	1 x 150 µl

## 5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Roche LightCycler® 1.5 or 2.0 real-time PCR system
  - \* CE conformity only exists if one of the mentioned devices is used.
- Roche LightCycler® capillaries
- Roche LightCycler® Cooling Block
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

## 6 TRANSPORT, STORAGE AND STABILITY

The MutaREAL® TPMT real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

## 7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

## 8 SAMPLE MATERIAL

Starting material for the MutaREAL® TPMT real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

## 9 REAL-TIME-PCR

### 9.1 Important points before starting

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

### 9.2 Procedure

For amplification, three reaction tubes (LightCycler® capillaries) are needed per sample plus two additional reaction tubes per master mix for the negative and the positive controls. The following tables show the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mixes are pipetted as described in tables 2, 3 and 4:

#### Master mix 1 - TPMT\*2

Table 2: Preparation of master mix 1 - TPMT\*2

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 1 (yellow)	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix (blue)	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (green).

- For the positive control add **2 µl** of the supplied positive control 1 (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

### Master mix 2 - TPMT\*3B

Table 3: Preparation of master mix 2 - TPMT\*3B

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 2 ( <b>brown</b> )	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix ( <b>blue</b> )	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 2 (**orange**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

### Master mix 3 - TPMT\*3C

Table 4: Preparation of master mix 3 - TPMT\*3C

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Detection mix 3 ( <b>purple</b> )	10.5 µl	10.5 µl * (N + (N * 0.1))
Enzyme mix ( <b>blue</b> )	12.5 µl	12.5 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 3 (**white**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

Close the capillaries with the lids, transfer them into the LightCycler® carousel and centrifuge them in the LightCycler® centrifuge (if a tabletop centrifuge is used, centrifuge the capillaries in the inserts of the cooling block at 3000 rpm for 15 s). Then transfer the carousel to the LightCycler® and start the PCR programme described in 9.3.

### 9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 5.

Table 4: Real-Time-PCR thermal profile

Description	Time	Temperature	Heating rate	Cycles	Acquisition
Initial Denaturation	120 s	94 °C	20 °C/s	1	none
Denaturation	10 s	94 °C	20 °C/s	45	none
Primer annealing	10 s	58 °C	20 °C/s		<b>single</b>
Elongation	8 s	72 °C	20 °C/s		none
Melting curve	20 s	95 °C	20 °C/s	1	none
	20 s	40 °C	20 °C/s	1	none
	0 s	75 °C	0.2 °C/s	1	<b>constant</b>
Cooling	30 s	40 °C	20 °C/s	1	-

## 10 DATA ANALYSIS

Add an analysis of the type „genotyping“ for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is 640 nm.

### TPMT \*2 - G238C

Temperature mutation allele: 56.0 °C (+/-2 °C)

Temperature wild type allele: 66.0 °C (+/-2 °C)

The following graph shows the typical results for TPMT\*2: **blue curve** - negative control, **red curve** - homozygous wild type, **green curve** - heterozygous mutation.


Fig. 1: Evaluation of TPMT \*2

### TPMT\*3B - G460A

Temperature mutation allele: 56.5°C (+/-2°C)

Temperature wild type allele: 64.0°C (+/-2°C)

The following graph shows the typical results for TPMT\*3B: **blue curve** - negative control, **green curve** - homozygous wild type, **black curve** - heterozygous mutation.


Fig. 2: Evaluation of TPMT\*3B.

### TPMT\*3C - A719G

Temperature wild type allele: 52.0°C (+/-2°C)

Temperature mutation allele: 58.0°C (+/-2°C)

The following graph shows the typical results for TPMT\*3C: **blue curve** - negative control, **beige curve** - homozygous wild type, **purple curve** - heterozygous mutation.


Fig. 3: Evaluation of TPMT\*3C.

### 10.1 Special cases - TPMT G460A and A719G

If the mutations are present on one allele, the allele is called TPMT\*3A. This results in the following special cases for the evaluation of the analysis:

TPMT\*3B and TPMT\*3C heterozygous: The genotype is either TPMT\*3B/\*3C or TPMT\*1/\*3A (if both mutations are present on one allele). Which genotype is present cannot be conclusively determined by this method.

TPMT\*3B heterozygous and TPMT\*3C homozygous mutation: The genotype is TPMT\*3A/\*3C.

TPMT\*3B homozygous mutant and TPMT\*3C heterozygous mutant: the genotype is TPMT\*3A/\*3B.

TPMT\*3B and TPMT\*3C homozygous: the genotype is TPMT\*3A/\*3A.

The positive control 1 (**red**) provided contains a template that is heterozygous for the point mutation G238C. Positive control 2 (**orange**) contains a template that is heterozygous for the mutation G460A. Positive control 3 (**white**) contains a template that is heterozygous for the mutation A719G.

## 11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on [info@immundiagnostik.com](mailto:info@immundiagnostik.com).

### No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8 °C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.


## 12 LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

## 13 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Enzyme mix		Contains sufficient for <n> test
	Detection mix 1		Upper limit of temperature
	Detection mix 2		Manufacturer
	Detection mix 3		Lot number
	Positive control 1		Consult instructions for use

<b>CONTROL 2</b>   +	Positive control 2	<b>CONTENT</b>	Content
<b>CONTROL 3</b>   +	Positive control 3		Use by YYYY-MM-DD
<b>CONTROL</b>   -	Negative control	<b>IVD</b>	<i>In vitro</i> diagnostic medical device

## 14 LITERATURE

1. Skrzypczak-Zielinska et al., Molecular Diagnosis & Therapy, 2016, 20:493-9
2. Yates et. al., Annals of Internal Medicine, 1997, 126(8):608-14
3. Tamm et al., Clinical Pharmacology & Therapeutics, 22. Oktober 2016, Accepted Article

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)