MutaREAL[®] DPD-2

Real-Time-PCR-Kit

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MutaREAL[®] DPD-2

Real-Time-PCR Kit

*For the analysis of the polymorphisms *2A, *13, *D949V and the haplotype B3 in the DPYD gene*

Valid from 2025-11-24

KF292032
KF292096

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	19
2. INTRODUCTION	19
3. PRINCIPLE OF THE TEST	19
4. PACKAGE CONTENTS	20
5. EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	20
6. TRANSPORT, STORAGE AND STABILITY	21
7. WARNINGS AND PRECAUTIONS	21
8. SAMPLE MATERIAL	22
9. REAL-TIME-PCR	22
9.1 <i>Important points before starting</i>	22
9.2 <i>Procedure</i>	23
9.3 <i>Instrument settings</i>	25
10. DATA ANALYSIS	25
10.1 <i>DPD*2A</i>	26
10.2 <i>DPD*Haplotyp B3</i>	27
10.3 <i>DPD*D949V</i>	28
10.4 <i>DPD*13</i>	29
11. TROUBLESHOOTING	30
12. LIMITATIONS OF THE METHOD	30
13. ABBREVIATIONS AND SYMBOLS	31
14. LITERATURE	31

1. INTENDED USE

The MutaREAL® DPD-2 Real-Time PCR Kit is a molecular biological test for the investigation of mutations and polymorphisms of the dihydropyrimidine dehydrogenase gene (DPYD) from genomic DNA. The test is based on FRET-based real-time PCR. The gene variants investigated are associated with dihydropyrimidine dehydrogenase deficiency, which has an influence on the metabolism of 5-Fluorouracil (5-FU). The following alleles and mutations are detected by the test:

Table 1: Detected alleles and mutations.

Allele / Mutation		rs number	Influence on the 5-FU metabolism
*2A	IVS14+1G>A	rs3918290	very high
*13	1679T>G	rs55886062	very high
D949V	2846A>T	rs67376798	very high
Haplotype B3	c.1129-5926C>G	rs75017182	high

2. INTRODUCTION

The enzyme dihydropyrimidine dehydrogenase (DPD) belongs to the group of oxidoreductases and catalyses the degradation of various cytostatic drugs such as 5-fluorouracil (5-FU). Mutations in the DPYD gene lead to reduced activity of the enzyme. This can lead to an accumulation of the cytostatic drug in the body during treatment with 5-FU. This can lead to severe side effects and life-threatening poisoning. [1]

3. PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® DPD-2 real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay contains specific primers that flank the respective target sequence and hybridisation probes that bind to the target sequence in the region of the SNP. One of the hybridisation probes is labelled with a donor fluorophore and, after appropriate excitation, transfers its energy to the acceptor fluorophore with which the other hybridisation probe is labelled if they are in close proximity. After the energy transfer, the acceptor dye emits light with a longer wavelength. Energy transfer can only occur when both hybridisation probes have bound to the target sequence. The amount of hybridised probe pairs and thus the fluorescence signal increases with the amount of amplified PCR product. Here, the fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of amplification by melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. One of the hybridisation probes binds to a part of the target sequence that is present in the wild type and the mutation. The second hybridisation probe spans the mutation site. As the temperature rises, the mismatched and thus less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and thus the fluorescence signal decreases only at a higher temperature.

4. PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF292032) or 96 (KF292096) reactions.

Table 2: Components of the MutaREAL® DPD-2 Real-Time-PCR Kit.

Label	Lid Colour	Content	
		32	96
Reaction mix 1 (DPD*2A)	yellow	1 x 805 µl	3 x 805 µl
Reaction mix 2 (DPD*Haplo B3)	white	1 x 805 µl	3 x 805 µl
Reaction mix 3 (DPD*D949V)	black	1 x 805 µl	3 x 805 µl
Reaction mix 4 (DPD*13)	colourless	1 x 805 µl	3 x 805 µl
Positive control 1 (DPD*2A)	red	1 x 30 µl	1 x 90 µl
Positive control 2 (DPD*Haplo B3)	orange	1 x 30 µl	1 x 90 µl
Positive control 3 (DPD*D949V)	purple	1 x 30 µl	1 x 90 µl
Positive control 4 (DPD*13)	brown	1 x 30 µl	1 x 90 µl
Negative control	green	1 x 150 µl	1 x 150 µl

5. EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Roche LightCycler® 1.5, 2.0, 480 II or PRO real-time PCR system
- DNase/RNase-free Roche LightCycler® capillaries or 96-well plates/strips (white)
- Roche LightCycler® Cooling Block
- Sterile reaction tubes
- Calibrated precision pipettes with variable volumes from 2 to 1000 µl

- Sterile single-use pipette tips with aerosol barrier
- (Table) centrifuge (without specific requirements)
- Vortex mixer (without specific requirements)
- Water, PCR grade
- EDTA plasma tubes for sample collection

6. TRANSPORT, STORAGE AND STABILITY

The MutaREAL® DPD-2 real-time PCR kit is transported frozen on dry ice. All components are to be stored protected from light at a minimum of -20°C immediately after receipt. Avoid more than five freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

After opening, the reaction mixes can be stored at +2 to +8°C for up to 6 months.

Be sure to protect the reaction mixes from direct sunlight during the entire test procedure.

7. WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.

- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8. SAMPLE MATERIAL

Starting material for the MutaREAL® DPD-2 real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9. REAL-TIME-PCR

9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.
- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the reaction mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, four reaction tubes per sample and two additional reaction tubes per master mix are required for the negative and the positive control. The following table shows the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mix is pipetted as described in tables 3, 4, 5 and 6:

Master mix 1 (*2A)

Table 3: Preparation of master mix 1 (*2A)

Reagent	Volume per 25 µl - reaction mix	Master mix volume
Reaction mix 1 (yellow)	23 µl	23 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 1 (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Master mix 2 (*Haplo B3)

Table 4: Preparation of master mix 2 (*Haplo B3)

Reagenz	Volume per 25 µl - reaction mix	Master mix volume
Reaction mix 2 (white)	23 µl	23 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 2 (**orange**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Master mix 3 (*D949V)

Table 5: Preparation of master mix 3 (*D949V)

Reagenz	Volume per 25 µl - reaction mix	Master mix volume
Reaction mix 3 (black)	23 µl	23 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 3 (**purple**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

Master mix 4 (*13)

Table 6: Preparation of master mix 3 (*13)

Reagenz	Volume per 25 µl - reaction mix	Master mix volume
Reaction mix 4 (black)	23 µl	23 µl * (N + (N * 0.1))

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each reaction tube.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control 4 (**brown**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding reaction tube.

LightCycler® 1.5 / 2.0: Close the capillaries with the lids, transfer them into the LightCycler® carousel and centrifuge them in the LightCycler® centrifuge (if a table-top centrifuge is used, centrifuge the capillaries in the inserts of the cooling block at 3000 rpm for 15 s). Then transfer the carousel to the LightCycler® device and start the PCR program described in 9.3.

LightCycler® 480 II / PRO: Seal the wells with a sealing foil and centrifuge the plate (2000 rpm for 15 s). Then transfer the plate into the LC 480 II / PRO and start the PCR program described in 9.3.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 7.

Table 7: Real-Time-PCR thermal profile

Description	Time	Temperature	Cycles	Acquisition
Initial Denaturation	2 min	95 °C	1	none
Denaturation	10 s	94 °C	50	none
Primer annealing	25 s	56 °C		single
Elongation	25 s	72 °C		none
Melting curve*/**	30 s	95 °C	1	none
	30 s	60 °C	1	none
	30 s	50 °C	1	none
	30 s	40 °C	1	none
	0 s	80 °C	1	constant
Cooling	30 s	40 °C	1	-

* For the Roche LightCycler® 1.5/2.0, select a ramp rate of 0.05 °C/s.

** For the Roche LightCycler® 480 II, 6 acquisitions (per °C) must be selected.

10. DATA ANALYSIS

Add an analysis of the type „genotyping“ for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is to be selected at 660 / 705 nm (F3).

10.1 DPD*2A

LightCycler® 1.5 / 2.0

Temperature mutation allele*: 50.0 °C (+/-2 °C)

Temperature wild type allele*: 60.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **black curve** - negative control, **green curve** - homozygous mutated, **blue curve** - heterozygous mutant, **red curve** - homozygous wild type.


Fig. 1: Evaluation of master mix 1 - LightCycler® 2.0.

LightCycler® 480 II / PRO

Temperature mutation allele*: 51.0 °C (+/-2 °C)

Temperature wild type allele*: 60.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **red curves** - homozygous wild type or homozygous mutant, **green curve** - heterozygous mutant.


Fig. 2: Evaluation of master mix 1 - LightCycler® 480 II.

10.2 DPD*Haplotyp B3

LightCycler® 1.5 / 2.0

Temperature mutation allele*: 61.0 °C (+/-2 °C)

Temperature wild type allele*: 63.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **green curve** - negative control, **pink curve** - heterozygous mutated, **blue curve** homozygous wild type.


Fig. 3: Evaluation of master mix 2 - LightCycler® 2.0.

LightCycler® 480 II / PRO

Temperature mutation allele*: 61.0 °C (+/-2 °C)

Temperature wild type allele*: 63.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **red curve** - homozygous wild type or heterozygous mutant.


Fig. 4: Evaluation of master mix 2 - LightCycler® 480 II.

10.3 DPD*D949V

LightCycler® 1.5 / 2.0

Temperature mutation allele*: 56.0 °C (+/-2 °C)

Temperature wild type allele*: 60.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **green curve** - negative control, **yellow curve** - homozygous wild type, **pink curve** - heterozygous mutant, **purple curve** - homozygous mutant.


Fig. 5: Evaluation of master mix 3 - LightCycler® 2.0.

LightCycler® 480 II / PRO

Temperature mutation allele*: 56.0 °C (+/-2 °C)

Temperature wild type allele*: 62.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **red curves** - homozygous wild type or homozygous mutant, **green curve** - heterozygous mutant.


Fig. 6: Evaluation of master mix 3 - LightCycler® 480 II.

10.4 DPD*13

LightCycler® 1.5 / 2.0

Temperature mutation allele*: 66.0 °C (+/-2 °C)

Temperature wild type allele*: 57.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **green curve** - negative control, **pink curve** - homozygous wild type, **brown curve** - heterozygous mutant, **grey curve** - homozygous mutant.


Fig. 7: Evaluation of master mix 4 - LightCycler® 2.0.

LightCycler® 480 II / PRO

Temperature mutation allele*: 66.5 °C (+/-2 °C)

Temperature wild type allele*: 58.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **blue curve** - negative control, **red curves** - homozygous wild type or homozygous mutant, **green curve** - heterozygous mutant.


Fig. 8: Evaluation of master mix 4 - LightCycler® 480 II.

** Depending on the PCR device used, the T_m values determined may deviate by more than $\pm 2^\circ\text{C}$. In this case use the positive controls (red / orange / purple and brown) as a reference and apply $\pm 2^\circ\text{C}$ to these T_m values.*

The positive controls 1, 2, 3 & 4 provided contain one DNA each which is heterozygous for the DPD*2A, DPD*Haplotype B3, DPD*D949V or DPD*13 mutation.

11. TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Reaction mixes have been subjected to more than five freeze cycles or have been stored at $2-8^\circ\text{C}$ for more than 6 months. Repeat the analysis with a fresh aliquot or new reaction mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The reaction mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

12. LIMITATIONS OF THE METHOD

The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. When alleles are detected, the polymorphism tested is indicated. Other rare alleles may be present and are not covered by this method. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13. ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleic acid		Catalog number
PCR	Polymerase chain reaction		To be used with
	Reaction mix 1		Contains sufficient for <n> test
	Reaction mix 2		Upper limit of temperature
	Reaction mix 3		Manufacturer
	Reaction mix 4		Lot number
	Positive control 1		Consult instructions for use
	Positive control 2		Content
	Positive control 3		<i>In vitro</i> diagnostic medical device
	Positive control 4		Use by YYYY-MM-DD
	Negative control		

14. LITERATURE

1. Del Re et al., EPMA J. 2010 Sep; 1(3): 495–502

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com