MutaPLEX[®] MRSA real time PCR kit

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MutaPLEX® MRSA real time PCR kit

*Test for the qualitative in vitro detection of MRSA DNA
in clinical specimens*

Valid from 2017-11-10

REF **KG190332**

Σ 32

REF **KG190396**

Σ 96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1	INTENDED USE	18
2	PATHOGEN INFORMATION	18
3	PRINCIPLE OF THE TEST	18
4	PACKAGE CONTENTS	19
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	19
6	TRANSPORT, STORAGE AND STABILITY	19
7	IMPORTANT NOTES	20
8	GENERAL PRECAUTIONS	20
9	SAMPLE MATERIAL	20
10	SAMPLE PREPARATION	20
11	CONTROL DNA	21
	<i>DNA isolation from clinical samples</i>	<i>21</i>
12	REAL TIME PCR	21
	<i>12.1 Important points before starting</i>	<i>21</i>
	<i>12.2 Procedure</i>	<i>22</i>
	<i>12.3 Instrument settings</i>	<i>23</i>
13	DATA ANALYSIS	24
14	ASSAY VALIDATION	26
15	LIMITATIONS OF THE METHOD	27
16	TROUBLESHOOTING	27
17	KIT PERFORMANCE	29
	<i>17.1 Diagnostic Sensitivity and Specificity</i>	<i>29</i>
	<i>17.2 Analytical Sensitivity</i>	<i>29</i>
	<i>17.3 Analytical Specificity</i>	<i>29</i>
18	ABBREVIATIONS AND SYMBOLS	30
19.	LITERATURE	31

1 INTENDED USE

The MutaPLEX® MRSA real time PCR™ is an assay for the detection of MRSA DNA in clinical specimens in open real time PCR machines.

2 PATHOGEN INFORMATION

Staphylococcus aureus are gram-positive coccal bacteria which are ubiquitously found in the environment. About 25–30% of the human population are long-term carriers of *Staphylococcus aureus* because the bacteria are frequently part of the flora found in the nose and on skin. *Staphylococcus aureus* can cause a range of illnesses such as minor skin infections, like furuncles and abscesses, pyomyositis, but also life-threatening diseases such as pneumonia, endocarditis, toxic shock syndrome (TSS), and sepsis.

Of increasing importance worldwide are methicillin-resistant *Staphylococcus aureus* (MRSA) strains. Especially in hospitals, MRSA present a danger, because they are resistant to all β -lactam antibiotics (e.g. penicillin) and often possess further resistances to other antibiotics.

3 PRINCIPLE OF THE TEST

The MutaPLEX® MRSA real time PCR kit contains specific primers and dual-labelled probes for the amplification and detection of MRSA DNA in clinical specimens. The PCR targets the *orfX/SSCmec* junction and allows for the detection of MRSA in clinical samples, even those containing coagulase-negative *Staphylococci*. Furthermore, MutaPLEX® MRSA real time PCR kit allows the detection of the methicillin resistance gene *mecA/mecC*, to eliminate false positive results through dropout mutants.

The presence of nucleic acid is detected by an increase in fluorescence due to hydrolysis of the probes during amplification. The fluorescence of the pathogen-specific probes is measured in the FAM channel. The fluorescence of the *mecA/mecC* gene-specific probes is measured in the Cy5 channel. For a positive MRSA result, both channels need to show an amplification.

Furthermore, MutaPLEX® MRSA real time PCR kit contains a control DNA, which is added during DNA extraction and detected in the same reaction by a differently labelled probe.

The control DNA allows the detection of PCR inhibition and acts as control for the isolation of the nucleic acid from the clinical specimen. The fluorescence of the control DNA is measured in the VIC®/HEX/JOE/TET channel.

4 PACKAGE CONTENTS

The reagents supplied are sufficient for 32 (KG190332) or 96 (KG190396) reactions, respectively.

Table 1: Components of the MutaPLEX® MRSA real time PCR kit.

Label	Lid Colour	Content	
		32	96
Reaction Mix	yellow	1 x 512 µl	2 x 768 µl
Positive control	red	1 x 50 µl	1 x 100 µl
Negative control	green	1 x 50 µl	1 x 100 µl
Control DNA	colourless	1 x 160 µl	2 x 240 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA isolation kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038, MutaCLEAN® Plus lysis buffer, KG1036, or the magnet particle based system MutaCLEAN® Mag RNA/DNA, KG1023)
- PCR grade water
- Sterile microtubes
- Pipets (adjustable volume)
- Sterile pipet tips with filter
- Table centrifuge
- Vortex mixer
- Real time PCR instrument
- Optical PCR reaction tubes with lid
- Optional: Liquid handling system for automation
- Optional: BLP-DNA (bacterium-like particles, KG7013, see chapter 11)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

6 TRANSPORT, STORAGE AND STABILITY

The MutaPLEX® MRSA real time PCR-Kit is shipped on dry ice. All components must be stored at maximum -20 °C in the dark immediately after receipt. Do not use reagents after the date of expiry printed on the package.

Up to 20 freeze and thaw cycles are possible.

For convenience, opened reagents can be stored at 2–8 °C for up to 6 months.

Protect kit components from direct sunlight during the complete test run.

7 IMPORTANT NOTES

- The MutaPLEX® MRSA real time PCR must be performed by qualified personnel only.
- Good Laboratory Practice (GLP) has to be applied.
- Clinical samples must always be regarded as potentially infectious material and all equipment used has to be treated as potentially contaminated.

8 GENERAL PRECAUTIONS

- Stick to the protocol described in the instructions for use.
- Set up different laboratory areas for the preparation of samples and for the set up of the PCR in order to avoid contaminations.
- Pipettes, tubes and other materials must not circulate between those different laboratory areas.
- Always use filter tips.
- Regularly decontaminate equipment and benches with ethanol-free decontaminant.
- Do not combine MutaPLEX® MRSA real time PCR-Kit components of different lot numbers.

9 SAMPLE MATERIAL

Starting material for the MutaPLEX® MRSA real time PCR is bacterial DNA isolated or released from clinical specimens (e.g. throat swabs, nasal swabs, bronchial lavage, wound swabs).

10 SAMPLE PREPARATION

The MutaPLEX® MRSA real time PCR is suitable for the detection of MRSA DNA isolated from clinical specimens with appropriate isolation methods.

Commercial kits for DNA isolation such as MutaCLEAN® Universal RNA/DNA (KG1038) are recommended.

Important: In addition to the samples, always run a water control in your extraction. Treat this water control analogous to a sample.

Comparing the amplification of the control DNA in the samples to the amplification of the internal control in the water control will give insights on possible inhibitions of the real time PCR. Furthermore, possible contaminations during DNA extraction will be detectable.

Please note chapter 11 “Control DNA”.

If the real time PCR is not performed immediately, store extracted DNA according to the instructions given by the DNA extraction kit’s manufacturer.

11 CONTROL DNA

A control DNA is supplied and can be used as extraction control or only as inhibition control. This allows the user to control the DNA isolation procedure and to check for possible real time PCR inhibition.

The bacterium-like particles (BLP-DNA, KG7013) are not supplied.

DNA isolation from clinical samples

a) Control DNA or BLP-DNA used as extraction control

MutaPLEX® MRSA control DNA or BLP-DNA is added to the DNA extraction.

Add 5 µl control DNA or BLP-DNA per extraction (5 µl x (N+1)). Mix well. Perform the DNA isolation according to the manufacturer’s instructions. Please follow protocol A.

The control DNA must be added to the lysis buffer of the extraction kit.

b) Control DNA used as internal control of the real time PCR

If only inhibition will be checked, please follow protocol B.

12 REAL TIME PCR

12.1 Important points before starting

- Please pay attention to chapter 7 “Important Notes”.
- Before setting up the real time PCR familiarise yourself with the real time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run, one positive control and one negative control should be included.

- Before each use, all reagents should be thawed completely at room temperature, thoroughly mixed, and centrifuged very briefly.
- We recommend to keep reagents and samples at 2–8°C (e.g. on ice or a cooling block) at all times.

12.2 Procedure

If the control DNA or BLP-DNA is used to control both, the real time PCR and the DNA isolation procedure, please follow protocol A. If the control DNA is solely used to detect possible inhibition of the real time PCR, please follow protocol B.

Protocol A

The control DNA or BLP-DNA was added during DNA extraction (see chapter 11 “Control DNA”). In this case, prepare the master mix according to table 2.

The master mix contains all of the components needed for PCR except the sample. Prepare a volume of master mix for at least one sample more than required, in order to compensate for pipetting inaccuracy.

Table 2: Preparation of the master mix (control DNA was added during DNA extraction)

Volume per reaction	Volume master mix
16 µl Reaction Mix	16 µl x (N+1)

Protocol B

The control DNA is used for the control of the real time PCR only (see chapter 11 “Control DNA”). In this case, prepare the master mix according to table 3.

The master mix contains all of the components needed for real PCR except the sample. Prepare a volume of master mix for at least one sample more than required, in order to compensate for pipetting inaccuracy.

Table 3: Preparation of the master mix (control DNA is added directly to the master mix)

Volume per reaction	Volume master mix
16 µl Reaction Mix	16 µl x (N+1)
0.5 µl Control DNA *	0.5 µl x (N+1)*

*The increase in volume caused by adding the control DNA is not taken into account when preparing the PCR assay. The sensitivity of the detection system is not impaired.

Protocol A and B: real time PCR set up

- Place the number of optical PCR reaction tubes needed into the respective tray of the real time PCR instrument.

- Pipet 16 µl of the master mix into each optical PCR reaction tube.
- Add 4 µl of the eluates from the DNA isolation (including the eluate of the water control), the positive control and the negative control to the corresponding optical PCR reaction tube (table 4).
- Close the optical PCR reaction tubes immediately after filling in order to reduce the risk of contamination.

Table 4: Preparation of the real time PCR

Component	Volume
Master mix	16.0 µl
Sample	4.0 µl
Total volume	20.0 µl

12.3 Instrument settings

For the real time PCR use the thermal profile shown in table 5.

Table 5: real time PCR thermal profile

Description	Time	Temperature	No of cycles
Reverse transcription	10 min	45 °C	1
Initial Denaturation	5 min	95 °C	1
Amplification of DNA			45
Denaturation	10 s	95 °C	
Annealing and extension	40 s	60 °C	
	Aquisition at the end of this step		

Dependent on the real time instrument used, further instrument settings have to be adjusted according to table 6.

Table 6: Overview of the instrument settings required for the MutaPLEX® MRSA real time PCR.

Real time PCR Instrument	Parameter	Detection Channel	Notes		
LightCycler 480II	MRSA SCCmec Control DNA mecA/mecC Mutation/ Deletion	FAM (465-510) HEX (533-580) CY5 (618-660)	Colour compensation kit MutaPLEX® CC-1 (KG19-5-CC) not required		
			Melt factor	Quant factor	Max integration time (s)
			1	10	1
			1	10	2
Stratagene Mx3000P/ Mx3005P	MRSA SCCmec Control DNA mecA/mecC Mutation/ Deletion	FAM HEX Cy5	Gain 8	Reference dye: none	
			Gain 1		
			Gain 4		
ABI 7500	MRSA SCCmec Control DNA mecA/mecC Mutation/ Deletion	FAM JOE Cy5	Option reference dye ROX: NO		
Rotor-Gene Q, Rotor-Gene 3000 Rotor-Gene 6000	MRSA SCCmec Control DNA mecA/mecC Mutation/ Deletion	Green Yellow Red	Gain 5		
			Gain 5		
			Gain 5		

13 DATA ANALYSIS

The MRSA-specific amplification is measured in the FAM channel. The fluorescence of the mecA/mecC-Gene specific probes is measured in the Cy5 channel. The amplification of the control DNA is measured in the VIC®/HEX/JOE/TET channel.

The following results can occur:

Ct values			Interpretation
FAM channel SCCmec	Cy5 channel resistance gene mecA/mecC	HEX channel	
pos	pos	pos or neg	Positive result, the sample contains MRSA DNA. The result for the control DNA is irrelevant.
pos	neg	pos or neg	Negative result, the sample contains MS-MRSA DNA. The result for the control DNA is irrelevant.
neg	pos	pos or neg	Negative result, the sample contains MR-CoNS DNA. The result for the control DNA is irrelevant.
neg	neg	27-35*	Negative result, the sample contains no MRSA/MS-MRSA and MR-CoNS DNA.
neg	neg	> 35*/ neg	No diagnostic statement can be made. The real time PCR is either inhibited or errors occurred while DNA extraction.

*Depending on the PCR instrument and/or the chosen extraction method, the Ct values might be shifted. The water control can be used as reference. If the HEX Ct value of a sample differs a lot from the water control, partial inhibition has occurred, leading to false negative results in case of weak positive samples.

Figure 1 and figure 2 show examples for positive and negative real time PCR results.


Figure 1: The positive sample shows amplification signal in the bacteria specific channel (FAM/Cy5), whereas no fluorescence signal is detected in the negative sample.


Figure 2: The positive sample as well as the negative sample show a signal in the control DNA-specific VIC®/HEX/JOE/TET channel. The amplification signal of the control DNA in the negative sample shows that the missing signal in the MRSA-specific FAM channel is not due to PCR inhibition or failure of DNA isolation, but that the sample is a true negative.

14 ASSAY VALIDATION

Set a threshold as follows:

Negative controls

All negative controls should be below the threshold. If there is a potential contamination (appearance of a curve in the negative control or a cluster of curves in speci-

mens at high CT – for example above 36), results obtained are not interpretable and the whole run (including extraction) has to be repeated.

Positive controls

All the positive controls must show a positive (i. e. exponential) amplification curve. The positive controls must fall below a CT of 30.

Internal controls

All internal controls must show a positive (i. e. exponential) amplification curve. The internal control must fall below a CT of 33. If the internal control is above CT 34, this points to a purification problem or a strong positive sample that can inhibit the IC. In the latter case, the assay is valid. If a water control run is performed, the IC must fall below a CT of 33.

15 LIMITATIONS OF THE METHOD

The results must always be considered in relation to the clinical symptoms. Therapeutical consequences should be made in consideration of clinical data.

A negative test result does not exclude a MRSA infection.

16 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a real time PCR.

No fluorescence signal in the FAM or Cy5 channel of the positive control

The selected channel for analysis does not comply with the protocol

Select the FAM channel for analysis of the MRSA-specific amplification, the Cy5 channel for the *mecA/mecC*-specific amplification and the VIC®/HEX/JOE/TET channel for the amplification of the control DNA .

Incorrect configuration of the real time PCR

Check your work steps and compare with chapter “Procedure”.

The programming of the thermal profile is incorrect

Compare the thermal profile with the protocol (table 5).

Incorrect storage conditions for one or more kit components or kit expired

Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in chapter “Transport, Storage and Stability”.

Weak or no signal of the control DNA and simultaneous absence of a signal in the MRSA-specific FAM or Cy5 channel***real time PCR conditions do not comply with the protocol***

Check the real time PCR conditions (chapter 12).

real time PCR inhibited

Make sure that you use an appropriate isolation method (see “Sample preparation”) and follow the manufacturer’s instructions. Make sure that the ethanol-containing washing buffer of the isolation kit has been completely removed. An additional centrifugation step at high speed is recommended before elution of the DNA.

DNA loss during isolation process

In case the control DNA was added before extraction, the lack of an amplification signal can indicate that the DNA isolation was not successful. Make sure that you use an appropriate isolation method (commercial kits are recommended) and stick to the manufacturer’s protocol.

Incorrect storage conditions for one or more components or kit expired

Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in chapter “Transport, Storage and Stability”.

Detection of a fluorescence signal in the FAM or Cy5 channel of the negative control***Contamination during preparation of the PCR***

Repeat the real time PCR in replicates. If the result is negative in the repetition, the contamination occurred when the samples were pipetted into the optical PCR reaction tubes. Make sure to pipet the positive control last and close the optical PCR reaction tube immediately after adding the sample. If the same result occurs, one or more of the kit components might be contaminated. Make sure that work space and instruments are decontaminated regularly. Use a new kit and repeat the real time PCR.

17 KIT PERFORMANCE

17.1 Diagnostic Sensitivity and Specificity

During the validation study of the MutaPLEX® MRSA real time PCR kit, 14 positive and 13 negative samples were tested. The diagnostic sensitivity was found to be 100% and the diagnostic specificity 100% (table 7).

The positive predictive value was found to be 100%, the negative predictive value showed to be 100%.

Table 7: Overview of the amount of samples tested and the resulting positive and negative predictive values

	positive samples	negative samples
MutaPLEX® MRSA positive	14	0
MutaPLEX® MRSA negative	0	13
Sensitivity	100%	
Specificity	100%	

17.2 Analytical Sensitivity

The limit of detection (LoD) of MutaPLEX® MRSA real time PCR kit was determined using serial dilutions of synthetic target DNA sequences in a Stratagene Mx3000 real time PCR instrument. The LoD of MutaPLEX® MRSA real time PCR kit is at least 1 cfu per reaction each.

17.3 Analytical Specificity

The specificity of the MutaPLEX® MRSA real time PCR kit was evaluated additionally with different other relevant viruses and bacteria found in clinical samples (table 8).


The primers and probes detected the MRSA sample but not any other pathogen tested.


Table 8: Bacterial and viral pathogens used for the determination of the specificity of primers and probes of the MutaPLEX® MRSA real time PCR kit.

Strain	Expected result	Result
<i>Streptococcus agalactiae</i>	negative	negative
<i>Coxsackievirus</i> Strain P.B.	negative	negative
<i>Coxsackievirus</i> Strain B.S.	negative	negative
<i>Herpes simplex virus</i>	negative	negative

Strain	Expected result	Result
<i>Borrelia burgdorferi</i>	negative	negative
Tick borne encephalitis	negative	negative
Influenza A	negative	negative
Influenza B	negative	negative
Respiratory syncytial virus A	negative	negative
Respiratory syncytial virus B	negative	negative
<i>Legionella pneumophila</i> serogroup 1	negative	negative
Cytomegalovirus	negative	negative
MRSA	positive	positive

18 ABBREVIATIONS AND SYMBOLS

DNA	Deoxyribonucleid Acid	BLP	Bacterium-Like Particles
CT	Cycle threshold		Control DNA
MRSA	Methicillin-resistant <i>Staphylococcus aureus</i>		To be used with
MS-MRSA	Methicillin-suceptible MRSA, <i>mecA</i> dropout mutant	OrfX/SCCmec	Junction for <i>S. aureus</i> DNA and SCCmec cassette
MSSA	Methicillin-suceptible <i>Staphylococcus aureus</i>	MR-ConS	Methicillin-resistant coagulase negative <i>Staphylococcus</i>
<i>mecA</i> / <i>mecC</i>	Two variants of the methicillin resistance gene		Manufacturer
PCR	Polymerase Chain Reaction		Catalog number
	Negative control		Contains sufficient for <n> test

	Positive control		Upper limit of temperature
	Reaction Mix		Use by
	Content		Lot number
	Consult instructions for use		In vitro diagnostic medical device

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