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MutaPLEX[®] MRSA seqc real time PCR kit

*Test for the qualitative in vitro detection of MRSA DNA
in biological specimens*

Valid from 2021-10-15

REF **KG194196**

IVD

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1	INTENDED USE	24
2	PATHOGEN INFORMATION	24
3	PRINCIPLE OF THE TEST	24
4	PACKAGE CONTENTS	25
5	EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER	25
6	TRANSPORT, STORAGE AND STABILITY	25
7	WARNINGS AND PRECAUTIONS	26
8	SAMPLE MATERIAL	27
9	SAMPLE PREPARATION	27
10	CONTROL DNA	27
11	REAL TIME PCR	27
	11.1 <i>Important points before starting</i>	27
	11.2 <i>Procedure</i>	28
	11.3 <i>Instrument settings</i>	29
12	DATA ANALYSIS	31
13	ASSAY VALIDATION	33
14	LIMITATIONS OF THE METHOD	34
15	TROUBLESHOOTING	35
16	KIT PERFORMANCE	36
	16.1 <i>Analytical Sensitivity</i>	36
	16.2 <i>Analytical Specificity</i>	37
	16.3 <i>Linear Range</i>	40
	16.4 <i>Precision</i>	41
	16.5 <i>Diagnostic Sensitivity</i>	42
17	ABBREVIATIONS AND SYMBOLS	43
18.	LITERATURE	44

1 INTENDED USE

The MutaPLEX® MRSA seqc real time PCR is an assay for the detection of DNA of MRSA extracted from biological specimens.

2 PATHOGEN INFORMATION

Staphylococcus aureus are gram-positive coccal bacteria which are ubiquitously found in the environment. About 25–30% of the human population are long-term carriers of *Staphylococcus aureus* because the bacteria are frequently part of the flora found in the nose and on skin. *Staphylococcus aureus* can cause a range of illnesses such as minor skin infections, like furuncles and abscesses, pyomyositis, but also life-threatening diseases such as pneumonia, endocarditis, toxic shock syndrome (TSS), and sepsis.

Of increasing importance worldwide are methicillin-resistant *Staphylococcus aureus* (MRSA) strains. Especially in hospitals, MRSA present a danger, because they are resistant to all β -lactam antibiotics (e.g. penicillin) and often possess further resistances to other antibiotics.

3 PRINCIPLE OF THE TEST

The MutaPLEX® MRSA seqc real time PCR kit contains specific primers and dual-labelled probes for the amplification and detection of MRSA DNA extracted from biological specimen. The PCR targets the SCCmec/orfX junction and allows for the detection of MRSA in biological samples, even those containing coagulase-negative *Staphylococci*. Furthermore, MutaPLEX® MRSA seqc real time PCR kit allows the detection of the methicillin resistance gene *mecA/mecC*, to eliminate false positive results through dropout mutants.

The presence of nucleic acid is detected by an increase in fluorescence due to hydrolysis of the probes during amplification. The fluorescence of the pathogen-specific probes is measured in the FAM channel. The fluorescence of the *mecA/mecC* gene-specific probes is measured in the Cy5 channel. For a positive MRSA result, both channels need to show an amplification.

Furthermore, MutaPLEX® MRSA seqc real time PCR kit contains a control DNA (internal process control, IPC), which is added during DNA extraction and detected in the same reaction by a HEX-labelled probe.

The Control DNA allows the detection of PCR inhibition and acts as control, that the nucleic acid was isolated from the biological specimen.

Additionally, MutaPLEX® MRSA seqc real time PCR Kit contains an Internal System Control (ISC). The ISC consists of primers and probes for the detection of a house

keeping gene (Succinate dehydrogenase) in the eluate from biological specimens. The ISC helps preventing false negative results due to insufficient sample drawing or transport. The amplification of the succinate dehydrogenase target sequence is measured in the ROX channel.

4 PACKAGE CONTENTS

The reagents supplied are sufficient for 96 reactions.

Table 1: Components of the MutaPLEX® MRSA seqc real time PCR kit .

Label	Lid Colour	Content
		96
Reaction Mix	yellow	1 x 1344 µl
Positive control	red	1 x 150 µl
Negative control	green	1 x 150 µl
Control DNA	colourless	1 x 480 µl

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA isolation kit (e.g. MutaCLEAN® Mag RNA/DNA, KG1023/KG1024)
- PCR grade water
- Sterile microtubes
- Pipets (adjustable volume)
- Sterile pipet tips with filter
- Table centrifuge
- Vortex mixer
- Real time PCR instrument
- Optical PCR reaction tubes with lid or optical PCR reaction plate with optical foil
- Optional: Liquid handling system for automation

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

6 TRANSPORT, STORAGE AND STABILITY

The MutaPLEX® MRSA seqc real time PCR-Kit is shipped on dry ice. All components must be stored at maximum -20 °C in the dark immediately after receipt. Do not use reagents after the date of expiry printed on the package.

Up to 20 freeze and thaw cycles are possible.

For convenience, opened reagents can be stored at 2–8 °C for up to 6 months.

Protect kit components from direct sunlight during the complete test run.

7 WARNINGS AND PRECAUTIONS

- Read the Instructions for Use carefully before using the product.
- Before first use check the product and its components for:
- Use of this product is limited to personnel specially instructed and trained in the techniques of real time PCR procedures.
- Specimens should always be treated as infectious and/or biohazardous in accordance with safe laboratory procedures.
- Avoid microbial and nuclease (DNase/RNase) contamination of the eluates and the components of the kit.
- Always use DNase/RNase-free disposable pipette tips with aerosol barriers.
- Always wear protective disposable powder-free gloves when handling kit components.
- Use separated and segregated working areas for (1) sample preparation, (2) reaction setup and (3) amplification/detection activities. The workflow in the laboratory should proceed in unidirectional manner.
- Always wear disposable gloves in each area and change them before entering a different area.
- Dedicate supplies and equipment to the separate working areas and do not move them from one area to another.
- Store positive and/or potentially positive material separated from all other components of the kit.
- Do not open the reaction tubes/plates post amplification to avoid contamination with amplicons.
- Additional controls may be tested according to guidelines or requirements of local, state and/or federal regulations or accrediting organizations.
- Do not autoclave reaction tubes after the PCR, since this will not degrade the amplified nucleic acid and will bear the risk to contaminate the laboratory area.

- Discard sample and assay waste according to your local safety regulations.
- Do not combine kit components of different lot numbers.

8 SAMPLE MATERIAL

Starting material for the MutaPLEX® MRSA seqc real time PCR is bacterial DNA isolated from biological specimens.

9 SAMPLE PREPARATION

The MutaPLEX® MRSA seqc real time PCR is suitable for the detection of MRSA DNA isolated from biological specimens with appropriate isolation methods.

Commercial kits for DNA isolation such as MutaCLEAN® Mag RNA/DNA (KG1023/KG1024) are recommended.

Important: In addition to the samples, always run a water control in your extraction. Treat this water control analogous to a sample.

Comparing the amplification of the control DNA in the samples to the amplification of the internal control in the water control will give insights on possible inhibitions of the real time PCR. Furthermore, possible contaminations during DNA extraction will be detectable.

Please note chapter 10 “Control DNA”.

If the real time PCR is not performed immediately, store extracted DNA according to the instructions given by the DNA extraction kit's manufacturer.

10 CONTROL DNA

A control DNA is supplied and can be used as extraction control. This allows the user to control the DNA isolation procedure and to check for possible real time PCR inhibition.

Add 5 µl control DNA per extraction (5 µl x (N+1)). Mix well. Perform the DNA isolation according to the manufacturer's instructions.

The control DNA must be added to the lysis buffer of the extraction kit.

11 REAL TIME PCR

11.1 *Important points before starting*

- Please pay attention to chapter 7 “Important Notes”.

- Before setting up the real time PCR familiarise yourself with the real time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run, one positive control and one negative control should be included.
- Before each use, all reagents should be thawed completely at room temperature, thoroughly mixed, and centrifuged very briefly.

11.2 Procedure

The master mix contains all of the components needed for PCR except the sample. Prepare a volume of master mix for at least one sample more than required, in order to compensate for pipetting inaccuracy.

Table 2: Preparation of the master mix.

Volume per reaction	Volume master mix
14 µl Reaction Mix	14 µl x (N+1)

Real time PCR set up

- Place the number of optical PCR reaction tubes needed into the respective tray of the real time PCR instrument / take an optical PCR reaction plate.
- Pipet 14 µl of the master mix into each optical PCR reaction tube / the optical PCR reaction plate.
- Add 6 µl of the eluates from the DNA isolation (including the eluate of the water control), the positive control and the negative control to the corresponding optical PCR reaction tube / the optical PCR reaction plate (table 3).
- Close the optical PCR reaction tubes / the optical PCR reaction plate immediately after filling in order to reduce the risk of contamination.

Table 3: Preparation of the real time PCR

Component	Volume
Master mix	14.0 µl
Sample	6.0 µl
Total volume	20.0 µl

11.3 Instrument settings

For the real time PCR use one of the thermal profiles shown in table 4 and table 5.

Table 4: real time PCR thermal profile

Description	Time	Temperature	No of cycles
Initial Denaturation	5 min	95 °C	1
Amplification of DNA			45
Denaturation	10 s	95 °C	
Annealing and extension	40 s	60 °C	
	Aquisition at the end of this step		

If in the same run samples should be tested for pathogens with RNA genome, use the thermal profile shown in Table 5.

Table 5: real time PCR thermal profile

Description	Time	Temperature	No of cycles
Reverse transcription	10 min	45 °C	1
Initial Denaturation	5 min	95 °C	1
Amplification of DNA			45
Denaturation	10 s	95 °C	
Annealing and extension	40 s	60 °C	
	Aquisition at the end of this step		

Dependent on the real time instrument used, further instrument settings have to be adjusted according to table 6.

Table 6: Overview of the instrument settings required for the MutaPLEX® MRSA real time PCR.

Real time PCR instrument	Parameter	Detection Channel	Notes		
LightCycler 480II	SCCmec/orfX Control DNA (IPC) ISC mecA/mecC Mutation/Deletion	465-510 533-580 533-610 618-660	Colour compensation Kit MutaPLEX® CC-1 (KG19-5-CC) required		
			Melt factor	Quant factor	Max integration time (s)
			1	10	1
			1	10	2
			1	10	3
Stratagene Mx3000P/ Mx3005P	SCCmec/orfX Control DNA (IPC) ISC mecA/mecC Mutation/Deletion	FAM HEX ROX Cy5	Gain 8 Gain 1 Gain 1 Gain 4	Reference dye: none	
QuantStudio 5 Bio-Rad CFX96 Bio-Rad CFX Opus 96	SCCmec/orfX Control DNA (IPC) ISC mecA/mecC Mutation/Deletion	FAM HEX ROX Cy5	Option reference dye ROX: NO		
Mic qPCR Cycler	SCCmec/orfX Control DNA (IPC) ISC mecA/mecC Mutation/Deletion	Green Yellow Red Orange	Gain 8 Gain 10 Gain 10 Gain 10		

12 DATA ANALYSIS

The following results can occur:

Signal/ C _T Values				Interpretation
FAM Channel SCCmec/orfX	Cy5 Channel resistance gene mecA/mecC	ROX Channel ISC	HEX Channel IPC	
positive	positive	positive or negative ¹	positive or negative ²	Positive result, the sample contains MRSA DNA.
positive	negative	positive or negative ¹	positive or negative ²	Negative result, the sample contains MS-MRSA DNA.
negative	positive	positive or negative ¹	positive or negative ²	Negative result, the sample contains MR-CoNS DNA.
negative	negative	positive	≤ 34 ³	Negative result, the sample contains no MRSA/ MS-MRSA and MR-CoNS DNA.
negative	negative	negative	negative oder > 34 ³	No diagnostic statement can be made. The real time PCR is either inhibited or errors occurred while DNA extraction.

¹ If the analysed samples originate from cultivation, the ROX channel may be negative.

² A strong positive signal in the FAM or the Cy5 can inhibit the IPC. In such cases the result for the Control DNA can be neglected.

³ In case of high C_T values, the IPC should be compared to the water extraction control as described in the chapter 'Assay Validation'.

Figures 1 - 4 show examples for positive and negative real time PCR results.


Fig. 1: The Positive Sample shows pathogen specific amplification in the FAM channel (Positive Sample and Positive Control), whereas no fluorescence signal is detected in the Negative Sample or the Negative Control (Mx3005P).


Fig. 2: The Positive Sample shows pathogen specific amplification in the Cy5 channel (Positive Sample and Positive Control), whereas no fluorescence signal is detected in the Negative Sample or the Negative Control (Mx3005P).


Fig. 3: Signals of the amplification of the ISC in the ROX channel. The amplification signal in the positive and the negative sample indicates a sufficient amount of DNA in the sample eluate and confirms the integrity of the sampling (Mx3005P).


Fig. 4: The Positive Sample and the Negative Sample show a signal in the Control DNA specific HEX channel (IPC). The amplification signal of the Control DNA in the Negative Sample shows that the missing signals in the pathogen specific channels FAM and Cy5 are not due to PCR inhibition or failure of DNA isolation, but that the sample is a true Negative Sample (Mx3005P).

13 ASSAY VALIDATION

Negative controls

The Negative Control must show no C_T in the FAM, Cy5, ROX and HEX channel.

Positive controls

The Positive Control must show a positive (i.e. exponential) amplification curve in the different channels FAM, Cy5 and ROX. The Positive Control must fall below C_T 30.

Internal controls

The following values for the amplification of the internal controls are valid using the nucleic acid extraction kit MutaCLEAN® Mag RNA/DNA. All internal controls (ISC and IPC, seqc – sample and extraction quality control) must show a positive (i.e. exponential) amplification curve.

The Control DNA (IPC) must fall below a C_T of 34. If the Control DNA is above C_T 34 this points to a purification problem or a strong positive sample that can inhibit the IPC. In the latter case, the assay is valid. It is recommended to perform the extraction of a water control in each run. The IPC in the water control must fall below a C_T of 34.

For accurately drawn respiratory swab samples, the ISC shows C_T values from app. 15 to app. 28. A heavily delayed signal of higher than a C_T of 34 indicates a low sample amount. Therefore, false negative results cannot be ruled out. In case of no amplifications neither in the FAM nor in the Cy5 channel, there must be an amplification curve in the ROX channel (ISC) and the HEX (IPC) channel when using eluates of primary samples from humans.

If other nucleic acid extraction kits are used, the customer must define own cutoffs. In this case the C_T value of the Control DNA (IPC) in an eluate from a sample should not be delayed for more than 4 C_T in comparison to an eluate from an extracted water control.

14 LIMITATIONS OF THE METHOD

- Strict compliance with the Instruction for Use is required for optimal results.
- Use of this product is limited to personnel specially instructed and trained in the techniques of real time PCR and in vitro diagnostic procedures.
- Good laboratory practice is essential for proper performance of this assay.
- All reagents should be closely monitored for impurity and contamination. Any suspicious reagents should be discarded.
- This assay must not be used on a biological specimen directly. Appropriate nucleic acid extraction methods have to be conducted prior to using this assay.
- The presence of PCR inhibitors may cause false negative or invalid results.
- As with any diagnostic test, results of the diarellaMRSA seqc real time PCR Kit need to be interpreted in consideration of all clinical and laboratory findings.

15 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a real time PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No fluorescence signal in the FAM, Cy5, ROX channel of the positive control

The selected channel for analysis does not comply with the protocol

Select the FAM channel for analysis of the MRSA specific amplification, the Cy5 channel for the *mecA/mecC* specific amplification, the ROX channel for the amplification of the ISC and the HEX channel for the amplification of the Control DNA (IPC).

Incorrect configuration of the real time PCR

Check your work steps and compare with chapter "Procedure".

The programming of the thermal profile is incorrect

Compare the thermal profile with the protocol (tables 4 and 5).

Incorrect storage conditions for one or more kit components or kit expired

Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in chapter "Transport, Storage and Stability".

Weak or no signal of the control DNA (IPC) and ISC and simultaneous absence of a signal in the MRSA-specific FAM or Cy5 channel

real time PCR conditions do not comply with the protocol

Check the real time PCR conditions (chapter 11).

real time PCR inhibited

Make sure that you use an appropriate isolation method (see "Sample preparation") and follow the manufacturer's instructions. Make sure that the ethanol-containing washing buffer of the isolation kit has been completely removed. An additional centrifugation step at high speed is recommended before elution of the DNA.

DNA loss during isolation process

In case the control DNA was added before extraction, the lack of an amplification signal can indicate that the DNA isolation was not successful. Make sure that you use an appropriate isolation method (commercial kits are recommended) and stick to the manufacturer's protocol.

Incorrect storage conditions for one or more components or kit expired

Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in chapter "Transport, Storage and Stability".

Detection of a fluorescence signal in the FAM or Cy5 channel of the negative control***Contamination during preparation of the PCR***

Repeat the real time PCR in replicates. If the result is negative in the repetition, the contamination occurred when the samples were pipetted into the optical PCR reaction tubes. Make sure to pipet the positive control last and close the optical PCR reaction tube immediately after adding the sample. If the same result occurs, one or more of the kit components might be contaminated. Make sure that work space and instruments are decontaminated regularly. Use a new kit and repeat the real time PCR.

Detection of a fluorescence signal in the ROX channel of the Negative Control***Contamination with human DNA during preparation of the real time PCR***

As long as the ROX channel shows very high C_T values, the contamination is negligible.

If the FAM and Cy5 channel are negative in the Negative Control, the PCR is still valid for the detection of MRSA.

16 KIT PERFORMANCE**16.1 Analytical Sensitivity**

For the FAM and Cy5 channel the limit of detection (LoD) of MutaPLEX® MRSA seqc real time PCR Kit was determined using serial dilutions of the AcroMetrix MRSA Positive Control, CE-IVD. Since there is no ISC within the AcroMetrix MRSA Positive Control, CE-IVD, the LoD for the ROX channel was determined using serial dilutions of a synthetic DNA-fragment containing the specific gene target sequence. The determination of the LoD was done on a Stratagene Mx3005P.

The LoD of MutaPLEX® MRSA seqc real time PCR Kit is < 2.5 genome copies per µl for the FAM, Cy5 and ROX channel.

16.2 Analytical Specificity

The specificity of the MutaPLEX® MRSA seqc real time PCR was evaluated with different ring trial samples of known status and different other relevant viruses and bacteria found in biological samples and basing on in silico analyses.

All ring trial samples were detected correctly. Results are shown in table 8 and table 9.

The results for the sample analysis are shown table 10 and table 11. The results for the in silico analysis are shown in table 12

Table 8: Ring trial samples tested for the validation of the sensitivity of the MutaPLEX® MRSA seqc real time PCR Kit. Results for the FAM channel.

sample	expected result	result
<i>cMRSA (S. aureus, oxaR, PVL-pos, spa:t 008)</i>	positive	positive
MRSA (<i>S. aureus, oxaR, PVL-neg</i>)	positive	positive
Escherichia coli K12	negative	negative
<i>cMSSA + CoNS (S. aureus, S. epidermidis oxaR, PVL-pos)</i>	negative	negative
<i>CoNS (S. epidermidis, oxaS)</i>	negative	negative
MSSA + CoNS (<i>S. aureus, S. hämolyticus oxaR, PVL-neg</i>)	negative	negative
MRSA SCCmec TypV (<i>S. aureus, oxaR, PVL-neg</i>)	positive	positive
<i>cMRSA (S. aureus, oxaR, PVL-pos)</i>	positive	positive
MRSA (<i>S. aureus, PVL-neg, pSA442 neg</i>)	positive	positive
CoNS <i>oxaS</i>	negative	negative
<i>cMRSA (S. aureus, oxaR, PVL-pos)</i>	positive	positive
<i>cMRSA (S. aureus, oxaR, PVL-pos)</i>	positive	positive

Table 9: Ring trial samples tested for the validation of the sensitivity of the MutaPLEX® MRSA seqc real time PCR Kit. Results for the Cy5 channel.

sample	expected result	result
<i>cMRSA (S. aureus, oxaR, PVL-pos, spa:t 008)</i>	positive	positive
MRSA (<i>S. aureus, oxaR, PVL-neg</i>)	positive	positive

sample	expected result	result
Escherichia coli K12	negative	negative
cMSSA + CoNS (<i>S. aureus</i> , <i>S. epidermidis oxaR</i> , PVL-pos)	positive	positive
CoNS (<i>S. epidermidis</i> , <i>oxaS</i>)	negative	negative
MSSA + CoNS (<i>S. aureus</i> , <i>S. hämolyticus oxaR</i> , PVL-neg)	positive	positive
MRSA SCCmec TypV (<i>S. aureus</i> , <i>oxaR</i> , PVL-neg)	positive	positive
cMRSA (<i>S. aureus</i> , <i>oxaR</i> , PVL-pos)	positive	positive
MRSA (<i>S. aureus</i> , PVL-neg, pSA442 neg)	positive	positive
CoNS <i>oxaS</i>	negative	negative
cMRSA (<i>S. aureus</i> , <i>oxaR</i> , PVL-pos)	positive	positive
cMRSA (<i>S. aureus</i> , <i>oxaR</i> , PVL-pos)	negative	negative

Tabelle 10: Eluted DNA/RNA from bacterial and viral pathogens tested for the determination of the analytical specificity of MutaPLEX® MRSA seqc real time PCR Kit, FAM channel.

sample	expected result FAM channel	MutaPLEX® MRSA seqc FAM channel
<i>Pneumocystis jirovecii</i>	negative	negative
<i>Streptococcus agalactiae</i>	negative	negative
Coxsackie A9 Strain P.B.	negative	negative
Herpes simplex Virus Type 2 Str. G	negative	negative
<i>Borrelia burgdorferi</i>	negative	negative
<i>Staphylococcus ueberis</i>	negative	negative
<i>Streptococcus dysagalactiae</i>	negative	negative
<i>Entrococcus faecalis</i>	negative	negative
<i>Klebsiella</i>	negative	negative
<i>Staphylococcus intermedius</i>	negative	negative
<i>Pseudomonas aeruginosa</i>	negative	negative
<i>Staphylococcus sciuri</i>	negative	negative
<i>Legionella pneumophila</i>	negative	negative
TBE Virus K617	negative	negative

sample	expected result FAM channel	MutaPLEX® MRSA seqc FAM channel
Influenza A Virus (H1N1)	negative	negative
Influenza B Virus	negative	negative
Respiratory Syncytial Virus A	negative	negative
Respiratory Syncytial Virus B	negative	negative
Cytomegalie Virus AD169	negative	negative

Tabelle 11: Eluted DNA/RNA from bacterial and viral pathogens tested for the determination of the analytical specificity of MutaPLEX® MRSA seqc real time PCR Kit, Cy5 channel.

sample	expected result Cy5 channel	MutaPLEX® MRSA seqc Cy5 channel
<i>Pneumocystis jirovecii</i>	negative	negative
<i>Streptococcus agalactiae</i>	negative	negative
Coxsackie A9 Strain P.B.	negative	negative
Herpes simplex Virus Type 2 Str. G	negative	negative
<i>Borrelia burgdorferi</i>	negative	negative
<i>Staphylococcus ueberis</i>	negative	negative
<i>Streptococcus dysagalactiae</i>	negative	negative
<i>Entrococcus faecalis</i>	negative	negative
<i>Klebsiella</i>	negative	negative
<i>Staphylococcus intermedius</i>	negative	negative
<i>Pseudomonas aeruginosa</i>	negative	negative
<i>Staphylococcus sciuri</i>	negative	negative
<i>Legionela pneumophila</i>	negative	negative
TBE Virus K617	negative	negative
Influenza A Virus (H1N1)	negative	negative
Influenza B Virus	negative	negative
Respiratory Syncytial Virus A	negative	negative
Respiratory Syncytial Virus B	negative	negative
Cytomegalie Virus AD169	negative	negative

Tabelle 12: Inclusivity of the MutaPLEX® MRSA seqc real time PCR Kit Primers and Probes (in silico analysis).

9 - 1000 whole genome sequences		Homology	Comment
SCCmec / orfX junction	forward Primer Mix	9 / 9 - 593 / 593	based on different SCCmec cassettes
	reverse Primer	1000 / 1000	only orfX
	Probe	1000 / 1000	only orfX
mec A	forward Primer	789 / 789	
	reverse Primer	789 / 789	
	Probe	789 / 789	
mec C	forward Primer	18 / 18	
	reverse Primer	18 / 18	
	Probe	18 / 18	
Succinate dehydrogebase	forward Primer	1000 / 1000	
	reverse Primer	1000 / 1000	
	Probe	1000 / 1000	

16.3 Linear Range

The linear range of the MutaPLEX® MRSA seqc real time PCR Kit was evaluated by analysing logarithmic dilution series of in vitro transcripts of the target sequences.


Fig. 5: Determination of the linear range of the MutaPLEX® MRSA seqc real time PCR Kit for SCCmec in the FAM channel.


Abb. 6: Determination of the linear range of the MutaPLEX® MRSA seqc real time PCR Kit for mecA/mecC in the Cy5 channel.


Abb. 7: Determination of the linear range of the MutaPLEX® MRSA seqc real time PCR Kit for the ISC in the ROX channel.

16.4 Precision

The precision of the MutaPLEX® MRSA seqc real time PCR Kit was determined as intra-assay variability, inter-assay variability and inter-lot variability.

Variability data are expressed by standard deviation and coefficient of variation. The data are based on quantification analyses of defined concentrations of SCCmec specific synthetic DNA, mecA/mecC specific synthetic DNA, ISC specific synthetic DNA and on the threshold cycle of the Control RNA (IPC). The results are shown in Table 13.

Table 13: Precision of the MutaPLEX® MRSA seqc real time PCR Kit.

SCCmec (FAM)	copies/ µl	Standard Deviation	Coefficient of Variation [%]
Intra-Assay-Variability	250	0.33	0.99

SCCmec (FAM)	copies/ µl	Standard Deviation	Coefficient of Variation [%]
Inter-Assay-Variability	250	0.47	1.39
Inter-Lot-Variability	250	0.05	0.14

mecA/mecC (Cy5)	copies/ µl	Standard Deviation	Coefficient of Variation [%]
Intra-Assay-Variability	25	0.25	0.78
Inter-Assay-Variability	25	0.23	0.72
Inter-Lot-Variability	25	0.17	0.52

ISC (ROX)	copies/ µl	Standard Deviation	Coefficient of Variation [%]
Intra-Assay-Variability	2.5	0.22	0.69
Inter-Assay-Variability	2.5	0.20	0.62
Inter-Lot-Variability	2.5	0.22	0.67

IPC (HEX)	copies/ µl	Standard Deviation	Coefficient of Variation [%]
Intra-Assay-Variability	250	0.34	1.20
Inter-Assay-Variability	250	0.85	2.96
Inter-Lot-Variability	250	0.26	0.91



16.5 Diagnostic Sensitivity

The diagnostic sensitivity of real time (RT-) PCR assays is mainly dependent on the DNA/RNA extraction method used to isolate DNA and RNA from various biological specimens. DNA/RNA extraction reagents are not part of the Immundiagnostik AG real time (RT-) PCR kits. Immundiagnostik AG real time (RT-) PCR kits include an extraction control and guidelines for the validation criteria of the extraction control in each reaction. The extraction control indicates inhibition of the real time (RT-) PCR and/or inefficient nucleic acid extraction. It cannot be used as a calibrator.

Therefore, Immundiagnostik AG guarantees the analytical sensitivities and specificities of the real time (RT-) PCR kits, performed with eluted DNA and RNA from reference materials and ring trial samples and with synthetic nucleic acid fragments. Immundiagnostik AG does not guarantee diagnostic sensitivities. If diagnostic sensitivities are mentioned in manuals of Immundiagnostik AG real time (RT-) PCR kits, the data are strictly correlated to a specific nucleic acid extraction method that has

been used during the validation of the respective kits and cannot be transferred to other extraction methods. It is the responsibility of the user to qualify the extraction methods used for DNA/RNA isolation from biological samples.

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