ID-Vit[®] Pantothersäure

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
ID-Vit[®] Pantothenic acid

Microbiological test kit for the determination of total free pantothenic acid (vitamin B₅) in serum using a Lactobacillus plantarum coated microtitre plate
For use in human and veterinary medicine and in research

Valid from 2022-06-20

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	15
2. INTRODUCTION	15
3. MATERIAL SUPPLIED	16
4. MATERIAL REQUIRED BUT NOT SUPPLIED	16
5. STORAGE AND PREPARATION OF REAGENTS	17
5.1 <i>Water</i>	17
5.2 <i>Preparation of the controls</i>	17
5.3 <i>Preparation of the standard curve</i>	17
5.4 <i>Preparation of the sterile assay medium</i>	18
5.5 <i>Microtiter plate [PLATE]</i>	18
6. SAMPLE STORAGE AND PREPARATION	19
6.1 <i>Sample dilution</i>	19
7. ASSAY PROCEDURE	19
7.1 <i>Principle of the test</i>	19
7.2 <i>Test preparations</i>	19
7.3 <i>Test procedure</i>	20
7.4 <i>Measurement</i>	20
8. EVALUATION OF RESULTS	20
8.1 <i>Calculation</i>	21
8.2 <i>Expected values</i>	21
8.3 <i>Quality control</i>	21
9. PERFORMANCE CHARACTERISTICS	21
9.1 <i>Precision and reproducibility</i>	22
9.2 <i>Recovery</i>	22
10. LIMITATIONS	23
11. PRECAUTIONS	23
12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24
13. REFERENCES	24

1. INTENDED USE

ID-Vit® Pantothenic acid is a microtiter plate test kit based on a microbiological method which measures the total free pantothenic acid content in serum. The test kit contains the standard and all reagents required to perform the test. It is sufficient for 96 determinations including standard curves. An ELISA reader is required for the evaluation of the results. For use in human and veterinary medicine and in research. For *in vitro* diagnostic use only.

2. INTRODUCTION

Pantothenic acid is the reactive thiol function of CoA and ACP

Pantothenic acid (vitamin B₅) is synthesised by most microorganisms and plants from pantoic acid. The vitamin is an integral part of 4'-phosphopantetheine, which is a component of coenzyme A (CoA). CoA plays a key role in the metabolism of numerous compounds, especially lipids and the ultimate catabolic disposition of carbohydrates and ketogenic amino acids. About 80% of the vitamin in animal tissues is in CoA form, and the rest exists mainly as phosphopantetheine and phosphopantethenate.

Another essential role of pantothenic acid is its participation in the 4'-phosphopantetheine moiety of acyl carrier protein (ACP), where the phosphodiester-linked prosthetic group uses the sulfhydryl terminus to exchange with malonyl-CoA to form an ACP-S malonyl thioester, which can chain elongate during fatty acid biosynthesis.

Pantothenic acid deficiency

Pantothenic acid deficiency is exceedingly rare. Because of its rarity, most information about pantothenic acid deficiency has been obtained from experiments: Pantothenic acid deficiency has been induced in humans by use of a metabolic antagonist, w-methyl pantothenic acid along with a pantothenic acid-deficient diet. Subjects became irascible and developed postural hypotension and rapid heart rate on exertion, epigastric distress with anorexia and constipation, numbness and tingling of the hands and feet. Because pantothenic acid is involved with so many vital processes in the body, it is not surprising that a broad number of complications might result from deficiency.

From recent research it is known that the pantothenic acid derivative, pantetheine (two molecules of pantetheine joined by a disulfide bond), has a hypocholesterolemic effect. A metabolic antagonist of pantothenic acid, pantoyl g-amino butyric acid (called pantoyl-GABA), is widely used in Japan as an antidementia drug for treating cognitive impairments in pathological states such as Alzheimer's disease, presumably through increasing cholinergic activity *in vivo*.

Indications

Suspicion of inadequate intake of pantothenic acid, e.g.

- dialysis patients
- alcohol abusius
- Crohn's disease, Colitis ulcerosa

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity	
			KIF004	KIF004.2
KIF004	PLATE	<i>Lactobacillus plantarum</i> -precoated microtiter plate	1 x	2 x
	DIL	Water	6 x 30 ml	7 x 30 ml
	ASYMED	Pantothenic acid assay medium	4 x	4 x
	STD	Pantothenic acid standard, lyophilized.	4 x	3 x
	FOL	Adhesive cover foil	1 x whole 3 x half	3 x whole
	FRA	Replacement holder for microtiter strips	1 x	1 x
	CTRL1	Pantothenic acid control 1, lyophilized.	4 x	3 x
	CTRL2	Pantothenic acid control 2, lyophilized.	4 x	3 x

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber, 37 °C
- Water bath (90 °C–100 °C)
- ELISA reader 610–630 nm (540–550 nm)
- Calibrated precision pipettors and sterile 20–1 000 µl tips
- 5 ml and 10 ml pipets
- 1.5–2 ml reaction vials, sterile
- 0.2 µm sterile polyethersulfone (PES) filter with a sterile disposable syringe
- 15 ml centrifuge tubes, sterile (e.g. Falcon tubes)
- Biocentrifuge (10 000 g)

5. STORAGE AND PREPARATION OF REAGENTS

- Store test kit and reagents at 2–8 °C.
- Prepare reagents freshly and use them immediately after preparation. Discard remaining unused reagents and waste in accordance with country, federal, state, and local regulations.
- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. Prepare only the appropriate amount necessary for each run. The kit can be used up to 4 times within the expiry date stated on the label.

5.1 Water

- Water [DIL] (for medium [ASYMED], standard [STD] and controls [CTRL1, CTRL2])
- Push the lid up and pull it back to the rim of the glass, then twist the whole cap off.

5.2 Preparation of the controls

- The lyophilised controls [CTRL1, CTRL2] have to be resuspended with each 1.25 ml water [DIL] from the test kit, then homogenise using a vortex.
- After reconstitution, the controls are treated like samples.
- The concentration of the controls changes from lot to lot and is stated in the product specification.

5.3 Preparation of the standard curve

- For the preparation of the standard curve, standard concentrate is needed. To prepare standard concentrate, resuspend the lyophilised standard [STD] with x ml (x = please see the enclosed quality control protocol for the needed volume) water [DIL] supplied with the test kit, then homogenise using a vortex.
- Prepare a standard curve in 6 sterile reaction tubes (1.5–2 ml volume) from standard concentrate and water [DIL] following the scheme depicted in the table below:

Pantothenic acid [µg/l]	Water [DIL] [µl]	+	Standard concentrate [µl]	=	Total volume [µl]
Blank: 0	975	+	0	=	975
Standard 1: 2.3	975	+	25	=	1 000
Standard 2: 4.6	950	+	50	=	1 000
Standard 3: 18.4	400	+	100	=	500
Standard 4: 27.6	350	+	150	=	500
Standard 5: 36.8	300	+	200	=	500

5.4 Preparation of the sterile assay medium

- Fresh sterile assay medium has to be prepared each time before performing a test.
- Remove lyophilised assay medium from the desiccant bag in the assay medium bottle by taking the bag with a forceps and shaking it whilst still inside the bottle. Then remove the clean desiccant bag and discard it.
- Add 10 ml water [DIL] to the assay medium bottle [ASYMED], close the bottle firmly and shake it. This amount is sufficient for 6 microtiter stripes.
- Heat the medium bottle in a water bath at 90–100 °C for 5 min, shake well at least 2 times during this incubation time. Take care that the medium bottle is always firmly closed.
- Quickly cool the medium bottle to < 30 °C (at 2–8 °C for 10 min).
- Filter the medium using a disposable syringe (10 ml) and the 0.2 µm PES filter into a sterile centrifuge tube (15 ml, e.g. Falcon).
- After this preparation, the sterile assay medium can be used in the test.

5.5 Microtiter plate [PLATE]

- Store the microtiter plate [PLATE] in the aluminium packaging containing the desiccant bag at 2–8 °C.
- The microtiter plate [PLATE] has to be protected from humidity and contamination.
- Take care that the aluminium packaging is not damaged.
- Carefully close the aluminium packaging after opening.

- Take only the microtiter stripes needed directly before usage to avoid contamination

6. SAMPLE STORAGE AND PREPARATION

- Use serum for analysis.
- Samples are stable at 2–8 °C for 3 days in the dark. For longer storage, samples should be frozen and kept at -20 °C.
- Hemolytic samples may give erroneous results and should not be used for analysis. Lipemic samples should be centrifuged at 13 000 *g* before assaying to obtain fat free serum as far as possible.
- Samples should be centrifuged (at least 5 min at 10 000 *g*) prior to measurement. Use the resulting supernatant in the test.

6.1 Sample dilution

Take 50 µl from the sample/control, add 350 µl water [DIL] and mix. The sample treatment and dilution result in a total dilution of 1:8 (= sample dilution factor).

7. ASSAY PROCEDURE

7.1 Principle of the test

The serum samples are diluted and then transferred into the wells of a microtiter plate coated with *Lactobacillus plantarum*. The addition of pantothenic acid in either standards or samples gives a pantothenic acid-dependent growth response until pantothenic acid is consumed. After incubation at **37 °C** for **24 h**, the growth of *Lactobacillus plantarum* is measured turbidimetrically at 610–630 nm (alternatively at 540–550 nm) in an ELISA reader and compared to a standard curve generated from the dilution series. The amount of pantothenic acid is directly proportional to the turbidity.

7.2 Test preparations

Take as many microtiter strips as needed from kit. Return unused strips and any unused test kit component to the original packaging, and put in the refrigerator. Bring all necessary reagents to room temperature.

7.3 Test procedure

- Take as many microtiter strips as needed from the kit and put them in the second microtiter strip holder [FRA].
- Put 150 µl sterile assay medium into the cavities.
- Add each 150 µl of the prepared standard curve, samples and controls into the respective cavities. Pre-rinse each pipet tip with standard, control or sample solution, respectively.
- Carefully seal the plate with adhesive cover foil [FOL]. Important: the cavities must be made airtight by pressing the foil down with the hand!
- Keep at 37 °C for 24 h in an incubator.

7.4 Measurement

- Press the adhesive cover foil [FOL] firmly down again with the hand.
- Upturn the microtiter plate [PLATE], put it onto a tabletop and shake the microbes well.
- Turn the microtiter plate [PLATE] over again and carefully remove the adhesive cover foil [FOL]. During this, fix the strips in the frame with your hand because the foil is highly adhesive.
- Remove air bubbles in the cavities using a pipet tip or a needle.
- Read turbidity in an ELISA reader at E 610–630 nm (alternatively at E 540–550 nm).

Please note

- After 24 h incubation time, the microtiter plate [PLATE] may be stored for a maximum of 48 h in the refrigerator before measuring the turbidity.
- To prevent time-loss through public holidays or weekends, the microtiter plate [PLATE] may also be evaluated after 65 h incubation.

8. EVALUATION OF RESULTS

We recommend to use the 4 parameter algorithm to calculate the results. The sample dilution factor has to be considered for data evaluation.

The blank should have an optical density < standard 1. It serves as optical control to exclude contaminations and is not included in the calculation of results.

8.1 Calculation

Pantothenic acid in $\mu\text{g/l}$ = value from the standard curve \times sample dilution factor (8)

8.2 Expected values

The concentration of pantothenic acid was determined in 74 samples of different blood donors. The median value was 91.4 (81.4) $\mu\text{g/l}$. The 2-SD area was 36–147 $\mu\text{g/l}$. Figure 1 shows the distribution of the values.

Number of samples	74
Mean	91.4
Median	81.4
SD	27.7
MW-2*SD	36.0
MW+2*SD	146.8


Fig. 1: Distribution of pantothenic acid values in blood donor samples

Please note

A concentration range of 18.4–294.4 $\mu\text{g/l}$ pantothenic acid is covered at a sample dilution of 1:8.

We recommend each laboratory to develop its own normal range as normal ranges strongly depend on the choice of the patient collective. The values mentioned above are only for orientation and can deviate from other published data.

8.3 Quality control

The extinction of the highest standard has to be > 0.6 .

Results, generated from the analysis of control samples, should be evaluated for acceptability. The results for the samples may not be valid if within the same assay one or more values of the quality control sample or the highest standard are outside the acceptable limits.

9. PERFORMANCE CHARACTERISTICS

The following performance characteristics have been collected using human serum samples.

9.1 Precision and reproducibility

Intraassay (n = 28)

	Pantothenic acid [$\mu\text{g/l}$]	CV [%]
Sample	81.0	3.0

Interassay (n = 5)

	Pantothenic acid [$\mu\text{g/l}$]	CV [%]
Sample	92.4	4.9

9.2 Recovery

Samples from 3 patients were spiked with pantothenic acid and analysed. The mean values are shown below.

Sample (n=5)	Mean value original sample [$\mu\text{g/l}$]	Spike [$\mu\text{g/l}$]	Pantothenic acid expected [$\mu\text{g/l}$]	Pantothenic acid measured [$\mu\text{g/l}$]	Recovery Rate [%]
A	112.5	18.4	130.9	131.9	105
		36.8	149.3	141.5	79
		55.2	167.7	158.1	83
Recovery rate in total [%]					89

Sample (n=5)	Mean value measured in original sample [$\mu\text{g/l}$]	Spike [$\mu\text{g/l}$]	Pantothenic acid expected [$\mu\text{g/l}$]	Pantothenic acid measured [$\mu\text{g/l}$]	Recovery Rate [%]
B	96.61	18.4	115.0	113.8	93
		36.8	133.4	133.8	101
		55.2	151.8	165.3	125
Recovery rate in total [%]					106

Sample (n=5)	Mean value measured in original sample [µg/l]	Spike [µg/l]	Pantothenic acid expected [µg/l]	Pantothenic acid measured [µg/l]	Recovery Rate [%]
C	106.21	18.4	124.6	122.5	88
		36.8	143.0	138.6	88
		55.2	161.4	176.1	127
Recovery rate in total [%]					101

10. LIMITATIONS

Whole blood cannot be used in the assay.

11. PRECAUTIONS

- As the test is based on a microbiological method, the general guidelines for sterile work should be observed as far as possible (preferably work in a sterile bench / PCR hood, use of sterile instruments or equipment).
- GLP (Good Laboratory Practice) guidelines have to be observed.
- Water quality is extremely important for the test. Only the water delivered with the test kit [DIL] should be used.
- For sterile filtration, only a sterile polyethersulfone filter must be used.
- It is essential to run a standard curve for each separate assay.
- Controls should be measured with each assay.
- We recommend measurements in duplicate.
- If a higher dilution results in a higher value measured, inhibitors like antibiotics might be present.
- Reagents should not be used beyond the expiration date shown on the label.
- Wear gloves during the test.
- Used microtiter stripes [PLATE] and materials that have been in contact with patient samples should be handled and disposed as potentially infectious.

12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- ID-Vit® is a trademark of Immundiagnostik AG.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Control samples should be analysed with each run.
- The assay should always be performed according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

13. REFERENCES

1. Burtis, C.A. & Ashwood, E.R., 1999. Tietz textbook of clinical chemistry 3rd ed., W.B. Saunders.
2. Coronel, F. et al., 1991. Treatment of hyperlipemia in diabetic patients on dialysis with a physiological substance. *American journal of nephrology*, **11**(1), pp.32–6.

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com