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Manual
For professional use only

ID-Vit[®] Folic acid in whole blood

*Microbiological test kit for the determination of folic acid
in whole blood using a Lactobacillus rhamnosus-coated
microtiter plate
For use in human and veterinary medicine and in research*


Valid from 2024-05-07

REF KIF005VB

REF KIF005VB.2

Σ 96

Σ 2 x 96

IVD

CE

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Safety information

These accessories are to be used exclusively in accordance with the enclosed instructions for use. Important safety information for this product can be found in chapter 6.

Table of Contents

1. INTENDED PURPOSE	19
2. INTRODUCTION	19
3. PRINCIPLE OF THE TEST	20
4. MATERIAL SUPPLIED	21
5. MATERIAL REQUIRED BUT NOT SUPPLIED	22
6. PRECAUTIONS	22
7. STORAGE AND PREPARATION OF REAGENTS	22
7.1 Water	23
7.2 Preparation of the sterile assay medium	23
7.3 Preparation of the lysis-reagent	23
7.4 Preparation of the controls	23
7.5 Preparation of the standard curve	24
7.6 Microtiter plate (PLATE)	24
8. SAMPLE STORAGE AND PREPARATION	24
8.1 Sample pretreatment	25
9. ASSAY PROCEDURE	25
9.1 Test preparations	25
9.2 Test procedure	25
9.3 Measurement	26
10. EVALUATION OF RESULTS	26
10.1 Calculation	26
10.2 Quality control	27
11. LIMITATIONS	27
12. PERFORMANCE CHARACTERISTICS	27
12.1 Precision and reproducibility	27
12.2 Recovery	27
12.3 Linearity	29
13. GENERAL NOTES ON THE TEST	30
14. REFERENCES	31
15. SYMBOLS	31

1. INTENDED PURPOSE

ID-Vit® Folic acid is a microtiter plate test kit based on a microbiological method for determining the total content of folic acid from whole blood. All required reagents, as well as the standard are included in the test. An ELISA reader is required for the evaluation. For use in human and veterinary medicine and in research. For *in vitro* diagnostic use only.

2. INTRODUCTION

Folic acid (vitamin B₉) belongs to the group of the vitamin B complex; it is water-soluble, sensitive to light, oxygen, increased temperatures and is essentially involved in all growth and development processes of the body. Folic acid is necessary for the formation of red blood cells, for optimal bone marrow function and healthy nerve activity. Folic acid is also essential for all cell divisions (hence its importance in foetal development).

Although folic acid is found in most plant and animal foods an undersupply of folic acid is the most common vitamin deficiency in Europe and North America. According to the German Nutrition Society only one in four Germans consumes enough folic acid - often the result of an unbalanced diet with little fresh fruit and vegetables. But age, various diseases and the intake of certain medicines such as cotrimoxazole can also lead to absorption disorders and an associated deficiency. Lowered folic acid levels occur because of:

- a decreased supply (e.g. folic acid antagonists).
- a disrupted resorption (e.g. in celiac disease. CED).
- an increased requirement (e.g. increased alcohol consumption. when taking contraceptives. during pregnancy. anaemia or cancer).

Symptoms of Deficiency

The early symptoms of deficiency are fatigue, irritability, depression, poor concentration and loss of appetite; other consequences are inflammation of the mucous membranes, anaemia and severe neurological damage. During pregnancy, when the need for folic acid doubles, a deficiency can lead to premature birth and severe malformations. Optimal folic acid supply during pregnancy can reduce the risk of neural tube defects in the foetus by 85%. Since both a folic acid deficiency and a vitamin B₁₂ deficiency can cause megaloblastic anaemia, the determination of both vitamins is particularly important in this clinical picture in order to supplement the right vitamin. Treatment of megaloblastic anaemia caused by vitamin B₁₂ deficiency with folic acid alone can lead to irreversible damage to the central nervous system.

Folic acid and arteriosclerosis

A folic acid deficiency is known to be the most common cause of hyperhomocysteinaemia. Hyperhomocysteinaemia has been recognised as an independent factor in arteriosclerosis. Therefore, the determination of folic acid can be carried out within the framework of a coronary disease risk analysis. Apart from the influence of folic acid on the homocysteine levels a further positive effect on the endothelial function in heart patients has been established – development of nitrate tolerance during continuous nitrate therapy, e. g. in such patients, an increased release of oxygen radicals occurs without folic acid supplementation (Verhaar et al. 2002).

Indications

- Hyperchrome, macrocytic anemia
- Long-term therapy with antiepileptic drugs or folic acid antagonists
- Long-term haemodialysis
- (Multiple) pregnancy / planned pregnancy
- Enhanced erythropoiesis
- Chronic liver diseases
- Hemoblastosis
- Psoriasis, dermatitis
- Stomatitis, glossitis
- Chronic alcohol abus

3. PRINCIPLE OF THE TEST

The EDTA whole blood samples are pre-treated and diluted with a buffer mixture and then transferred into the wells of a microtiter plate precoated with *Lactobacillus rhamnosus*. The addition of folic acid in either standards or samples gives a folic acid-dependent growth response until folic acid is consumed. After incubation at **37 °C** for **46–50 h** the growth of *Lactobacillus rhamnosus* is measured turbidimetrically at 610–630 nm (alternatively at 540–550 nm) in an ELISA reader and compared to a standard curve generated from the dilution series. The amount of folic acid is directly proportional to the turbidity.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity	
			KIF005VB	KIF005VB.2
KIF005VB/ KIF005VB.2	DIL	Water	4 x 30 ml	7 x 30 ml
	PLATE	microtiter plate, precoated with <i>Lactobacillus rhamnosus</i>	1 x	2 x
	SOL	Sample treatment solution, ready-to-use	5 x 4.6 ml	10 x 4.6 ml
	PAF	Folic acid buffer, lyophilized	5 x	10 x
	ASYMED	Folic acid assay medium, ready-to-use	4 x	4 x
	ASYBUF	Folic acid medium treat- ment buffer, ready-to-use	4 x 1.5 ml	4 x 1.5 ml
	STD	Folic acid standard, lyophilized	4 x	3 x
	FOL	Adhesive cover foil	1 x whole 3 x half	3 x whole
	CTRL1	Folic acid control 1, lyophilized	4 x	3 x
	CTRL2	Folic acid control 2, lyophilized	4 x	3 x
	FRA	Replacement holder for microtiter strips	1 x	1 x

For reorders of single components use the catalogue number followed by the label as product number.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber, 37°C
- ELISA reader 610–630 nm (540–550 nm)
- Calibrated precision pipettors and sterile single use 20–1 000 µl tips
- 5 ml and 10 ml pipets
- 1.5–2 ml reaction vials
- 0.2 µm sterile polyethersulfone (PES) filter with a disposable syringe (10 ml)
- 15 ml centrifuge tubes (e.g. Falcon tubes)
- Biocentrifuge (10 000 g)
- Vortex

6. PRECAUTIONS

- The test is based on a microbiological method. Contaminations lead to erroneous results.
- Water quality is extremely important for the test. Use only the water delivered with the test kit (**DIL**).
- For sterile filtration, only a sterile polyethersulfone filter must be used.
- It is essential to run a standard curve for each separate assay.
- Measure controls with each assay.
- We recommend measurements in duplicate.
- Do not use reagents beyond the expiration date shown on the label.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- Used microtiter stripes and materials that have been in contact with patient samples must be handled and disposed of as potentially infectious

7. STORAGE AND PREPARATION OF REAGENTS

- Store test kit and reagents at 2–8°C.
- Prepare reagents freshly and use them immediately after preparation.
- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. Prepare only the appropriate amount necessary for each run. The kit can be used up to 3x (KIF005VB.2) or 4x (KIF005VB) within the expiry date stated on the label.

7.1 Water

- Water (**DIL**) for medium (**ASYMED**), standard (**STD**), controls (**CTRL1**, **CTRL2**) and dilutions.
- Push the lid up and pull it back to the rim of the glass, then twist the whole cap off.

7.2 Preparation of the sterile assay medium

- Fresh sterile assay medium has to be prepared each time before performing a test.
- Remove the desiccant bag from the lyophilised assay medium bottle by taking the bag with a forceps and shaking it whilst still inside the bottle. Then remove the clean desiccant bag and discard it.
- Add 10 ml water (**DIL**) to the assay medium bottle (**ASYMED**), close the bottle firmly and vortex well. This amount is sufficient for 6 microtiter stripes.
- Filter the medium using a disposable syringe (10 ml) and the 0.2 µm PES filter into a centrifuge tube (15 ml, e.g. Falcon).
- After this preparation, the sterile assay medium can be used in the test.

7.3 Preparation of the lysis-reagent

- Transfer 4.5 ml sample preparation buffer (**SOL**) into vials containing lyophilised buffer (**PAF**), then homogenise using a vortex.
- The lysis-reagent cannot be stored.

7.4 Preparation of the controls

- The lyophilised controls (**CTRL1**, **CTRL2**) have to be resuspended each with **125 µl** lysis-reagent (corresponds to a 1:10 dilution), then homogenise using a vortex.
- After reconstitution, the controls are treated like samples.
- The concentration of the controls changes from batch to batch and can be found in the product specification.

7.5 Preparation of the standard curve

- For the preparation of the standard curve, standard concentrate is needed. To prepare standard concentrate, resuspend the lyophilised standard (**STD**) with x ml water (**DIL**) (x = see quality control protocol) supplied with the test kit, then homogenise using a vortex.
- Prepare a standard curve in 6 sterile reaction tubes (1.5–2 ml volume) from standard concentrate and water (**DIL**) following the scheme depicted in the table below:

Folic acid [µg/l]	Water (DIL) [µl]	+	Standard concentrate [µl]	=	Total volume [µl]
Blank: 0	500	+	0	=	500
Standard 1: 0,04	450	+	50	=	500
Standard 2: 0,12	350	+	150	=	500
Standard 3: 0,20	250	+	250	=	500
Standard 4: 0,28	150	+	350	=	500
Standard 5: 0,36	50	+	450	=	500

7.6 Microtiter plate (PLATE)

- Store the microtiter plate (**PLATE**) in the aluminium packaging containing the desiccant bag at 2–8°C.
- The microtiter plate (**PLATE**) has to be protected from humidity and contamination.
- Take care that the aluminium packaging is not damaged.
- Carefully close the aluminium packaging after opening.
- Take only the microtiter stripes needed directly before usage to avoid contamination.

8. SAMPLE STORAGE AND PREPARATION

- Use EDTA whole blood for analysis.
- Samples are stable at 2–8°C for 8 hours in the dark. For longer storage, samples can be frozen and kept at -20°C for up to 5 months.

8.1 Sample pretreatment

Add 25 µl of EDTA whole blood sample to 225 µl of the prepared lysis-reagent (ratio 1:10). Transfer the controls into a 1.5 ml reaction tube. Incubate samples and controls at 37°C for 30 min. Cool at 2–8°C for 10 min and then centrifuge (10 000 g) for 10 min.

8.2 Sample dilution

Take 10 µl from the supernatant of the pre-treated EDTA whole blood and add 740 µl water [DIL] and mix. The sample pretreatment and dilution corresponds in total to a **1:750 dilution** (= sample dilution factor).

Take 10 µl of the supernatant of the pre-treated **controls** and add 740 µl water [DIL] and mix. The sample pre-treatment and dilution corresponds in total to a **1:750 dilution** (= sample dilution factor).

9. ASSAY PROCEDURE

9.1 Test preparations

Take as many microtiter strips as needed from the kit. Return unused strips and any unused test kit components to the original packaging, and store in the refrigerator. Bring all necessary reagents to room temperature.

9.2 Test procedure

- Take as many microtiter strips as needed from the kit and put them in the second microtiter strip holder (**FRA**).
- Put 150 µl sterile assay medium into each cavity.
- Add 150 µl of the prepared standard dilutions (blank, standard 1–5), samples and controls into the respective cavities. Pre-rinse each pipet tip with standard, control or sample solution, respectively.
- Carefully seal the plate with adhesive cover foil (**FOL**). Important: the cavities must be made airtight by pressing the foil down with the hand!
- Keep at **37°C** for **46–50 h** in an incubator.

9.3 Measurement

- Press the adhesive cover foil (**FOL**) firmly down again with the hand.
- Turn the microtiter plate (**PLATE**) upside down, place it onto a tabletop and shake the microbes well.
- Turn the microtiter plate (**PLATE**) over again and carefully remove the adhesive cover foil (**FOL**). During this, fix the strips in the frame with your hand because the foil is highly adhesive.
- Remove air bubbles in the cavities using a pipet tip or a needle.
- Read turbidity in an ELISA reader at E 610–630 nm (alternatively at E 540–550 nm).

Please note

- After 46–50 h incubation time the microtiter plate (**PLATE**) may be stored for a maximum of 48 h in the refrigerator before measuring the turbidity.

10. EVALUATION OF RESULTS

We recommend to use the 4 parameter algorithm to calculate the results. The sample dilution factor has to be considered for data evaluation.

The blank serves as a visual control to exclude contamination and is not taken into account in the calculation. The optical density must be < standard 1. If this is not the case, the analysis must be carried out again.

10.1 Calculation

Folic acid in $\mu\text{g/l}$ = value from the standard curve \times sample dilution factor (750).

Reference value for EDTA whole blood sample

Based on internal studies of EDTA whole blood samples of apparently healthy persons (n = 74) the following values were estimated.

Folic acid: 79.5–231.4 $\mu\text{g/l}$

Please note

With a sample dilution factor of 750 a range from 30–270 $\mu\text{g/l}$ folic acid is covered.

We recommend each laboratory to establish its own reference range as reference ranges are strongly dependent on the selection of the patient collective. The reference range for folic acid is given for guidance only and may differ from other published data.

10.2 Quality control

The extinction of the highest standard has to be > 0.6.

Results, generated from the analysis of control samples, should be evaluated for acceptability. The results for the samples may not be valid if within the same assay one or more values of the quality control sample or the highest standard are outside the acceptable limits.

11. LIMITATIONS

Only EDTA whole blood can be used for the test.

12. PERFORMANCE CHARACTERISTICS

The following performance characteristics have been collected using human EDTA whole blood samples.

12.1 Precision and reproducibility

Intraassay (n = 22)

Sample	Folic acid [$\mu\text{g/l}$]	CV [%]
1	70.24	7.2
2	146.03	8.5
3	239.31	5.6

Interassay (n = 33)

Sample	Folic acid [$\mu\text{g/l}$]	CV [%]
1	84.53	7.6
2	149.08	7.7
3	227.9	8.2

12.2 Recovery

For the determination of recovery 4 EDTA whole blood samples were spiked with folic acid in different concentrations and measured. The result is obtained by the ratio between the measurement result and the true concentration of a sample.

Sample (n=9)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
A	99.2	40.8	132.98	129.86	98
		76.8	168.98	186.44	110
		100.3	192.48	216.34	112
Recovery rate in total [%]					107

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
B	107.1	40.8	147.89	167.37	113
		76.8	183.89	195.19	106
		100.3	207.39	243.72	118
Recovery rate in total [%]					112

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
C	108.4	40.8	149.21	135.01	90
		76.8	185.21	168.53	91
		100.3	208.71	215.19	103
Recovery rate in total [%]					95

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
D	115.6	40.8	156.38	170.08	109
		76.8	192.38	211.48	110
		100.3	215.88	215.58	100
Recovery rate in total [%]					106

12.3 Linearity

EDTA whole blood samples from four patients were serially diluted and analysed by two operators. The results are shown below.

Sample	Dilution		Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
A	1:425	1	316.7	303.2	96
	1:500		271.7	259.5	95
	1:750		194.1	194.1	100
	1:1 000		145.6	153.8	106
	1:1 250		116.5	113.2	97
B	1:425	2	257.1	266.6	104
	1:500		218.6	226.5	104
	1:750		145.7	145.7	100
	1:1 000		109.3	117.8	108
	1:1 250		87.4	94.1	108
C	1:425	1	279.9	275.4	98
	1:500		237.9	219.4	92
	1:750		158.6	158.6	100
	1:1 000		119.0	120.3	101
	1:1 250		95.2	96.9	102

Sample	Dilution		Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
D	1:425	2	265.8	267.6	101
	1:500		225.9	227.9	101
	1:750		150.6	150.6	100
	1:1 000		113.0	118.1	105
	1:1 250		90.4	95.5	106

13. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- ID-Vit® is a trademark of Immundiagnostik AG.
- Do not use reagents beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Follow the guidelines for medical laboratories.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be made within 14 days after reception of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Analyse controls with each run.
- Always perform assay according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

14. REFERENCES

1. Obeid. R. & Herrmann. W.. 2006. Mechanisms of homocysteine neurotoxicity in neurodegenerative diseases with special reference to dementia. *FEBS letters*. **580**(13). pp.2994–3005.
2. Strohecker. R. & Henning. H.. 1963. Vitamin-Bestimmungen. Erprobte Methoden. E. Merck AG. ed.. Weinheim/Bergstraße: Verlag Chemie GmbH.
3. Verhaar. M.C.. Stroes. E. & Rabelink. T.J.. 2002. Folates and cardiovascular disease. *Arteriosclerosis, thrombosis, and vascular biology*. **22**(1). pp.6–13.

15. SYMBOLS

	Temperature limitation		Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device		To be used with
	Manufacturer		Content sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com