ID-Vit[®] Vitamin B₁₂

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For use in human and veterinary medicine and in research


Valid from 2023-01-26

REF KIF012

REF KIF012.2

Σ 96

Σ 2 x 96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	16
2. INTRODUCTION	16
3. PRINCIPLE OF THE TEST	17
4. MATERIAL SUPPLIED	18
5. MATERIAL REQUIRED BUT NOT SUPPLIED	18
6. STORAGE AND PREPARATION OF REAGENTS	19
6.1 <i>Water</i>	19
6.2 <i>Preparation of the stabilisation reagent</i>	19
6.3 <i>Preparation of the controls</i>	19
6.4 <i>Preparation of the standard curve</i>	20
6.5 <i>Preparation of the sterile assay medium</i>	20
6.6 <i>Microtiter plate</i>	21
7. STORAGE AND PREPARATION OF SAMPLES	21
7.1 <i>Sample pretreatment</i>	21
7.2 <i>Sample dilution</i>	21
8. ASSAY PROCEDURE	22
8.1 <i>Test preparations</i>	22
8.2 <i>Test procedure</i>	22
8.3 <i>Measurement</i>	22
9. EVALUATION OF RESULTS	23
9.1 <i>Calculation</i>	23
9.2 <i>Quality control</i>	23
10. LIMITATIONS	24
11. PERFORMANCE CHARACTERISTICS	24
11.1 <i>Precision and reproducibility</i>	24
11.2 <i>Recovery</i>	24
12. PRECAUTIONS	25
13. GENERAL NOTES ON THE TEST	26
14. REFERENCES	27

1. INTENDED USE

ID-Vit® Vitamin B₁₂ is a microtiter plate test kit based on a microbiological method which measures the total vitamin B₁₂ content in serum. The test kit contains the standard and all reagents required to perform the test. It is sufficient for 96 determinations including standard curves. An ELISA reader is required for the evaluation of the results. For use in human and veterinary medicine and in research. For *in vitro* diagnostic use only.

2. INTRODUCTION

Vitamin B₁₂ (cobalamin), a collective term for a group of various substituted corrinoids with cobalt as the central atom, is found free and also protein-bound in food. The protein-bound form is degraded by pancreatic protease, releasing free B₁₂ which binds to intrinsic factor, a protein secreted by gastric parietal cells of the stomach mucosa. The cobalamin-intrinsic factor complex is bound to mucous membrane cells of the stomach in the ileum and absorbed by the cells. In the case of high doses, a diffusion of the complex also takes place. Vitamin B₁₂ is bound to the protein transcobalamin II (TC-II) within the cells. TC-II serves as a transport protein for vitamin B₁₂ in the circulation system.

Vitamin B₁₂ is involved as a co-enzyme in metabolic processes and plays an important role in the formation of the blood, the development of the nervous system and the regeneration of the mucous membranes. In addition, there is a direct relationship to the formation of folic acid because methylcobalamin is involved in the transfer of methyl groups for the synthesis of methionin from homocystein.

Vitamin B₁₂ deficiency

Vitamin B₁₂ deficiency is rarely caused by dietary factors. In most cases, it results from a resorption disorder of the intestines or defective development of intrinsic factor. Since vitamin B₁₂ resorption can be reduced up to 50% in the elderly, an increased supplement is recommended. Pregnant women with a lacto-vegetarian diet are also recommended to increase their intake because their vitamin B₁₂ stores in the liver could be exhausted.

The classical vitamin B₁₂ deficiency disease is pernicious anemia. In the early stages of the disease, vitamin B₁₂ deficiency symptoms are manifested as weariness, palpitations, pallor or jaundice.

Indications for vitamin B₁₂ determination

- Megaloblastic (pernicious) anemia
- Hyperhomocysteinaemia (patients on dialysis, old people)
- Homocysteinuria
- Peripheral neuropathy
- Patients with CED, gastritis, gastrectomy, gluten intolerance or intestinal re-sorption disorders
- Pancreatic insufficiency
- Patients with thrombosis
- Alcoholism
- Chronic liver and kidney disease
- Vitamin B₁₂ deficiency from diet (vegan vegetarians)
- Pregnancy and lactation

3. PRINCIPLE OF THE TEST

The serum samples are pre-treated and diluted with a buffer mixture, and then transferred into the wells of a microtiter plate coated with *Lactobacillus delbrueckii* subsp. *lactis*. The addition of vitamin B₁₂ in either standards or samples gives a vitamin B₁₂-dependent growth response until vitamin B₁₂ is consumed. After incubation at **37°C** for **48 h**, the growth of *Lactobacillus delbrueckii* subsp. *lactis* is measured turbidimetrically at 610–630 nm (alternatively at 540–550 nm) in an ELISA reader and a standard curve is generated from the dilution series. The amount of vitamin B₁₂ is directly proportional to the turbidity.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity	
			KIF012	KIF012.2
KIF012	PLATE	<i>Lactobacillus delbrueckii</i> subsp. <i>lactis</i> -precoated microtiter plate	1 plate	2 plate
	SOL	Sample preparation buffer	7 x 5 ml	14 x 5 ml
	STAB	Stabilizer; lyophilised	7 x	14 x
	DIL	Water	4 x 30 ml	5 x 30 ml
	ASYMED	Vitamin B ₁₂ assay medium	4 x	4 x
	STD	Vitamin B ₁₂ standard; lyophilised	4 x	3 x
	FOL	Adhesive cover foil	1 x whole 3 x half	3 x whole
	FRA	Replacement holder for 96 well plates	1 x	1 x
	CTRL1	Vitamin B ₁₂ control 1; lyophilised	4 x	3 x
CTRL2	Vitamin B ₁₂ control 2; lyophilised	4 x	3 x	

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber, 37 °C
- Water bath (90 °C–100 °C)
- optional for sample preparation: thermoblock (95 °C)
- ELISA reader 610–630 nm (540–550 nm)
- Calibrated precision pipettors and 20–1 000 µl tips
- Pipettes of 5 and 10 ml
- 1.5–2 ml reaction vials, sterile
- 0.2 µm sterile polyethersulfone (PES) filter with a sterile disposable syringe
- 15 ml centrifuge tubes, sterile (e.g. Falcon tubes)
- Biocentrifuge (10 000 g)

6. STORAGE AND PREPARATION OF REAGENTS

- Store test kit and reagents at 2–8 °C.
- Prepare reagents freshly and use them immediately after preparation. Discard remaining unused reagents and waste in accordance with country, federal, state, and local regulations.
- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. Prepare only the appropriate amount necessary for each run. The kit can be used up to 4 times within the expiry date stated on the label.

6.1 Water

- Water [DIL] (for medium [ASYMED], standard [STD] and controls [CTRL1, CTRL2])
- Push the cap up and pull it back to the rim of the glass, then twist the whole cap off.

6.2 Preparation of the stabilisation reagent

- Add 4 ml sample preparation buffer [SOL] to a vial of lyophilised stabilisation reagent [STAB], then homogenise using a vortex.
- Stabilisation reagent cannot be stored.

6.3 Preparation of the controls

- The lyophilised controls [CTRL1, CTRL2] have to be resuspended with each 300 µl water [DIL] from the test kit, then homogenise using a vortex.
- After reconstitution, the controls are treated like samples.
- The concentration of the controls changes from lot to lot and is stated in the product specification.

6.4 Preparation of the standard curve

- For the preparation of the standard curve, standard concentrate is needed. To prepare standard concentrate, resuspend the lyophilised standard [STD] with x ml (x = please see the enclosed quality control protocol for the volume needed) water [DIL] supplied with the test kit, then homogenise using a vortex.
- In 6 sterile reaction tubes (1.5–2 ml volume), prepare a standard curve from standard concentrate and water [DIL] following the scheme depicted in the table below:

Vitamin B ₁₂ [ng/l]	Water [DIL] [µl]	+	Standard concentrate [µl]	=	Total volume [µl]
Blank: 0	700	+	0	=	700
Standard 1: 6	700	+	50	=	750
Standard 2: 18	400	+	100	=	500
Standard 3: 27	350	+	150	=	500
Standard 4: 36	300	+	200	=	500
Standard 5: 54	200	+	300	=	500

6.5 Preparation of the sterile assay medium

- Fresh sterile assay medium has to be prepared each time before performing a test.
- Remove lyophilised assay medium from the desiccant bag in the assay medium bottle by taking the bag with a forceps and shaking it whilst still inside the bottle. Then remove the clean desiccant bag and discard it.
- Add 10 ml water [DIL] to the assay medium bottle [ASYMED], close the bottle firmly and shake it. This amount is sufficient for 6 microtiter stripes.
- Heat the medium bottle in a water bath at 90–100 °C for 5 min, shake well at least 2 times during this incubation time. Take care that the medium bottle is always firmly closed.
- Quickly cool the medium bottle to < 30 °C (at 2–8 °C for 10 min).
- Filter the medium using a disposable syringe (10 ml) and the 0,2 µm PES filter into a sterile centrifuge tube (15 ml, e.g. Falcon).
- After this preparation, the assay medium can be used in the test.

6.6 Microtiter plate [PLATE]

- Store the microtiter plate [PLATE] in the aluminium packaging containing the desiccant bag at 2–8 °C.
- The microtiter plate [PLATE] has to be protected from humidity and contamination.
- Take care that the aluminium packaging is not damaged.
- Carefully close the aluminium packaging after opening.
- Take only the microtiter stripes needed directly before usage to avoid contamination

7. STORAGE AND PREPARATION OF SAMPLES

- Use serum for analysis.
- Samples are stable at 2–8 °C for one day in the dark. For longer storage, samples should be frozen and kept at -20 °C.
- Hemolytic samples may give erroneous results and should not be used for analysis. Lipemic samples should be centrifuged at 13 000 *g* before assaying to obtain fat free serum as far as possible.
- Samples should be centrifuged (at least 5 min at 10 000 *g*) prior to measurement. Use the resulting supernatant in the test.

7.1 Sample pretreatment

Add 75 µl serum/control to 300 µl of prepared stabilising solution (ratio 1:5), mix, heat to 95 °C for 30 min and then cool quickly (10 min at 2–8 °C). Afterwards, centrifuge for 10 min at 10 000 *g*.

7.2 Sample dilution

Take 100 µl from the supernatant of the prepared serum/control, add 400 µl water [DIL] and mix. The sample treatment and dilution result in a total dilution of 1:25 (= sample dilution factor).

8. ASSAY PROCEDURE

8.1 Test preparations

Take as many microtiter strips as needed from kit. Return unused strips and any unused test kit component to the original packaging, and put in the refrigerator. Bring all necessary reagents to room temperature.

8.2 Test procedure

- Take as many microtiter strips as needed from the kit and put them in the second microtiter strip holder [FRA].
- Put 150 µl sterile assay medium into the cavities.
- Add each 150 µl of the prepared standard curve, samples and controls into the respective cavities. Pre-rinse each pipet tip with standard, control or sample solution, respectively.
- Carefully seal the plate with adhesive cover foil [FOL]. Important: the cavities must be made airtight by pressing the foil down with the hand!
- Keep at 37 °C for 48 h in an incubator.

8.3 Measurement

- Press the adhesive cover foil [FOL] firmly down again with the hand.
- Upturn the microtiter plate [PLATE], put it onto a tabletop and shake the microbes well.
- Turn the microtiter plate [PLATE] over again and carefully remove the adhesive cover foil [FOL]. During this, fix the strips in the frame with your hand because the foil is highly adhesive.
- Remove air bubbles in the cavities using a pipet tip or a needle.
- Read turbidity in an ELISA reader at E 610–630 nm (alternatively at 540–550 nm).

Please note

- After 48 h incubation time, the microtiter plate [PLATE] may be stored for a maximum of 48 h in the refrigerator before measuring the turbidity.
- To prevent time-loss through public holidays or weekends, the microtiter plate [PLATE] may also be measured after 60 h incubation.

9. EVALUATION OF RESULTS

We recommend to use the 4 parameter algorithm to calculate the results. The sample dilution factor has to be considered for data evaluation.

The blank should have an optical density < standard 1. It serves as optical control to exclude contaminations and is not included in the calculation of results.

9.1 Calculation

Vitamin B₁₂ in ng/l = value from the standard curve × sample dilution factor (25)

Reference value for human serum

Based on studies of matrix samples of apparently healthy persons (n = 83) the following values were estimated.

Vitamin B₁₂:

200–830 ng/l (≅ 147.2–626.5 pmol/l)

Please note

A concentration range of 150–1 350 ng/l vitamin B₁₂ is covered at a sample dilution of 1:25.

We recommend each laboratory to develop its own normal range as normal ranges strongly depend on the choice of the patient collective. The values mentioned above are only for orientation and can deviate from other published data.

9.2 Quality control

The extinction of the highest standard has to be > 0.6.

Results, generated from the analysis of control samples, should be evaluated for acceptability. The results for the samples may not be valid if within the same assay one or more values of the quality control sample or the highest standard are outside the acceptable limits.

10. LIMITATIONS

Whole blood, EDTA plasma, and heparin plasma cannot be used in the assay.

11. PERFORMANCE CHARACTERISTICS

The following performance characteristics have been collected using human serum samples.

11.1 Precision and reproducibility

Intraassay (n = 21)

	Vitamin B ₁₂ [ng/l]	CV [%]
Sample	294	5.38

Interassay (n = 3)

	Vitamin B ₁₂ [ng/l]	CV [%]
Sample	285	8.0

11.2 Recovery

Samples from 3 patients were spiked with vitamin B₁₂ and analysed. The mean values are shown below.

Sample (n=5)	Mean value original sample [ng/l]	Spike [ng/l]	Vitamin B ₁₂ expected [ng/l]	Vitamin B ₁₂ measured [ng/l]	Recovery Rate [%]
A	566.58	187.5	754.08	726.36	85
		375.0	941.58	908.21	91
Recovery rate in total [%]					88

Sample (n=5)	Mean value measured in original sample [ng/l]	Spike [ng/l]	Vitamin B ₁₂ expected [ng/l]	Vitamin B ₁₂ measured [ng/l]	Recovery Rate [%]
B	481.3	187.5	668.8	681.45	107
		375.0	856.3	929.50	120
Recovery rate in total [%]					114

Sample (n=5)	Mean value measured in original sample [ng/l]	Spike [ng/l]	Vitamin B ₁₂ expected [ng/l]	Vitamin B ₁₂ measured [ng/l]	Recovery Rate [%]
C	526.44	187.5	713.94	762.23	126
		375.0	901.44	845.02	85
Recovery rate in total [%]					105

12. PRECAUTIONS

- As the test is based on a microbiological method, the general guidelines for sterile work should be observed as far as possible, (preferably work in a sterile bench / PCR hood, use of sterile instruments or equipment).
- GLP (Good Laboratory Practice) guidelines have to be observed.
- Water quality is extremely important for the test. Only the water delivered with the test kit [DIL] should be used.
- For sterile filtration, only a sterile polyethersulfone filter must be used.
- It is essential to run a standard curve for each separate assay.
- Controls should be measured with each assay.
- We recommend measurements in duplicate.
- If a higher dilution results in a higher measured value, inhibitors like antibiotics might be present.

- Wear gloves during the test.
- Avoid foaming when mixing reagents.
- Used microtiter stripes [PLATE] and materials that have been in contact with patient's samples should be handled and disposed as potentially infectious.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.

13. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- ID-Vit® is a trademark of Immundiagnostik AG.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Control samples should be analysed with each run.
- The assay should always be performed according to the enclosed manual.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

14. REFERENCES

1. Strohecker R, Henning H (1963) Vitamin-Bestimmungen. Hrsg. E. Merck AG, Darmstadt. *Verlag Chemie*, Weinheim
2. Kräutler, B., 2005. Vitamin B12: chemistry and biochemistry. *Biochemical Society transactions*, **33**, pp.806–810.
3. Obeid, R. & Herrmann, W., 2006. Mechanisms of homocysteine neurotoxicity in neurodegenerative diseases with special reference to dementia. *FEBS letters*, **580**(13), pp.2994–3005.
4. Partearroyo, T. et al., 2013. Vitamin B12 and Folic Acid Imbalance Modifies NK Cytotoxicity, Lymphocytes B and Lymphoproliferation in Aged Rats. *Nutrients*, **5**(12), pp.4836–4848.
5. Schwarz, J. et al., 2014. The influence of a whole food vegan diet with Nori algae and wild mushrooms on selected blood parameters. *Clinical laboratory*, **60**(12), pp.2039–2050.
6. Schwarz, J. et al., 2015. Biochemical Identification of Vitamin B12 Deficiency in a Medical Office. *Clinical laboratory*, **61**(7), pp.687–692.

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com