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

*For the determination of 1,25-(OH)<sub>2</sub> vitamin D<sub>3</sub>/D<sub>2</sub>  
in plasma and serum*

Valid from 2023-11-08

**KM1022**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED PURPOSE</b>	<b>20</b>
<b>2. INTRODUCTION</b>	<b>20</b>
<b>3. MATERIAL PROVIDED</b>	<b>22</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>23</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>23</b>
<i>Storage</i>	23
<i>Preparation of the mobile phases and test reagents</i>	24
<i>Preparation of the internal standard</i>	24
<i>Preparation of the calibrators and controls</i>	24
<b>6. SAMPLE PREPARATION</b>	<b>25</b>
<b>7. LC-MS/MS METHOD</b>	<b>26</b>
<b>8. CALCULATION</b>	<b>26</b>
<b>9. EXAMPLES OF CHROMATOGRAMS</b>	<b>26</b>
<i>Calibrator</i>	27
<i>Patient sample</i>	28
<i>Blank sample</i>	29
<b>10. QUALITY CONTROL</b>	<b>30</b>
<i>Normal range (plasma or serum)</i>	30
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>30</b>
<i>Trueness and linearity</i>	30
<i>Precision</i>	30
<i>Lower limit of quantification (LLOQ)</i>	31
<b>12. PRECAUTIONS</b>	<b>31</b>
<b>13. DISPOSAL</b>	<b>32</b>
<b>14. TECHNICAL HINTS</b>	<b>32</b>
<b>15. GENERAL NOTES ON THE TEST</b>	<b>32</b>
<b>16. SYMBOL EXPLANATION</b>	<b>34</b>

## 1. INTENDED PURPOSE

1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub>/D<sub>2</sub> IDKenhance® ImmuPlate® DIY LC-MS/MS kit is for the LC-MS/MS based quantification of 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub> and 1,25-(OH)<sub>2</sub>-vitamin D<sub>2</sub> in serum and plasma after immunoaffinity enrichment (ImmuPlate®) and derivatization (IDKenhance®). The assay is an *in vitro* diagnostic medical device for manual use by professional laboratory staff in a laboratory environment. It is used for the determination of the vitamin D metabolite composition in patients with suspected impaired vitamin D metabolism. It can thus be used as a support in the differentiated assessment of the vitamin D status, also in the context of a therapy of vitamin D deficiency by supplementation of vitamin D.

## 2. INTRODUCTION

Vitamin D is either produced in the skin (under the influence of UV light) or taken up from nourishment. The storage type of vitamin D, namely 25-OH-vitamin D, is formed in the liver. The hormone 1,25-(OH)<sub>2</sub>-vitamin D (D hormone) is formed in a second hydroxylation step in the kidney. The responsible enzyme, the kidney 1 $\alpha$ -hydroxylase, is subjected to a rigid control through hormones (especially parathyroid hormone) and its activity is influenced by the serum concentrations of calcium and phosphate.

The serum concentration of 1,25-(OH)<sub>2</sub>-vitamin D normally re-adjusts itself to the demands of metabolism. Deviations from the normal range of 1,25-(OH)<sub>2</sub>-vitamin D must therefore always be interpreted in the context of the remaining parameters of the calcium metabolism. The serum concentration of 1,25-(OH)<sub>2</sub>-vitamin D decreases only in seldom cases of vitamin D deficiency. For the diagnosis of vitamin D deficiency the precursor metabolite, 25-OH-vitamin D, should be measured.

The reason for a deficiency of 1,25-(OH)<sub>2</sub>-vitamin D can be found in metabolic disturbances, caused either by genetic defects of the enzyme 1 $\alpha$ -hydroxylase (rare) or kidney malfunctions (more common). Even a slightly impaired kidney function can lead to a decrease of the 1,25-(OH)<sub>2</sub>-vitamin D concentration.

Since 1,25-(OH)<sub>2</sub>-vitamin D has important functions in calcium metabolism as well as supplementing secretion of parathyroid hormone from the parathyroid glands, increasing kidney malfunctioning leads to development of renal osteopathy, which is characterised by osteomalacia and osteitis fibrosa.

Treatment of renal osteopathy consists of the administration of 1,25-(OH)<sub>2</sub>-vitamin D (calcitriol) or the prohormone 1α-hydroxy vitamin D. In renal tubules malfunctions decreased or relatively low levels of 1,25-(OH)<sub>2</sub>-vitamin D (e.g. diabetes insipidus, Fanconi syndrome) are found. A non-physiological over-production of 1,25-(OH)<sub>2</sub>-vitamin D arises in granulomatosis (e.g. sarcoidosis), where extra-renal synthesis of 1,25-(OH)<sub>2</sub>-vitamin D occurs. This can lead to hypercalcaemia. Also in idiopathic hypercalciuria a relatively high level of 1,25-(OH)<sub>2</sub>-vitamin D is found. Increased concentrations of 1,25-(OH)<sub>2</sub>-vitamin D can be seen in case of non-functional vitamin D receptors (rare), during calcium deficient nutrition, as well as a result from overproduction of parathyroid hormone (primary hyperthyroidism).

Supplemental vitamin D is available in two distinct forms: ergocalciferol (vitamin D<sub>2</sub>) and cholecalciferol (vitamin D<sub>3</sub>). Pharmacopoeias have officially regarded these two forms as equivalent and interchangeable, based on studies of rickets prevention in infants. The determination of 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>/D<sub>2</sub> as a measure of 1,25-(OH)<sub>2</sub>-vitamin D status provides an objective, quantitative measure of the biological response to vitamin D administration.

### Indications

- Defect of kidney functions
  - Chronic kidney failure
  - Haemodialysis following kidney transplantation
- Renal osteopathy
- Osteomalacia from various types of vitamin D metabolism disturbances
- Kidney tubules function disturbances (diabetes insipidus, Fanconi-Syndrom)
- Monitoring of therapy with active vitamin D metabolites
- Ideopathic hypercalciuria
- Hypercalcaemia

### 3. MATERIAL PROVIDED

Cat. No.	Identifier	Kit components	Quantity
KM0002	RECSOL	Reconstitution solution	2 x 15 ml
KM0003	WASHSOL	Wash solution	2 x 80 ml
KM0005	RMMAT	Sealing mat for ELUPLATE	1 piece
KM0006	IPMAT	Sealing mat for PLATE	1 piece
KM1000	INTSTD	Internal standard	2 x 600 µl
	MOPHAB	Mobile phase B	1 x 500 ml
	SOLA	Solution A	1 x 25 ml
KM1002	DER	Derivatization reagent	2 x 3 ml
	CAL	Calibrators 1–6; lyophilised (for concentration, see products pecification)	2 vials (à 600 µl) per level
	CTRL	Controls 1–3; lyophilised (for concentration, see product specification)	3 vials (à 600 µl) per level
	STOP	Stop solution	2 x 1.5 ml
	ACTREAG	Activation reagent	1 x 2.5 ml
	MOPHAA	Mobile phase A	1 x 500 ml
KM1100	ELUSOL	Elution solution	2 x 15 ml
KM1110	ELUPLATE	Elution plate with glass tubes	1 piece
KM1120	PLATE	96-well filter plate	1 piece
	SLU	Antibody slurry for extraction of 1,25-(OH) <sub>2</sub> -Vitamin D <sub>3</sub> /D <sub>2</sub>	1 x 10.5 ml

For reorders of single components, please use the catalogue number followed by the identifier without space as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

The following material is required for the 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub>/D<sub>2</sub> IDKenhance® ImmuPlate® DIY LC-MS/MS application (not included in the kit):

- Precision pipettors and disposable tips to deliver 10–1 000 µl
- Centrifuge capable of 550 g for 96-well plate
- Vortex mixer
- Vacuum centrifuge
- Overhead shaker or orbital shaker
- LC-MS/MS system
- LC-MS vials or 96-well plate suitable for the autosampler (if ELUPLATE is not compatible with the autosampler)

The following accessories, among others, can be ordered separately for the 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub>/D<sub>2</sub> IDKenhance® ImmuPlate® LC-MS/MS application at Immundiagnostik AG:

- tuning solution for 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>/D<sub>2</sub> (KM1002DERTU)
- tuning solution for the internal standard (KM1002DERTS)
- 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub>/D<sub>2</sub> IDKenhance® ImmuPlate® LC-MS/MS extraction kit (KM1122)
- UPLC column (KM1002SP)
- in-line filter (KM1000IF)
- in-line filter holder (KM1000IH)

Please ask for our single component price list.

Immundiagnostik AG offers further options for sample preparation for the 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>/D<sub>2</sub> LC-MS/MS application. Please contact us for your individual solution.

#### 5. PREPARATION AND STORAGE OF REAGENTS

##### *Storage*

The test reagents should be stored protected from light, dry and their specified storage temperature (CAL1–6, CTRL1–3, INTSTD, DER: -20 °C; all others 2–8 °C). The test reagents stored in this way are usable until the indicated expiry date.

**Note:** After preparation of the test reagents for the test procedure other stabilities might apply (see respective preparation step).

### *Preparation of the mobile phases and test reagents*

Before use, the mobile phases (MOPHAA and MOPHAB) and solution A (SOLA) must be activated by adding activation reagent (ACTREAG) according to the following chart:

Components			ACTREAG
Name	[ml]		[µl]
Mobile phases (MOPHAA and MOPHAB)	500	+	500
Solution A (SOLA)	25		25

Prior use mobile phases should be degassed.

**Note:** After activation with activation reagent (ACTREAG), the components mobile phase A (MOPHAA), mobile phase B (MOPHAB) and solution A (SOLA) can be stored up to 2 weeks. It is therefore recommended to prepare only as much as is needed for the test approach.

**Attention:** The activation reagent (ACTREAG) must be added under the fume cupboard. All vessels to be used must be absolutely clean, free of detergents and preferably made of LC-MS/MS suitable glass.

### *Preparation of the internal standard*

Immediately before use, the concentrate of the internal standard (INTSTD) must be diluted with the reconstitution solution (RECSOL) in a ratio of 1:5 (e.g., 0.6 ml INTSTD with 2.4 ml RECSOL). It is recommended to prepare only as much as needed for the test preparation (50 µl per sample).

### *Preparation of the calibrators and controls*

Dissolve calibrators (CAL1–6) and controls (CTRL1–3) in 600 µl of reconstitution solution (RECSOL) each while 30 s vortexing.

The concentration of the calibrators (CAL1–6) and controls (CTRL1–3) changes slightly from batch to batch; the exact content is given in the product specification.

## 6. SAMPLE PREPARATION

Serum and plasma samples are suited for the assay.

The samples must be centrifuged before use (minimum 5 min at 10 000 *g*).

Control samples should be analysed with each run.

Prior to use in the assay, allow all samples and reagents to come to room temperature (18–26 °C).

Mix samples and reagents well before use.

1.	Cut off approximately 5 min of the outlet side of a 200 µl pipette tip and pipette 100 µl SLU into each PLATE well to be used. Swirl the bottle SLU regularly to prevent the slurry from settling.
2.	Add 500 µl CAL, CTRL or sample into the corresponding wells of the filled PLATE.
3.	Add 50 µl diluted INTSTD (INTSTD diluted with RECSOL) into each filled well of the PLATE.
4.	Seal the PLATE securely with IPMAT (use roller), invert a few times and incubate for 1 h at room temperature in an orbital shaker (at least 1 200 rpm; make sure that all wells are effectively mixed) or an overhead shaker (15–20 rpm).
5.	Put the PLATE on a 96-well plate (wash plate) and centrifuge at 550 <i>g</i> for 2 min for removal of any SLU from the IPMAT.
6.	Remove and discard the IPMAT. Leave the PLATE on wash plate, centrifuge 2 min at 550 <i>g</i> until dryness; discard the flow-through.
7.	Pipette 500 µl WASHSOL into each used well of the PLATE and centrifuge at 550 <i>g</i> for 2 min until dryness; discard the flow-through. Perform this washing step a total of three times.
8.	Place the PLATE on the ELUPLATE.
9.	Add 250 µl ELUSOL into each well used, centrifuge at 550 <i>g</i> for 2 min, and collect the eluate containing the analytes in the glass tubes.
10.	Evaporate the eluate under a nitrogen stream at 37 °C or in a vacuum centrifuge.

11.	Add 50 µl DER into each glass tube and incubate for 1 h at 1 200 rpm on a shaker in the dark (attention: reaction is light sensitive).
12.	Add 25 µl STOP into each glass tube and incubate for 2 min at 1 200 rpm on a shaker.
13.	Evaporate the eluate under a nitrogen stream at 37 °C or in a vacuum centrifuge.
14.	Add 175 µl of activated SOLA (SOLA activated with ACTREAG) in every used well, seal the ELUPLATE with the RMMAT (or alternatively transfer the solution into suitable LC-MS vials for the autosampler) and shake at 1 200 rpm for 1 min.
15.	Injection into the LC-MS/MS system (see application note).

## 7. LC-MS/MS METHOD

Please refer to the application note or contact [lcms@immundiagnostik.com](mailto:lcms@immundiagnostik.com) for the parameters for setting the LC-MS/MS method.

## 8. CALCULATION

Linear regression is used as model for evaluation of the results. The six calibrator concentration points are connected by a straight line (weighting factor: 1/x). This is used to calculate the concentration of the samples and controls via the peak area.

## 9. EXAMPLES OF CHROMATOGRAMS

It must be noted that the retention time and signal intensity may vary depending on the device. The chromatograms shown here were recorded by using a Waters TQ-XS LC-MS/MS system.

## Calibrator


Fig. 1: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub> for calibrator 3.


Fig. 2: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>2</sub> for calibrator 3.

Patient sample


Fig. 3: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub> for a patient sample (citrate plasma).


Fig. 4: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>2</sub> for a patient sample (citrate plasma). The signal of the analyte is highlighted in the box due to low endogenous concentration.

*Blank sample*


Fig. 5: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>3</sub> for a blank sample.


Abb. 6: Chromatogram for the internal standard (dashed line) and the analyte (continuous line) of 1,25-(OH)<sub>2</sub>-Vitamin D<sub>2</sub> for a blank sample.

## 10. QUALITY CONTROL

Control samples should be analysed with each run. The results of the control samples are used to confirm the accuracy of the method. The test results may not be valid, if one or more values of the quality control sample are outside the acceptable range (see product specification).

### *Normal range (plasma or serum)*

Healthy adults (age 20–50) :	17–53 pg/ml
Children up to 12:	ca. 40 % higher values
Pregnant women (week 8–42):	ca. 60 % higher values
Persons older than 70:	ca. 40 % lower values

Vitamin D production in the skin is highly variable and depends, among other things, on the season and time of day, latitude, age, sun protection, etc.

We recommend each laboratory to establish its own norm concentration range.

## 11. PERFORMANCE CHARACTERISTICS

### *Trueness and linearity*

Trueness and linearity was assessed in accordance with the guidelines of CLSI EP-06-Ed2 using a certified reference material.

For 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>, a recovery between 94.0 % – 98.2 % was shown within the linear range. Linearity was assessed with a slope between 0.964 and non-linearity between -3.7 % and 1.8 %.

For 1,25-(OH)<sub>2</sub>-vitamin D<sub>2</sub>, a recovery between 96.5 % – 100.3 % was shown within the linear range. Linearity was assessed with a slope between 0.973 and non-linearity between -0.9 % and 2.9 %.

The assay has been verified as linear in a range up to 500 pg/ml (correlation = 1.0; 95 % confidence interval (CI) of the slope is between 0.95 and 0.98; 95 % CI of the y-intercept brackets the origin).

### *Precision*

Precision was determined according to the guidelines of CLSI EP05-A3. For repeatability and intermediate precision, 3 levels were measured on 20 measurement days with one run in triplicates by one user. For reproducibility, 3 levels were measured in fivefold values by 3 users on 5 measurement days.

Sample		Repeatability [CV]	Intermediate precision [CV]	Reproducibility [CV]
low	1,25-(OH) <sub>2</sub> - vitamin D <sub>3</sub> (15.6 pg/ml)	1.5 %	1.7 %	2.8 %
	1,25-(OH) <sub>2</sub> - vitamin D <sub>2</sub> (16.7 pg/ml)	1.6 %	1.7 %	2.3 %
middle	1,25-(OH) <sub>2</sub> - vitamin D <sub>3</sub> (36.7 pg/ml)	1.7 %	2.3 %	2.2 %
	1,25-(OH) <sub>2</sub> - vitamin D <sub>2</sub> (36.6 pg/ml)	1.4 %	2.2 %	2.8 %
high	1,25-(OH) <sub>2</sub> - vitamin D <sub>3</sub> (89.5 pg/ml)	1.8 %	2.4 %	3.6 %
	1,25-(OH) <sub>2</sub> - vitamin D <sub>2</sub> (87.8 pg/ml)	1.8 %	2.0 %	3.6 %

### *Lower limit of quantification (LLOQ)*

For 1,25-(OH)<sub>2</sub>-vitamin D<sub>3</sub>, the lower limit of quantification (LLOQ) is 1.41 pg/ml and for 1,25-(OH)<sub>2</sub>-vitamin D<sub>2</sub>, the LLOQ is 0.85 pg/ml. This was determined in accordance with the guidelines of CLSI EP17-A2 and is based on 120 determinations with 2 kit lots on 3 measurement days and a CV target of 20% using a precision profile.

It must be noted that the LLOQ may vary depending on the device.

## 12. PRECAUTIONS

- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious. As a precaution, it is recommended that the human material used is always considered potentially infectious.

- The GHS symbols indicated on the individual components and specifications of the material safety data sheets (available on request from Immundiagnostik AG) must be noted. When working with these reagents, the legal protective precautions must be adhered to.

### 13. DISPOSAL

Mobile phases (MOPHAA, MOPHAB), solution A (SOLA), activation reagent (ACTREAG), derivatization reagent (DER), stop solution (STOP) and elution solution (ELUSOL) must be disposed as non-halogenated solvents.

The calibrators (CAL1–6) and controls (CTRL1–3) should be disposed due to their treatment as potentially infectious material in accordance with local regulations.

### 14. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiry date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Plugs and caps of different reagents should not be swapped.
- The individual components of the kit are designed for a maximum of the specified number of test runs. Any part of the components that has already been used must not be reused, but must be disposed of properly in accordance with local regulations.

### 15. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- IDKenhance® and ImmuPlate® are brands of Immundiagnostik AG.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- Please contact Immundiagnostik AG if one or more components of the kit are damaged, missing (see material supplied) or precipitates are visible in the ready-to-use solutions.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

## 16. SYMBOL EXPLANATION

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Contains plasma derivatives or human blood		Consult instructions for use
	Consult specification data sheet		Do not re-use
	Unique Device Identification		Contains material of animal origin
	Medicinal substance		Contains material of human origin

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)