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14. LITERATUR

1. Cavia-Saiz M, Muñiz Rodríguez P, Llorente Ayala B, García-González M, Coma-Del Corral MJ, García Girón C. The role of plasma IDO activity as a diagnostic marker of patients with colorectal cancer. *Mol Biol Rep*. 2014 Apr; **41**(4):2275-9.
2. Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev*. 2014 Mar; **23**(3):461-8
3. Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology*. 2013 Mar 1; **2**(3):e23428
4. Eussen SJPM, Ueland PM, Vollset SE, Nygård O, Midttun Ø, Sulo G, Tell GS: Kynurenines as predictors of acute coronary events in the Hordaland Health Study. *Int J Cardiol*. 2015 Jun 15; **189**:18-24
5. Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis*. 2012 Jul; **18**(7):1214-20.
6. Murray C, Griffin ÉW, O'Loughlin E, Lyons A, Sherwin E, Ahmed S, Stevenson N J, et al. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity*. 2015, **48**: 274-286.
7. Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Igland J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J*. 2013 Sep; **34**(34):2689-96.

8. Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc.* 2014; **3**:e001094
9. Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol.* 2013 Sep 30; **168**(2):1435-40
10. Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol.* 2012 March; **19**(3): 436–442
11. Yan EB, Frugier T, Lim C K, Heng B, Sundaram G, Tan M, Rosenfeld J V, et al. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation.* 2015, **12**(1), 110.

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Manual

IDK® Kynurenine ELISA high sensitive

For the determination of L-kynurenine in rodents (plasma, serum and brain tissue), in cell culture supernatant and CSF

For research use only

Valid from 2021-07-14

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel: + 49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. MATERIAL SUPPLIED	17
3. MATERIAL REQUIRED BUT NOT SUPPLIED	18
4. STORAGE AND PREPARATION OF REAGENTS	18
5. STORAGE AND PREPARATION OF SAMPLES	19
6. ASSAY PROCEDURE	20
<i>Principle of the test</i>	20
<i>Derivatisation procedure</i>	20
<i>Test procedure</i>	21
7. RESULTS	22
8. LIMITATIONS	24
9. QUALITY CONTROL	24
<i>Reference Range</i>	24
10. PERFORMANCE CHARACTERISTICS	25
<i>Precision and reproducibility</i>	25
<i>Spiking recovery</i>	25
<i>Dilution recovery</i>	25
<i>Analytical sensitivity</i>	26
<i>Specificity</i>	26
11. PRECAUTIONS	26
12. TECHNICAL HINTS	26
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	26
14. REFERENCES	27

1. INTENDED USE

This Immundiagnostik AG assay is intended for the quantitative determination of L-kynurenine in EDTA plasma, serum and brain tissue of rodents (mouse, rat), and in cell culture supernatant and CSF. For research use only. Not for use in diagnostic procedures.

For human plasma, serum, and urine samples we recommend our IDK® Kynurenine ELISA K 7728.

2. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
KR3728	PLATE	Microtiter plate, pre-coated	12 x 8 wells
KR3728	STD	Standards, ready-to-use (0, 0.04, 0.16, 0.64, 2.5, 10 µM)	6 x 200 µl
KR3728	CTRL 1	Control, ready-to-use (see specification for range)	1 x 200 µl
KR3728	CTRL 2	Control, ready-to-use (see specification for range)	1 x 200 µl
KR0001.C.100	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
KR3728	AB	L-kynurenine antibody, lyophilised	2 x 1 vial
KR3728	CONJ	Conjugate, ready-to-use	1 x 12 ml
KR3728	REABUF	Reaction buffer, ready-to-use	1 x 45 ml
KR3728	DER	Derivatisation reagent	2 x 25 mg
KR0008.04	DMSO	Dimethylsulfoxide (DMSO)	1 x 4 ml
KR0011.28	ASYBUF	Assay buffer, ready-to-use	1 x 28 ml
KR0002.15	SUB	TMB substrate (tetramethylbenzidine), ready-to-use	1 x 15 ml
KR0003.15	STOP	Stop solution, ready-to-use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipets and 10-1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 14000 g
- Laboratory balance
- Ultrasonic homogeniser
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 6)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

4. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2-8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20 °C**, thaw before use in the test and mix well. Re-freeze immediately after use. The standards and controls are stable at -20 °C until the expiry date stated on the label.
- **DMSO** crystallises at 2-8 °C. Dissolve the crystals at room temperature or in a water bath.
- The **lyophilised derivatisation reagent (DER)** is stable at **2-8 °C** until the expiry date stated on the label. Reconstitute one vial of DER (25 mg) with **1.5 ml DMSO**.

Allow to dissolve for 10 minutes and mix thoroughly with a vortex-mixer. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. Please note: DMSO attacks all plastics but not polypropylene products and laboratory glass.

- The **lyophilised L-kynurenine antibody (AB)** is stable at **2-8 °C** until the expiry date stated on the label. Reconstitute one vial of AB with **3 ml of wash buffer. L-kynurenine antibody** (reconstituted AB) **can be stored at 2-8 °C for 1 month.**
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2-8 °C.**

5. STORAGE AND PREPARATION OF SAMPLES

EDTA plasma, serum, cell culture supernatant, CSF

- EDTA-Plasma and serum of rodents (mouse, rat), cell culture supernatant and CSF (cerebrospinal fluid) is suited for this test system. The samples are stable for 3 days at 2-8 °C. For longer storage, samples must be frozen at -20 °C or -80 °C.
- Samples with visible amounts of precipitates should be centrifuged.
- Samples are used **undiluted.**

If the sample volume of serum or plasma of rodents is less than 20 µl, dilute these samples 1:2 with reaction buffer (10 µl sample + 10 µl REABUF). This dilution factor must be considered in data evaluation.

Brain tissue

- Brain tissue of rodents (mouse, rat) can be tested in this assay. Place the tissue sample in a microcentrifuge tube and determine the exact wet mass. Add 40 µl of reaction buffer (REABUF) per milligram of tissue to the tube, e.g. **25 mg tissue + 1000 µl REABUF.** Homogenise the sample by sonication and centrifuge at 14,000 *g* for 10 min at 4 °C. Transfer the supernatant to a clean tube. 200 µl of the supernatant are used in the test (see sample preparation procedure). Store samples and supernatants at -20 °C or -80 °C.

For sample preparation, a derivatisation reagent for derivatisation of L-kynurenine is added to all samples (details are given in the derivatisation procedure).

6. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatisation reagent for L-kynurenine derivatisation. Afterwards, the derivatised samples and a polyclonal L-kynurenine-antiserum are incubated in wells of a microtiter plate coated with L-kynurenine-derivative (tracer).

During the incubation period, the target L-kynurenine in the sample competes with the tracer, immobilised on the wall of the microtiter wells, for the binding of the polyclonal antibodies. The L-kynurenine in the sample displaces the antibodies out of the binding to the tracer.

During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-kynurenine antibodies. After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The colour changes from blue to yellow and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow colour is inverse proportional to the L-kynurenine concentration in the sample; this means, high L-kynurenine concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standards. L-kynurenine, present in the samples, is determined directly from this curve.

Derivatisation procedure

Bring **all reagents and samples to room temperature** (15-30 °C) and mix well.

Derivatisation of standards, controls and samples is carried out in vials (e.g. 1.5 ml vials).

The reagents provided with this kit are sufficient for up to 48 derivatisations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1.	<p><i>Standards/controls, plasma, serum, cell culture supernatant, CSF:</i> Add 200 µl of standards (STD), controls (CTRL), samples in the corresponding vials.</p> <p>Add 200 µl reaction buffer (REABUF) into each vial (STD, CTRL, sample).</p> <p><i>Brain tissue:</i> Add 200 µl of supernatant after centrifugation into the corresponding vials.</p>
----	--

2.	Add 50 µl of freshly prepared derivatisation reagent into each vial (STD, CTRL, sample) and mix thoroughly by repeated inversion or several seconds on a vortex mixer.
3.	Incubate for 45 min at room temperature (15-30 °C) on a horizontal shaker .
4.	Add 500 µl assay buffer (ASYBUF) into each vial, mix thoroughly by repeated inversion or on a vortex mixer.
5.	Incubate for 5 min at room temperature (15-30 °C) on a horizontal shaker .

2 x 50 µl of the derivatised standards, controls and samples are used in the ELISA as duplicates.

Test procedure

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips covered at 2-8 °C. Strips are stable until expiry date stated on the label.

6.	Before use , wash the wells 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
7.	For the analysis in duplicate take 2 x 50 µl of the derivatised standards/controls/samples out of the vials and add into the respective wells of the microtiter plate.
8.	Add 50 µl L-kynurenine antibody into each well.
9.	Cover the strips tightly with foil and incubate overnight at 2-8°C .
10.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
11.	Add 100 µl conjugate (CONJ) into each well.
12.	Cover the strips and incubate for 1 hour at room temperature (15-30 °C) on a horizontal shaker.

13.	Discard the content of each well and wash 5 times with 250 µl wash buffer . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
14.	Add 100 µl substrate (SUB) into each well.
15.	Incubate for 8-12 min* at room temperature (15-30 °C) in the dark .
16.	Add 100 µl stop solution (STOP) into each well and mix well.
17.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference.

* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

CSF

The concentrations can be determined directly from the standard curve in $\mu\text{mol/l}$. **No factor** is required.

EDTA plasma, serum

The concentrations of undiluted samples can be determined directly from the standard curve, **no factor** is needed. Results from 1:2 diluted samples must be multiplied by the **factor of 2**.

Cell culture supernatant


Cell culture media show matrix effects when very low levels have to be determined. We suggest you proceed as follows: High L-kynurenine levels from 1 to 10 $\mu\text{mol/l}$ in the sample can be determined with the standard curve provided with the kit. For samples with lower L-kynurenine levels we recommend preparing an own standard curve by diluting the highest standard (STD6) with the cell culture media used in your laboratory.

The concentrations can be determined directly from the standard curve. **No factor** is required.

Brain tissue

The obtained L-kynurenine levels of brain tissue samples have to be multiplied by the **factor of 3.73**.

In the following, an example of a standard curve is given; do not use it for the calculation of your results.

8. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be diluted with reaction buffer (REABUF) and re-assayed. Please consider this dilution factor when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

Analytical sensitivity × sample dilution factor to be used

Analytical sensitivity see chapter "Performance Characteristics".

9. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference Range

CSF: Yan, E.B. et al., 2015. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*, 12(1):110.

Plasma (mice): Murray, C. et al., 2015. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity*, 48:274-286.

We recommend each laboratory to establish its own reference range.

10. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 10)

Sample	L-kynurenine [$\mu\text{mol/l}$]	CV [%]
1	1.04	8.1
2	2.94	4.9

Inter-Assay (n = 10)

Sample	L-kynurenine [$\mu\text{mol/l}$]	CV [%]
1	1.07	9.3
2	2.83	11.6

Spiking recovery

Two samples were spiked with different L-kynurenine concentrations and measured in this assay. The mean recovery rate was 98.7 % (n = 2).

sample [$\mu\text{mol/l}$]	spike [$\mu\text{mol/l}$]	expected [$\mu\text{mol/l}$]	measured [$\mu\text{mol/l}$]	recovery [%]
0.98	1.5	2.48	2.25	90.1
	3.0	3.98	3.95	99.2
1.13	1.5	2.63	2.51	95.4
	3.0	4.13	4.21	102.0

Dilution recovery

Two samples were diluted with reaction buffer and measured in this assay. The mean recovery rate was 97.4 % (n = 2).

sample [$\mu\text{mol/l}$]	dilution [$\mu\text{mol/l}$]	expected [$\mu\text{mol/l}$]	measured [$\mu\text{mol/l}$]	recovery [%]
2.25	1:2	1.13	1.19	105.7
	1:4	0.56	0.55	97.2
3.95	1:2	1.97	1.85	94.0
	1:4	0.99	0.92	92.8

Analytical sensitivity

The zero-standard was measured 48 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be 0.025 [$\mu\text{mol/l}$].

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against compounds with structural similarity to L-kynurenine. The specificity is calculated in percent in relation to the L-kynurenine binding activity.

3-HK (3-hydroxy-DL-kynurenine)	< 0.5 %
L-tryptophan	< 0.08 %

11. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact.

Warning: Causes serious eye irritation

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.

- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The guidelines for laboratories should be followed.
- IDK® is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

14. REFERENCES

1. Cavia-Saiz M, Muñiz Rodríguez P, Llorente Ayala B, García-González M, Coma-Del Corral MJ, García Girón C. The role of plasma IDO activity as a diagnostic marker of patients with colorectal cancer. *Mol Biol Rep.* 2014 Apr; **41**(4):2275-9.
2. Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev.* 2014 Mar; **23**(3):461-8
3. Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology.* 2013 Mar 1; **2**(3):e23428
4. Eussen S J PM, Ueland P M, Vollset S E, Nygård O, Midttun Ø, Sulo G, Tell GS: Kynurenines as predictors of acute coronary events in the Hordaland Health Study. *Int J Cardiol.* 2015 Jun 15; **189**:18-24

5. Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis*. 2012 Jul; **18**(7):1214-20.
6. Murray C, Griffin ÉW, O'Loughlin E, Lyons A, Sherwin E, Ahmed S, Stevenson N J, et al. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity*. 2015, **48**: 274-286.
7. Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Igland J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J*. 2013 Sep; **34**(34):2689-96.
8. Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc*. 2014; **3**:e001094
9. Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol*. 2013 Sep 30; **168**(2):1435-40
10. Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol*. 2012 March; **19**(3): 436–442
11. Yan EB, Frugier T, Lim CK, Heng B, Sundaram G, Tan M, Rosenfeld JV, et al. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*. 2015, **12**(1), 110.

Used symbols:

	Temperature limitation		Catalogue Number
	For research use only		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		