# BSP ELISA

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
# BSP ELISA

*For the in vitro determination of underglycosylated  
BSP (bone sialoprotein) in serum and EDTA plasma*

Valid from 2021-06-25

**KR4222**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

|   |           |
|---|-----------|
| <b>1. INTENDED USE</b>                                  | <b>13</b> |
| <b>2. INTRODUCTION</b>                                  | <b>13</b> |
| <b>3. MATERIAL SUPPLIED</b>                             | <b>13</b> |
| <b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>14</b> |
| <b>5. STORAGE AND PREPARATION OF REAGENTS</b>           | <b>14</b> |
| <b>6. STORAGE AND PREPARATION OF SAMPLES</b>            | <b>15</b> |
| <i>Sample storage</i>                                   | 15        |
| <i>Sample dilution</i>                                  | 15        |
| <b>7. ASSAY PROCEDURE</b>                               | <b>15</b> |
| <i>Principle of the test</i>                            | 15        |
| <i>Test procedure</i>                                   | 16        |
| <b>8. RESULTS</b>                                       | <b>17</b> |
| <b>9. LIMITATIONS</b>                                   | <b>18</b> |
| <b>10. QUALITY CONTROL</b>                              | <b>18</b> |
| <i>Reference range</i>                                  | 18        |
| <b>11. PERFORMANCE CHARACTERISTICS</b>                  | <b>18</b> |
| <b>12. PRECAUTIONS</b>                                  | <b>18</b> |
| <b>13. TECHNICAL HINTS</b>                              | <b>19</b> |
| <b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>19</b> |

## 1. INTENDED USE

This Immundiagnostik AG assay is an enzyme immunoassay intended for the quantitative determination of underglycosylated bone sialoprotein (BSP) in serum and EDTA-plasma. For research use only. Not for use in diagnostic procedures.

## 2. INTRODUCTION

Bone sialoprotein (BSP) is a phosphorylated glycoprotein that comprises 12 % of the total non-collagenous proteins in bone and is thought to be involved in bone mineralisation and remodelling. The polypeptide chain of BSP has a molecular weight of 35 kDa. Native BSP has an apparent molecular weight of 70–80 kDa due to complex posttranslational modifications involving glycosylation. BSP is synthesised mainly by skeletal-associated cell types, including fibroblasts, hypertrophic chondrocytes, osteoblasts, osteocytes, osteoclasts, as well as trophoblasts.

BSP is overexpressed and underglycosylated in human primary breast and prostate cancers that metastasise to the bone, whereby BSP expression correlates with the development of microcalcifications. Underglycosylated BSP is also expressed by other tumours, such as lung cancer, thyroid cancer, multiple myeloma and neuroblastoma, that predominantly metastasise to bones.

## 3. MATERIAL SUPPLIED

| Cat. No.     | Label     | Kit components  | Quantity     |
|--------------|-----------|---|--------------|
| KR4222       | PLATE     | Microtiter plate, pre-coated                                  | 12 x 8 wells |
| KR0001.C.100 | WASHBUF   | Wash buffer concentrate, 50x                                  | 1 x 100 ml   |
| KR4222       | CONJ      | Conjugate concentrate, peroxidase-labelled                    | 1 x 150 µl   |
| KR4222       | STD       | Standards, lyophilised (see specification for concentrations) | 2 x 6 vials  |
| KR4222       | CTRL 1    | Control, lyophilised (see specification for range)            | 2 x 1 vial   |
| KR4222       | CTRL 2    | Control, lyophilised (see specification for range)            | 2 x 1 vial   |
| KR4222       | SAMPLEBUF | Sample dilution buffer, ready-to-use                          | 1 x 15 ml    |

| Cat. No.  | Label | Kit components                                 | Quantity  |
|-----------|-------|--|-----------|
| KR0002.15 | SUB   | Substrate (tetramethylbenzidine), ready-to-use | 1 x 15 ml |
| KR0003.15 | STOP  | Stop solution, ready-to-use                    | 1 x 15 ml |

For reorders of single components, use the catalogue number followed by the label as product number.

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Calibrated precision pipettors and 10–1000 µl single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles >0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

#### 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASH-BUF)** has to be diluted with ultrapure water **1:50** before use (20 ml WASH-BUF + 980 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be re-dissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2–8 °C** until the expiry date stated on the label. **Wash buffer** (1:50 diluted WASHBUF) can be stored in a closed flask at **2–8 °C for 1 month**.
- The **lyophilised standards (STD)** and **controls (CTRL)** are stable at **2–8 °C** until the expiry date stated on the label. Before use, the STD and CTRL have

to be reconstituted with **500 µl of ultrapure water**. Allow the vial content to dissolve for 10 minutes and mix thoroughly to ensure complete reconstitution. **Standards and controls** (reconstituted STD and CTRL) **can be stored at -20 °C for 3 weeks. Avoid repeated thawing and freezing.**

- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in **wash buffer** (100 µl CONJ + 10 ml wash buffer). The CONJ is stable at **2–8 °C** until the expiry date stated on the label. **Conjugate (1:101 diluted CONJ) is not stable and cannot be stored.**
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8 °C**.

## 6. STORAGE AND PREPARATION OF SAMPLES

### *Sample storage*

EDTA-plasma or serum can be stored at -20 °C.

### *Sample dilution*

EDTA plasma or serum samples must be diluted **1:5** before performing the assay, e.g. **50 µl** sample + **200 µl** dilution buffer (SAMPLEBUF), mix well.

**100 µl** of the dilution are used in the test per well.

## 7. ASSAY PROCEDURE

### *Principle of the test*

This ELISA is designed for the quantitative determination of underglycosylated bone sialoprotein in serum and EDTA plasma.

The test principle is based on interaction of the antigen in the samples or standards and with the antibody coated on the wells of the microtiter plate. A peroxidase-conjugated secondary antibody is used for detection and quantification, and tetramethylbenzidine (TMB) as a peroxidase substrate. The enzymatic reaction is terminated by an acidic stop solution. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Bone sialoprotein, present in the samples, is determined directly from this curve.

### Test procedure

Bring all **reagents and samples to room temperature** (15–30°C) and mix well.

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8°C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

|     |   |
|-----|---|
| 1.  | <b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.        |
| 2.  | Add each <b>100 µl standards/controls/diluted samples</b> into the respective wells.  |
| 3.  | Cover the strips and incubate for <b>2 h</b> at room temperature (15–30°C) <b>on a horizontal shaker*</b> .   |
| 4.  | Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper. |
| 5.  | Add <b>100 µl conjugate</b> (diluted CONJ) into each well.  |
| 6.  | Cover the strips and incubate for <b>1 h</b> at room temperature (15–30°C) <b>on a horizontal shaker*</b> .   |
| 7.  | Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper. |
| 8.  | Add <b>100 µl substrate</b> (SUB) into each well.   |
| 9.  | Incubate for <b>10–20 minutes**</b> at room temperature (15–30°C) in the dark.  |
| 10. | Add <b>50 µl stop solution</b> (STOP) and mix well.   |

- |     |   |
|-----|---|
| 11. | Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference. |
|-----|---|

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 2. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

### Serum and EDTA samples

The obtained results have to be multiplied with the **dilution factor of 5** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result with the dilution factor used.

## 9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted and re-assayed. For the following analysis, the changed dilution factor has to be taken into consideration.

## 10. QUALITY CONTROL

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

We recommend each laboratory to establish its own reference range.

## 11. PERFORMANCE CHARACTERISTICS

Will follow.

## 12. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

### **13. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- The guidelines for laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue number

For research use only

To be used with

Manufacturer

Contains sufficient for &lt;n&gt; tests

Lot number

Use by

Attention

Consult instructions for use

Consult specification data sheet