**Anleitung / Manual**

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**Manual**

# **OmicSnap® Tube (pre-filled)**

***For the collection of stool samples***

Valid from Valid from 2024-01-29

REV001

**REF** **KV10021CP**

**IVD**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

**SAFETY INFORMATION**

These accessories are to be used exclusively in accordance with the enclosed instructions for use. Important safety information for this product can be found in chapter 6.

## 1. INTENDED PURPOSE

The *OmicSnap*® Tube is an all-in-one system for manual collection of a defined amount of stool by professionals and laypersons. It allows stabilization and homogenization of the stool sample for further use in nucleic acid-based procedures in the laboratory.

## 2. MATERIAL SUPPLIED

Cat. No.	Label	Kit-Komponents	Amount
KV10021CP	<i>OmicSnap</i> ® Tube (pre-filled)	Stool sample tube filled with ceramic granulate and stabilization buffer	100 units

## 3. STORAGE AND PREPARATION OF REAGENTS

The stabilization buffer in the tube is stable until the expiry date when stored at 2–30 °C.

## 4. PREPARATIONS FOR STOOL SAMPLING

- The stool sample has to be thawed. For particularly heterogeneous samples, we recommend a mechanical homogenisation using an inoculation loop or similar device.
- Place the *OmicSnap*® tube with the **white lid** facing down. Afterwards, unscrew the **yellow lid** of the tube counterclockwise (to the left) to open the tube and remove the attached dipstick. Insert the dipstick several times into the stool sample until the specimen chamber is filled. Insert the dipstick several times into the stool sample until the specimen chamber is filled. Wipe off excess stool from the front end of the dipstick (e.g. on a piece of toilet paper).
- Stick the dipstick with the attached stool **only once** back into the *OmicSnap*® Tube. Excess stool will be stripped off by the cone shaped insert, letting 100 mg of sample into the buffer.
- To close the tube, turn the **yellow lid** clockwise (to the right) until it clicks into place and the tube can no longer be opened. Tilt the *OmicSnap*® Tube several times until the specimen chamber of the dipstick is free of stool and the stool sample is completely suspended in the buffer. Please do not shake. The suspension now contains 100 mg of stool.
- **Attention:** Repeated insertion of sample into the tube will falsify the results!
- Place the tube with the **yellow lid** facing down. Afterwards, open the **white lid** of the tube. The stool-buffer solution is now ready for use according to the corresponding manual.

## 5. TECHNICAL HINTS

- The *OmicSnap*® Tube are intended for single use only.
- Do not interchange lids from different *OmicSnap*® Tubes.
- The *OmicSnap*® Tube must not be used after the expiration date.
- The included ceramic granulate ensures a simplified cell disruption for the subsequent nucleic acid purification.

## 6. WARNING AND PRECAUTIONS

- The *OmicSnap*® tube is an accessory for an *in vitro* diagnostic device. Any *in vitro* diagnostic use of the tube in combination with a downstream sample preparation and nucleic acid assay must be evaluated for the parameter in question.
- The product is to be used exclusively according to the enclosed working instructions.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.
- The buffer contained in the *OmicSnap*® Tube can cause skin irritation as well as severe eye irritation. Protective gloves, clothing and goggles should be worn. Avoid inhalation of dust. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses if possible. Continue rinsing. If eye irritation persists: Seek medical advice/attention. **IF SWALLOWED:** Call a poison control center/doctor. **IF ON SKIN:** Wash with plenty of water. **IF INHALED:** Remove to fresh air and allow person to breathe freely. **IF YOU FEEL UNWELL:** Call a poison control center/doctor.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- This product has been produced and placed on the market in accordance with the regulation (EU) 2017/746 (IVDR).

## 7. DISPOSAL

Specimens and other potentially infectious materials must be disposed of in accordance with regulatory requirements.

## 8. CHANGE HISTORY

REV	Change	Date
001	Newly generated according to IVDR	2024-01-29

## 9. SYMBOL EXPLANATION

	Temperature limitation		Catalogue Number		Lot number
	<i>In Vitro</i> Diagnostic Medical Device		Manufacturer		Unique Device Identification
	Date of manufacture		Do not re-use		Consult instructions for use
	Contains sufficient for <n> tests		Skin and Eye Irritation		

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)