# anti-htTG IgA ELISA

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
# anti-htTG IgA ELISA

*For the in vitro determination of  
anti human tissue transglutaminase IgA antibodies  
in serum and plasma*

Valid from 2021-11-11

**REF** K 9399

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. PRINCIPLE OF THE TEST</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>17</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>17</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>18</b>
<b>6. SPECIMEN COLLECTION AND PREPARATION</b>	<b>18</b>
<b>7. ASSAY PROCEDURE</b>	<b>19</b>
<i>Test procedure</i>	19
<b>8. RESULTS</b>	<b>20</b>
<b>9. LIMITATIONS</b>	<b>21</b>
<b>10. QUALITY CONTROL</b>	<b>21</b>
<i>Reference range</i>	21
<b>11. PERFORMANCE CHARACTERISTICS</b>	<b>22</b>
<i>Accuracy – Precision</i>	22
<i>Linearity</i>	22
<i>Analytical sensitivity</i>	22
<b>12. PRECAUTIONS</b>	<b>23</b>
<b>13. TECHNICAL HINTS</b>	<b>23</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>24</b>

## 1. INTENDED USE

The Immundiagnostik assay is intended for the determination of human anti-tissue-transglutaminase IgA antibodies in serum and plasma. For *in vitro* diagnostic use only.

## 2. PRINCIPLE OF THE TEST

This enzyme immunoassay is a sandwich assay for the quantitative determination of anti-tissue transglutaminase antibodies (IgA anti-htTG).

The wells of the microtiter plate are coated with the antigen. In a first incubation step, the anti-htTG IgA antibodies are bound to the coated antigen. To remove all unbound substances, a washing step is carried out.


In a further incubation step, a peroxidase-labeled antibody is added. After another washing step to remove all unbound substances, the solid phase is incubated with the peroxidase substrate, tetramethylbenzidine (TMB). Finally, an acidic stop solution is added to terminate the enzymatic reaction, whereby the colour changes from blue to yellow.

The intensity of the yellow colour is directly proportional to the IgA anti-tissue-transglutaminase antibody concentration. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. IgA anti-tissue transglutaminase antibodies present in the samples are determined directly from this curve.

### Indications

- See our proposal for stepwise diagnosis of gluten sensitivity on page 16

### 3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
K 9399	PLATE	Microtiter plate, pre-coated	12 x 8 wells
K 0001.C.100	WASHBUF	Wash buffer concentrate, 10 x	2 x 100 ml
K 9399	STD	htTG-standard lyophilised (see specification for concentrations)	4 x 1 vial
K 9399	CTRL 1	Control, lyophilised (see specification for range)	4 x 1 vial
K 9399	CTR 2L	Control, lyophilised (see specification for range)	4 x 1 vial
K 9399	CONJ	Conjugate (rabbit anti-IgA antibody, peroxidase-labelled)	1 x 200 µl
K 0002.15	SUB	Substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 0003.15	STOP	Stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water\*
- Laboratory balance
- Calibrated precision pipettors and 10–1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker with 37°C incubator
- Multi-channel pipets or repeater pipets
- Centrifuge, 3000 g
- Vortex
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader

\* Immundiagnostik AG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

## 5. STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 ml WASHBUF + 900 ml ultrapure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37°C. The **WASHBUF** is stable at **2–8°C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2–8°C for 1 month**.
- **Preparation of the conjugate:** Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in **wash buffer** (100 µl CONJ + 10 ml wash buffer). The **CONJ** is stable at **2–8°C** until the expiry date stated on the label. **Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored**.
- The **lyophilised standards (STD)** and **controls (CTRL)** are stable at **2–8°C** until the expiry date stated on the label. **Reconstitution** details are given in the **specification data sheet**. **Standards and controls** (reconstituted STD and CTRL) **are not stable and cannot be stored**.
- The preparation of the **standard curve** (volumes and concentrations) is described in the product specification.
- All other test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at **2–8°C**.

## 6. SPECIMEN COLLECTION AND PREPARATION

### Serum and plasma

Dilute serum and plasma **1:1 000** in ELISA wash buffer in two dilution steps. For example:

- **10 µl** serum + **90 µl** wash buffer, mix well = **dilution I** (1:10)
- **10 µl** dilution I + **990 µl** wash buffer, mix well = **dilution II** (1:100)  
This results in a final dilution of 1:1000.

For analysis, pipet **100 µl** of **dilution II** per well.

## 7. ASSAY PROCEDURE

### *Test procedure*

Bring all **reagents and samples to room temperature** (15–30 °C) and mix well.

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2–8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	<b>Before use</b> , wash the wells <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2.	Add each <b>100 µl standards/controls/prepared samples</b> into the respective wells.
3.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal shaker*</b> .
4.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5.	Add <b>100 µl conjugate</b> (diluted CONJ) into each well.
6.	Cover the strips and incubate for <b>1 hour</b> at room temperature (15–30 °C) on a <b>horizontal shaker*</b> .
7.	Discard the content of each well and wash <b>5 times</b> with <b>250 µl wash buffer</b> . After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8.	Add <b>100 µl substrate</b> (SUB) into each well.
9.	Incubate for <b>10–20 min**</b> at room temperature (15–30 °C) in the <b>dark</b> .
10.	Add <b>100 µl stop solution</b> (STOP) into each well and mix well.

- |     |   |
|-----|---|
| 11. | Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm as a reference. |
|-----|---|

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the 4 parameter algorithm.

### 1. 4 parameter algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

### 3. Spline algorithm

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the paired values should be evaluated manually.

### Serum / plasma samples

The sample dilution of 1:1000 is already considered in the calibration curve.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## 9. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be further diluted and re-assayed.

## 10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

htTG IgA Serum cut off	> 22 AU/ml	positive
htTG IgA Serum cut off	< 16 AU/ml	negative

Results between 16–22 AU/ml are in the grey area.

We recommend each laboratory to establish its own reference range.

In order to avoid unnecessary biopsies, tTG-IgA titers within the grey area should be controlled after 3–6 months. Biopsies should be performed only in case of increased tTG-IgA titers. Until now, there is no information on correlation between age and titer level in children with negative tTG-IgA values. Age-normalized threshold values could be useful.\*

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### *Accuracy – Precision*

#### **Repeatability (Intra-Assay); n = 64**

The repeatability was assessed with one serum sample under **constant** parameters (same operator, instrument, day and kit lot).

Sample	Mean value [AU/ml]	CV [%]
1	83.80	4.1

#### **Reproducibility (Inter-Assay); n = 34**

The reproducibility was assessed with 2 serum samples under **varying** parameters (different operators, instruments, days and kit lots).

Sample	Mean value [AU/ml]	CV [%]
1	69.55	14.4
2	15.51	9.7

### *Linearity*

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP6-A by serial dilution of two serum samples. For anti-htTG IgA in serum and plasma, the method has been demonstrated to be linear from 5.55 to 266.67 AU/ml, showing a non-linear behaviour of less than  $\pm 20\%$  and a correlation of  $> 95\%$  in this interval.

### *Analytical sensitivity*

The zero-standard was measured 28 times. The following value was obtained without considering possibly used sample dilution factors.

Limit of blank, LoB

0.654 AU/ml

## 12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Control samples should be analyzed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The 10x Wash buffer concentrate (WASHBUF) contains surfactants which may cause severe eye irritation in case of eye contact

**Warning:** Causes serious eye irritation

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical Advice/attention.

- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch as wells from already opened microtiter plates are exposed to different conditions than sealed ones.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according to the enclosed manual.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

### Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		