MutaREAL[®] ApoE

Real-Time-PCR-Kit

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MutaREAL[®] ApoE

Real-Time-PCR Kit

*For the analysis of the amino acid substitution
Cys130Arg and Cys176Arg in the ApoE gene*

Valid from 2021-10-28

KF290732
KF290796

32/96

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

| | | |
|-----------|--|-----------|
| 1 | INTENDED USE | 13 |
| 2 | INTRODUCTION | 13 |
| 3 | PRINCIPLE OF THE TEST | 13 |
| 4 | PACKAGE CONTENTS | 14 |
| 5 | EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER | 14 |
| 6 | TRANSPORT, STORAGE AND STABILITY | 14 |
| 7 | WARNINGS AND PRECAUTIONS | 14 |
| 8 | SAMPLE MATERIAL | 15 |
| 9 | REAL-TIME-PCR | 15 |
| | 9.1 <i>Important points before starting</i> | 15 |
| | 9.2 <i>Procedure</i> | 16 |
| | 9.3 <i>Instrument settings</i> | 17 |
| 10 | DATA ANALYSIS | 18 |
| | 10.1 <i>Evaluation and classification in ApoE isoforms</i> | 19 |
| 11 | TROUBLESHOOTING | 20 |
| 12 | LIMITATIONS OF THE METHOD | 20 |
| 13 | ABBREVIATIONS AND SYMBOLS | 20 |
| 14 | LITERATURE | 21 |

1 INTENDED USE

The MutaREAL® ApoE Real-Time PCR Kit is a FRET-based molecular biology test for the investigation of point mutations in the ApoE gene, which are responsible for the exchange of the amino acids Cys130Arg (formerly pos. 112) and Cys176Arg (formerly pos. 158).

2 INTRODUCTION

ApolipoproteinE (ApoE) belongs to the family of apolipoproteins, which play an important role in lipid metabolism and influences the formation of very-low-density lipoproteins (VLDL). The most common ApoE alleles are ApoE2 (E2), ApoE3 (E3) and ApoE4 (E4). These are defined by the variations Cys130Arg and Cys176Arg. E3 carries cysteine at position 130 and arginine at 176, E2 carries cysteine at both positions, whereas E4 carries arginine at both positions. The different ApoE isoforms have an influence on cholesterol levels. [1]

3 PRINCIPLE OF THE TEST

The sequence-specific MutaREAL® ApoE real-time PCR kit is based on fluorescence resonance energy transfer (FRET). The assay includes two specific primers flanking the target sequence and two hybridisation probes that bind adjacent to the target sequence. One of the hybridisation probes is labelled with a donor fluorophore and, after appropriate excitation, transfers its energy to the acceptor fluorophore with which the other hybridisation probe is labelled when they are in close proximity. After the energy transfer, the acceptor dye emits light with a longer wavelength. Energy transfer can only occur when both hybridisation probes have bound to the target sequence. The amount of hybridised probe pairs and thus the fluorescence signal increases with the amount of amplified PCR product. Here, the fluorescence signal is proportional to the amount of PCR product.

Genotyping is carried out after completion of amplification by melting curve analysis. For this purpose, the temperature is slowly increased after a denaturation step and the dissociation behaviour of the hybridisation probes is recorded while continuously measuring the fluorescence. One of the hybridisation probes binds to a part of the target sequence that is present in the wild type and the mutation. The second hybridisation probe spans the mutation site. As the temperature rises, the mismatched and thus less stable probes dissociate first and the fluorescence decreases. The perfectly paired hybridisation probes dissociate later due to their higher binding energy and thus the fluorescence signal decreases only at a higher temperature.

4 PACKAGE CONTENTS

The components supplied are sufficient for the preparation of 32 (KF290732) or 96 (KF290796) reactions.

Table 1: Components of the MutaREAL® ApoE Real-Time-PCR Kit.

| Label | Lid Colour | Content | |
|-----------------------------|------------|------------|-------------|
| | | 32 | 96 |
| Enzyme mix | blue | 1 x 876 µl | 2 x 1313 µl |
| Detection mix 1 (Cys130Arg) | yellow | 1 x 368 µl | 3 x 368 µl |
| Detection mix 2 (Cys176Arg) | brown | 1 x 368 µl | 3 x 368 µl |
| Positive control | red | 1 x 30 µl | 1 x 90 µl |
| Negative control | green | 1 x 150 µl | 1 x 150 µl |

5 EQUIPMENT AND REAGENTS TO BE SUPPLIED BY USER

- DNA extraction kit (e.g. MutaCLEAN® Universal RNA/DNA, KG1038)
- Roche LightCycler® 1.5 or 2.0 real-time PCR system
 - * CE conformity only exists if one of the mentioned devices is used.
- Roche LightCycler® capillaries
- Roche LightCycler® Cooling Block
- Sterile reaction tubes
- Calibrated pipettes (variable volumes) and sterile disposable tips with filter
- Optional: Liquid handling system for automation

6 TRANSPORT, STORAGE AND STABILITY

The MutaREAL® ApoE real-time PCR kit is transported frozen on dry ice or cool packs. All components are to be stored protected from light at a minimum of -20 °C immediately after receipt. Avoid multiple freeze-thaw cycles (make aliquots if necessary). Do not use after the expiry date indicated on the package.

After opening, the detection mixes and the enzyme / buffer mix can be stored at +4 to +8 °C for up to 6 months.

Be sure to protect the detection mixes from direct sunlight during the entire test procedure.

7 WARNINGS AND PRECAUTIONS

Read the instructions for use carefully before using the product.

- All samples must be considered potentially infectious and/or biohazardous and all items that come into contact with the specimens must be considered potentially contaminated.
- Real-time PCR must be performed in laboratories suitable for this purpose and by specially trained personnel.
- The assay must always be carried out according to the instructions supplied with the kit.
- Areas for sample preparation and preparation of the PCR master mix should be strictly separated.
- Pipettes, tubes and other working materials must not circulate from one area to the other.
- Always use pipette tips with filters.
- Always wear powder-free disposable gloves when using the kit
- Clean pipettes and work surfaces regularly with suitable decontamination solution (no ethanol-containing agents).
- Contamination of eluates and kit components with microbes or nucleases (RNAs and DNAses) should be avoided.
- Positive and potentially positive material must be kept separate from all other kit components at all times.
- Do not open reaction tubes/plates after amplification in order to avoid contamination.
- In accordance with guidelines or requirements of local, state or federal regulations or authorised organisations, additional controls may be tested.
- Do not autoclave reaction tubes after PCR as this will not degrade the amplified nucleic acid and risks contaminating the laboratory area.
- Dispose of samples and test waste according to your local safety regulations.
- Refrigerate all PCR reagents while working.
- The purity (A260/A280) of the genomic DNA should be between 1.8 and 2.0.

8 SAMPLE MATERIAL

Starting material for the MutaREAL® ApoE real-time PCR kit is genomic DNA isolated from clinical samples (blood) using a suitable extraction kit.

9 REAL-TIME-PCR

9.1 *Important points before starting*

- Please pay attention to chapter 7 “Warnings and precautions”.

- Before setting up the Real-Time-PCR familiarise yourself with the Real-Time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run all Positive Controls and one Negative Control should be included.
- Before each use, all reagents must be gently thawed, thoroughly mixed (do not vortex) and briefly centrifuged.
- Protect the detection mixes from exposure to light.
- We recommend always cooling the reagents and the preparation in a cooling block (+4 to +8 °C) or on ice while working.

9.2 Procedure

For amplification, two reaction tubes (LightCycler® capillaries) are needed per sample plus two additional reaction tubes per master mix for the negative and the positive controls. The following table shows the volumes to be pipetted per sample. For the analysis it is recommended to prepare a master mix for the number of samples (incl. negative and positive controls) (N) plus 10 % to compensate for inaccuracies. The master mix is pipetted as described in tables 3 and 4:

Master mix 1 (Cys130Arg)

Table 2: Preparation of master mix 1 (Cys130Arg)

| Reagent | Volume per 25 µl - reaction mix | Master mix volume |
|-----------------------------------|---------------------------------|---------------------------|
| Detection mix 1 (yellow) | 10.5 µl | 10.5 µl * (N + (N * 0.1)) |
| Enzyme mix (blue) | 12.5 µl | 12.5 µl * (N + (N * 0.1)) |

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

Master mix 2 (Cys176Arg)

Table 3 Preparation of master mix 2 (Cys176Arg)

| Reagenz | Volume per 25 µl - reaction mix | Master mix volume |
|----------------------------------|---------------------------------|---------------------------|
| Detection mix 2 (brown) | 10.5 µl | 10.5 µl * (N + (N * 0.1)) |
| Enzyme mix (blue) | 12.5 µl | 12.5 µl * (N + (N * 0.1)) |

- Mix the Master Mix carefully by pipetting up and down or by inverting and centrifuge briefly.
- Add **23 µl** of the Master Mix to each capillary.
- For the negative control add **2 µl** of the supplied negative control (**green**).
- For the positive control add **2 µl** of the supplied positive control (**red**).
- For each sample to be analysed, add **2 µl** of the extracted genomic DNA to the corresponding capillary.

Close the capillaries with the lids, transfer them into the LightCycler® carousel and centrifuge them in the LightCycler® centrifuge (if a tabletop centrifuge is used, centrifuge the capillaries in the inserts of the cooling block at 3000 rpm for 15 s). Then transfer the carousel to the LightCycler® and start the PCR programme described in 9.3.

9.3 Instrument settings

For the Real-Time-PCR use the thermal profile shown in table 5.

Table 4: Real-Time-PCR thermal profile

| Description | Time | Temperature | Heating rate | Cycles | Acquisition |
|----------------------|-------|-------------|--------------|--------|-----------------|
| Initial Denaturation | 120 s | 94 °C | 20 °C/s | 1 | none |
| Denaturation | 10 s | 94 °C | 20 °C/s | 45 | none |
| Primer annealing | 25 s | 60 °C | 20 °C/s | | single |
| Elongation | 25 s | 72 °C | 20 °C/s | | none |
| Melting curve | 20 s | 94 °C | 20 °C/s | 1 | none |
| | 20 s | 40 °C | 20 °C/s | 1 | none |
| | 0 s | 80 °C | 0.2 °C/s | 1 | constant |
| Cooling | 30 s | 40 °C | 20 °C/s | 1 | - |

10 DATA ANALYSIS

Add an analysis of the type „genotyping“ for the evaluation of the melting curves. This forms the derivation of the fluorescence curve. The detection wavelength is between 670 nm and 705 nm.

ApoE Cys130Arg

Temperature wild type allele (Cys): 51.5 °C (+/-2 °C)

Temperature mutation allele (Arg): 61.0 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **green curve** - negative control, **red curve** - homozygous Cys / Cys, **blue curve** - heterozygous Cys / Arg.


Fig. 1: Evaluation of master mix 1 - Cys130Arg.

ApoE Cys176Arg

Temperature wild type allele (Cys): 52.5 °C (+/-2 °C)

Temperature mutation allele (Arg): 61.5 °C (+/-2 °C)

The following graph shows the typical results for the possible genotypes: **pink curve** - negative control, **black curve** - heterozygous Cys / Arg, **blue curve** - homozygous Arg / Arg.


Fig. 2: Evaluation of master mix 2 - Cys176Arg.

The provided positive control contains a template which is **heterozygous** for both mutations (ApoE E2/E4).

10.1 Evaluation and classification in ApoE isoforms

Table 5: Evaluation and classification in ApoE isoforms

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11 TROUBLESHOOTING

The following troubleshooting guide is included to help you with possible problems that may arise when performing a Real-Time-PCR. If you have further questions, please do not hesitate to contact our scientists on info@immundiagnostik.com.

No or weak fluorescence in the positive control or samples.

Check the PCR programme of the real-time PCR system and repeat the analysis with the corrected protocol.

Detection mixes have been subjected to more than two freeze cycles or have been stored at 2-8 °C for more than four days. Repeat the analysis with a fresh aliquot or new detection mix.

The quality of the starting DNA is not sufficient. Use freshly extracted DNA and determine the concentration/purity before use.

The detection mixes were not protected from light exposure. Repeat the analysis with a fresh aliquot or new PCR reagents.

12 LIMITATIONS OF THE METHOD


The result is provided to the treating physician as supporting material and should never be used exclusively for diagnosis or treatment recommendations. The diagnosis as well as the treatment decisions to be taken remain the full responsibility of the physician.

The accuracy of genetic tests is not 100%. However, it has been found to be over 98% accurate based on validation data for this test. Furthermore, genetic test results must be considered in the context of the patient's clinical representation and known familial risks in the patient's environment.

The test only analyses a selection of markers. Therefore, a negative test result of the patient does not completely exclude a risk of any kind.

13 ABBREVIATIONS AND SYMBOLS

| | | | |
|---|---------------------------|---|----------------------------------|
| DNA | Deoxyribonucleic acid |  | Catalog number |
| PCR | Polymerase chain reaction |  | To be used with |
|  | Enzyme mix |  | Contains sufficient for <n> test |

| | | | |
|---|---|---|------------------------------|
|  | Detection mix 1 |  | Upper limit of temperature |
|  | Detection mix 2 |  | Manufacturer |
|  | Positive control |  | Lot number |
|  | Negative control |  | Consult instructions for use |
|  | <i>In vitro</i> diagnostic medical device |  | Content |
|  | Use by YYYY-MM-DD | | |

14 LITERATURE

1. Eichner et al., Am J Epidemiol, 2002, 155(6):487-95

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com