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
ID-Vit[®] **Folic acid in whole blood**

***Microbiological test kit for the determination of folic acid
in whole blood using a Lactobacillus rhamnosus-coated
microtitre plate
For use in human and veterinary medicine and in research***

Valid from 2022-04-22

KIF005VB

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: info@immundiagnostik.com www.immundiagnostik.com

Table of Contents

1. INTENDED USE	19
2. INTRODUCTION	19
3. MATERIAL SUPPLIED	21
4. MATERIAL REQUIRED BUT NOT SUPPLIED	21
5. PRINCIPLE OF THE TEST	22
6. STORAGE AND PREPARATION OF REAGENTS	22
6.1 Water	22
6.2 Preparation of the lysis-reagent	22
6.3 Preparation of the controls	23
6.4 Preparation of the standard curve	23
6.5 Preparation of the sterile assay medium	23
6.6 Microtiter plate [PLATE]	24
7. SAMPLE STORAGE AND PREPARATION	24
7.1 Sample pretreatment	24
7.2 Sample dilution	25
8. ASSAY PROCEDURE	25
8.1 Test preparations	25
8.2 Test procedure	25
8.3 Measurement	25
9. EVALUATION OF RESULTS	26
9.1 Calculation	26
9.2 Quality control	26
10. LIMITATIONS	27
11. PERFORMANCE CHARACTERISTICS	27
11.1 Precision and reproducibility	27
11.2 Recovery	27
11.3 Linearity	29
12. PRECAUTIONS	30
13. GENERAL NOTES ON THE TEST	31
14. REFERENCES	31

1. INTENDED USE

ID-Vit® Folic acid is a microtiter plate test kit based on a microbiological method for determining the total content of folic acid from whole blood. All required reagents, as well as the standard are included in the test. With the test 96 determinations including the standards. An ELISA reader is required for the evaluation. For use in human and veterinary medicine and in research. For *in vitro* diagnostic use only.

2. INTRODUCTION

Folic acid (vitamin B9) belongs to the group of the vitamin B complex; it is water-soluble, sensitive to light, oxygen, increased temperatures and is essentially involved in all growth and development processes of the body. Folic acid is necessary for the formation of red blood cells, for optimal bone marrow function and healthy nerve activity. Folic acid is also essential for all cell divisions (hence its importance in foetal development).

Although folic acid is found in most plant and animal foods an undersupply of folic acid is the most common vitamin deficiency in Europe and North America. According to the German Nutrition Society only one in four Germans consumes enough folic acid - often the result of an unbalanced diet with little fresh fruit and vegetables. But age, various diseases and the intake of certain medicines such as cotrimoxazole can also lead to absorption disorders and an associated deficiency. Lowered folic acid levels occur because of:

- a decreased supply (e.g. folic acid antagonists).
- a disrupted resorption (e.g. in celiac disease. CED).
- an increased requirement (e.g. increased alcohol consumption. when taking contraceptives. during pregnancy. anaemia or cancer).

Symptoms of Deficiency

The early symptoms of deficiency are fatigue irritability, depression, poor concentration and loss of appetite; other consequences are inflammation of the mucous membranes, anaemia and severe neurological damage. During pregnancy, when the need for folic acid doubles, a deficiency can lead to premature birth and severe malformations. Optimal folic acid supply during pregnancy can reduce the risk of neural tube defects in the foetus by 85%. Since both a folic acid deficiency and a vitamin B12 deficiency can cause megaloblastic anaemia, the determination of both vitamins is particularly important in this clinical picture in order to supplement the right vitamin. Treatment of megaloblastic anaemia caused by vitamin B12 deficiency with folic acid alone can lead to irreversible damage to the central nervous system.

Folic acid and arteriosclerosis

A folic acid deficiency is known to be the most common cause of hyperhomocysteinaemia. Meanwhile, the hyperhomocysteinaemia has been recognised as an independent factor in arteriosclerosis. Therefore, the determination of folic acid can be carried out within the framework of a coronary disease risk analysis. Beside of the influence of folic acid on the homocysteine levels a further positive effect on the endothelial function in heart patients has been established – development of nitrate tolerance during continuous nitrate therapy, e.g. in such patients, an increased release of oxygen radicals occurs without folic acid supplementation (Verhaar et al. 2002).

Indications

- Hyperchrome, macrocytic anemia
- Long-term therapy with antiepileptic drugs or folic acid antagonists
- Long-term haemodialysis
- Multiple birth pregnancy / planned pregnancy
- Enhanced erythropoiesis
- Chronic liver diseases
- Hemoblastosis
- Psoriasis, dermatitis
- Stomatitis, glossitis
- Chronic alcohol abus

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity	
			KIF005VB	KIF005VB.2
KIF005VB	PLATE	microtiter plate, precoated with <i>Lactobacillus rhamnosus</i>	1 x	2 x
	SOL	Sample treatment solution, ready-to-use	5 x 4.6 ml	10 x 4.6 ml
	DIL	Water	4 x 30 ml	7 x 30 ml
	PAF	Folic acid buffer, lyophilized	5 x	10 x
	ASYMED	Folic acid assay medium, ready-to-use	4 x	4 x
	ASYBUF	Folic acid medium treatment buffer, ready-to-use	4 x 1.5 ml	4 x 1.5 ml
	STD	Folic acid standard, lyophilized	4 x	3 x
	FOL	Adhesive cover foil	1 x whole 3 x half	3 x whole
	CTRL1	Folic acid control 1, lyophilized	4 x	3 x
	CTRL2	Folic acid control 2, lyophilized	4 x	3 x
	FRA	Replacement holder for microtiter strips	1 x	1 x

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Incubator with a dark incubation chamber (37 °C)
- ELISA reader 610–630 nm (540–550 nm)
- Calibrated precision pipettors and sterile 2–1 000 µl tips
- 1.5–2 ml reaction vials, sterile
- 0.2 µm sterile polyethersulfone (PES) filter with a sterile disposable syringe (10 ml)

- 15 ml centrifuge tubes, sterile (e.g. Falcon tubes)
- Centrifuge (10000g)

5. PRINCIPLE OF THE TEST

The EDTA whole blood samples are pre-treated and diluted with a buffer mixture and then transferred into the wells of a microtiter plate precoated with *Lactobacillus rhamnosus*. The addition of folic acid in either standards or samples gives a folic acid-dependent growth response until folic acid is consumed. After incubation at **37°C** for **48 h** the growth of *Lactobacillus rhamnosus* is measured turbidimetrically at 610–630 nm (alternatively at 540–550 nm) in an ELISA reader and compared to a standard curve generated from the dilution series. The amount of folic acid is directly proportional to the turbidity.

6. STORAGE AND PREPARATION OF REAGENTS

- Store test kit and reagents at 2–8°C.
- Prepare reagents freshly and use them immediately after preparation. Discard remaining unused reagents and waste in accordance with country, federal, state and local regulations.
- To run the assay more than once ensure that reagents are stored at the conditions stated on the label. Prepare only the appropriate amount necessary for each run. The kit can be used up to 4 times within the expiry date stated on the label.

6.1 Water

- Water [DIL] (for medium [ASYMED], standard [STD] and controls [CTRL1, CTRL2]).
- Push the lid up and pull it back to the rim of the glass, then twist the whole cap off.

6.2 Preparation of the lysis-reagent

- Transfer 4.5 ml sample preparation buffer [SOL] into vials containing lyophilised buffer [PAF], then homogenise using a vortex.
- The lysis-reagent cannot be stored.

6.3 Preparation of the controls

- Each lyophilised control [CTRL1, CTRL2] have to be resuspended with 125 µl lysis-reagent (corresponds to a 1:10 dilution), then homogenise using a vortex.
- After reconstitution, the controls are treated like samples.
- The concentration of the controls changes from batch to batch and can be found in the product specification.

6.4 Preparation of the standard curve

- For the preparation of the standard curve, standard concentrate is needed. To prepare standard concentrate, resuspend the lyophilised standard [STD] with x ml (x = please see the enclosed quality control protocol for the volume needed) water [DIL] supplied with the test kit. Then homogenise using a vortex.
- Prepare a standard curve in 6 sterile reaction tubes (1.5–2 ml volume) from standard concentrate and water [DIL] following the scheme depicted in the table below:

Folic acid [µg/l]	Water [DIL] [µl]	+	Standard concentrate [µl]	=	Total volume [µl]
Blank: 0	500	+	0	=	500
Standard 1: 0,04	450	+	50	=	500
Standard 2: 0,12	350	+	150	=	500
Standard 3: 0,20	250	+	250	=	500
Standard 4: 0,28	150	+	350	=	500
Standard 5: 0,36	50	+	450	=	500

6.5 Preparation of the sterile assay medium

- Fresh sterile assay medium has to be prepared each time before performing a test.
- Remove lyophilised assay medium from the desiccant bag in the assay medium bottle by taking the bag with a forceps and shaking it whilst still inside the bottle. Then remove the clean desiccant bag and discard it.

- Add 10ml water [DIL] and 1 ml medium treatment buffer [ASYMED] to the assay medium bottle [ASYMED], close the bottle firmly and shake it. This amount is sufficient for 6 microtiter stripes.
- Heat the medium bottle in a water bath at 90–100°C for 5 min. shake well at least 2 times during this incubation time. Take care that the medium bottle is always firmly closed.
- Quickly cool the medium bottle to < 30°C (at 2–8°C for 10 min).
- Filter the medium using a disposable syringe (10 ml) and the 0.2 µm PES filter into a sterile centrifuge tube (15 ml. e. g. Falcon).
- After this preparation the sterile assay medium can be used in the test.

6.6 Microtiter plate [PLATE]

- Store the microtiter plate [PLATE] in the aluminium packaging containing the desiccant bag at 2–8°C.
- The microtiter plate [PLATE] has to be protected from humidity and contamination.
- Take care that the aluminium packaging is not damaged.
- Carefully close the aluminium packaging after opening.
- Take only the microtiter stripes needed directly before usage to avoid contamination.

7. SAMPLE STORAGE AND PREPARATION

- Samples are stable at 2–8°C for 8 hours in the dark. For longer storage samples should be frozen and kept at -20°C.

7.1 Sample pretreatment

Add 25 µl of EDTA whole blood sample to 225 µl of the prepared lysis-reagent (ratio 1:10). Transfer the controls into a 1.5 ml reaction tube. Incubate samples and controls at 37°C for 30 min. cool rapidly (at 2–8°C for 10 min) and centrifuge (10 000 g) for 10 min.

7.2 Sample dilution

Take 10 µl from the supernatant of the pre-treated **EDTA whole blood** and add 740 µl water [DIL] and mix. The sample pretreatment and dilution corresponds in total to a **1:750 dilution** (= sample dilution factor).

Take 10 µl of the supernatant of the pre-treated **controls** and add 740 µl water [DIL] and mix. The sample pre-treatment and dilution corresponds in total to a **1:750 dilution** (= sample dilution factor).

8. ASSAY PROCEDURE

8.1 Test preparations

Remove only the reagents and materials as needed to performed the test and put the rest of the test kit back in the refrigerator. Bring all the required reagents to room temperature before start.

8.2 Test procedure

- Take as many microtiter strips as needed from the kit and place them in the second microtiter strip holder [FRA].
- Add 150 µl sterile assay medium (ASYMED) into the cavities.
- Add each 150 µl of the prepared standard solutions, samples and controls into the respective cavities. Pre-rinse each pipet tip with standard, control or sample solution, respectively.
- Carefully seal the plate with adhesive cover foil [FOL]. **Important:** the cavities must be made airtight, e.g. by pressing the foil down with the hand!
- Keep at 37 °C for 48 h in an incubator.

8.3 Measurement

- Press the adhesive cover foil [FOL] firmly down again with the hand.
- Turn the microtiter plate [PLATE] upsidedown put it onto a tabletop and shake the microbes well.
- Turn the microtiter plate [PLATE] over again and carefully remove the adhesive cover foil [FOL]. During this fix the strips in the frame with your hand because the foil is highly adhesive.

- Remove air bubbles in the cavities using a pipet tip or a needle.
- Read turbidity in an ELISA reader at E 610–630 nm (alternatively at E 540–550 nm).

Please note

- After 48 h incubation time the microtiter plate [PLATE] may be stored for a maximum of 48 h in the refrigerator before measuring the turbidity.
- To prevent time-loss through public holidays or weekends the microtiter plate [PLATE] may also be evaluated after 60 h incubation.

9. EVALUATION OF RESULTS

We recommend to use the 4 parameter algorithm to calculate the results. The sample dilution factor has to be considered for data evaluation.

The blank should have an optical density < standard 1. It serves as optical control to exclude contaminations and is not included in the calculation of results.

9.1 Calculation

Folic acid in $\mu\text{g/l}$ = value from the standard curve \times sample dilution factor (750)

Reference value for EDTA whole blood sample

Based on internal studies of EDTA whole blood samples of apparently healthy persons ($n = 74$) the following values were estimated.

Folic acid: 79.5–231.4 $\mu\text{g/l}$

Please note

With a sample dilution factor of 750 a range from 30–270 $\mu\text{g/l}$ folic acid is covered.

We recommend each laboratory to establish its own reference range as reference ranges are strongly dependent on the selection of the sample collective. The reference range for folic acid is given for guidance only and may differ from other published data.

9.2 Quality control

The extinction of the highest standard has to be > 0.6 .

Immundiagnostik AG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. If one of the controls or the highest standard are outside the specified range Immundiagnostik AG cannot guarantee the accuracy of the measurement results.

10. LIMITATIONS

Serum and plasma cannot be used in the test. The test was validated with EDTA whole blood only. If another form of whole blood is desired it must first be validated.

11. PERFORMANCE CHARACTERISTICS

The following performance characteristics have been collected using human with EDTA whole blood samples.

11.1 Precision and reproducibility

Intraassay (n = 22)

Sample	Folic acid [$\mu\text{g/l}$]	CV [%]
1	70.24	7.2
2	146.03	8.5
3	239.31	5.6

Interassay (n = 33)

Sample	Folic acid [$\mu\text{g/l}$]	CV [%]
1	84.53	7.6
2	149.08	7.7
3	227.9	8.2

11.2 Recovery

For the determination of recovery samples of 4 EDTA whole blood samples were spiked with folic acid concentrations in different concentrations and measured. The result is obtained by the ratio between the measurement result and the true concentration of a sample.

Sample (n=9)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
A	99.2	40.8	132.98	129.86	98
		76.8	168.98	186.44	110
		100.3	192.48	216.34	112
Recovery rate in total [%]					107

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
B	107.1	40.8	147.89	167.37	113
		76.8	183.89	195.19	106
		100.3	207.39	243.72	118
Recovery rate in total [%]					112

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
C	108.4	40.8	149.21	135.01	90
		76.8	185.21	168.53	91
		100.3	208.71	215.19	103
Recovery rate in total [%]					95

Sample (n=8)	Mean value original sample [µg/l]	Spike [µg/l]	Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
D	115.6	40.8	156.38	170.08	109
		76.8	192.38	211.48	110
		100.3	215.88	215.58	100
Recovery rate in total [%]					106

11.3 Linearity

EDTA whole blood samples from four patients were serially diluted and analysed by two operators. The results are shown below.

Sample	Dilution		Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
A	1:425	1	316.7	303.2	96
	1:500		271.7	259.5	95
	1:750		194.1	194.1	100
	1:1 000		145.6	153.8	106
	1:1 250		116.5	113.2	97
B	1:425	2	257.1	266.6	104
	1:500		218.6	226.5	104
	1:750		145.7	145.7	100
	1:1 000		109.3	117.8	108
	1:1 250		87.4	94.1	108
C	1:425	1	279.9	275.4	98
	1:500		237.9	219.4	92
	1:750		158.6	158.6	100
	1:1 000		119.0	120.3	101
	1:1 250		95.2	96.9	102

Sample	Dilution		Folic acid expected [µg/l]	Folic acid measured [µg/l]	Recovery Rate [%]
D	1:425	2	265.8	267.6	101
	1:500		225.9	227.9	101
	1:750		150.6	150.6	100
	1:1 000		113.0	118.1	105
	1:1 250		90.4	95.5	106

12. PRECAUTIONS

- As the test is based on a microbiological method, the general guidelines for sterile work should be observed as far as possible (preferably work in a sterile bench / PCR hood, use of sterile instruments or equipment).
- GLP (Good Laboratory Practice) guidelines have to be observed.
- Water quality is extremely important for the test. Only the water delivered with the test kit [DIL] should be used.
- For sterile filtration a **sterile polyethersulfone** filter is mandatory.
- It is essential to run a standard curve for each separate assay.
- We recommend measurements in duplicate.
- If a higher dilution results in a higher value measured, inhibitors like antibiotics might be present.
- Wear gloves during the test.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.
- Used microtiter stripes [PLATE] and materials that have been in contact with patient samples should be handled and disposed as potentially infectious.

13. GENERAL NOTES ON THE TEST

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- ID-Vit® is a trademark of Immundiagnostik AG.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be followed.
- Incubation time and temperature as well as pipetting volumes of the components are defined by the producer. Any variation of the test procedure which is not coordinated with the producer may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- The assay should always be performed according to the enclosed manual.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.
- Control samples should be analysed with each run.
- Serious incidents are to be reported to Immundiagnostik AG and the national regulatory authorities.

14. REFERENCES

1. Obeid, R. & Herrmann, W., 2006. Mechanisms of homocysteine neurotoxicity in neurodegenerative diseases with special reference to dementia. *FEBS letters*, **580**(13), pp.2994–3005.
2. Strohecker, R. & Henning, H., 1963. Vitamin-Bestimmungen. Erprobte Methoden. E. Merck AG, ed., Weinheim/Bergstraße: Verlag Chemie GmbH.
3. Verhaar, M.C., Stroes, E. & Rabelink, T.J., 2002. Folates and cardiovascular disease. *Arteriosclerosis, thrombosis, and vascular biology*, **22**(1), pp.6–13.

Used symbols:

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use
	Consult specification data sheet		

Immundiagnostik AG

Stubenwald-Allee 8a
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

info@immundiagnostik.com

www.immundiagnostik.com