# IDK® TurbiFIT®

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## 14. LITERATUR


1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2)CD001216.
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. Inflamm Bowel Dis, Vol 15 No 5; May 2009, p. 756–759.
3. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811–822.
4. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808–810.

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## Manual

# IDK® TurbiFIT®

*For the quantitative in vitro determination  
of faecal occult blood in stool*


Valid from 2022-03-24

TU1031.1  
TU1031.2  
TU1031.8

100/200/800

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 70190-363

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

# Table of Contents

<b>1. INTENDED USE</b>	<b>15</b>
<b>2. INTRODUCTION</b>	<b>15</b>
<b>3. MATERIAL SUPPLIED</b>	<b>15</b>
<b>4. MATERIAL REQUIRED BUT NOT PROVIDED</b>	<b>16</b>
<b>5. STORAGE AND PREPARATION OF REAGENTS</b>	<b>16</b>
<b>6. STORAGE AND PREPARATION OF SAMPLES</b>	<b>16</b>
<i>Sample storage and stability</i>	16
<i>Extraction of the stool samples</i>	17
<b>7. ASSAY PROCEDURE</b>	<b>18</b>
<i>Principle of the test</i>	18
<i>Test Procedure</i>	18
<b>8. QUALITY CONTROL</b>	<b>19</b>
<i>Reference Range</i>	20
<b>9. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Analytical sensitivity</i>	20
<i>Prozone</i>	20
<i>Clinical sensitivity and specificity</i>	20
<i>Interferences</i>	21
<i>Cross reactivity</i>	21
<b>10. LIMITATIONS</b>	<b>21</b>
<b>11. PRECAUTIONS</b>	<b>22</b>
<b>12. DISPOSAL</b>	<b>22</b>
<b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>22</b>
<b>14. REFERENCES</b>	<b>23</b>
<b>15. ANNEX - COMPATIBLE TURBIDIMETRIC ANALYSERS</b>	<b>23</b>

## 1. INTENDED USE

IDK® TurbiFIT® assay is a latex turbidimetric assay for the quantitative detection of haemoglobin (faecal occult blood) in human stool samples.

This assay is simple and widely applicable. This product is optimised for several automated analyser. For professional *in vitro* diagnostic use only.

## 2. INTRODUCTION

Colorectal cancer is the second leading cause of illness and death in the Western world. The screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps can be detected. Colorectal cancer is also associated with local acute inflammatory reaction being visualised, in some cases, by white cell neutrophil scanning.

Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases.

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. IDK® TurbiFIT® detects intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine.

## 3. MATERIAL SUPPLIED

Label	Kit components	Quantity for cat. no.		
		TU1031.1	TU1031.2	TU1031.8
Reagent 1	Reagent 1, ready-to-use	1 x 22 ml	2 x 22 ml	8 x 22 ml
Reagent 2	Reagent 2, ready-to-use	1 x 8 ml	1 x 13 ml	4 x 13 ml
CAL 0–5	Calibrators 0–5, ready-to-use	6 x 1 ml	6 x 1 ml	6 x 1 ml
CTRL 1	Control 1, ready-to-use	2 x 1 ml	2 x 1 ml	2 x 1 ml
CTRL 2	Control 2, ready-to-use	2 x 1 ml	2 x 1 ml	2 x 1 ml

The supplied reagents are in neutral tubes that are not compatible with the used analyzer. Therefore, it is necessary to fill the corresponding reagents into compatible reagent vessels.

#### 4. MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Centrifuge, 3 000 g
- Vortex
- Calibrated precision pipettors (adjustable volume) and 10–1 000 µl single-use tips
- Disposable gloves and laboratory equipment.
- Automated analyser (list of compatible analysers see annex information)
- Reagent vessels that are approved for the used analyser
- **IDK® TurbiTUBE®** for sample collection (cat. no.: TU0091.A)

#### 5. STORAGE AND PREPARATION OF REAGENTS

- All test reagents are ready-to-use. Test reagents are stable until the expiry date (see label) when stored at 2–8 °C. Do not freeze any component of the kit. Keep away from direct sunlight.
- Always consider that the reagent containers must be properly closed to avoid any contamination.

#### 6. STORAGE AND PREPARATION OF SAMPLES

##### *Sample storage and stability*

Collect a sufficient amount of human stool in a clean and dry container (no preservatives or transport media is necessary).

The raw stool samples can be stored at 2–8 °C for one day. If not immediately tested, stool samples can be frozen and stored at -20 °C for a maximum of 6 months. In this case, the sample must be totally thawed, and brought to room temperature (15–30 °C) before testing.

Stool samples must be homogenised as thoroughly as possible prior to preparation.

Extracted stool samples are stable at 2–8 °C for 7 days, as well as 7 days at room temperature.

## Extraction of the stool samples

Stool Sample Application System (IDK® TurbiTUBE®) (cat. no.: TU0091.A)

### **Stool sample tube – Instructions for use**

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

**IDK® TurbiTUBE®** with 1.5 ml sample extraction buffer:

Applied amount of stool: 15 mg

Buffer Volume: 1.5 ml

Dilution Factor: 1:100

Please follow the instructions for the preparation of stool samples using the **IDK® TurbiTUBE®** as follows:

- a) The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b) Unscrew the tube (yellow part of cap) to open. Insert the yellow dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.  
**Hint:** In case of liquid stool samples we recommend to pipette 15 µl of liquid sample into the tube.
- c) Vortex the tube well until no stool sample remains in the notches.  
**Important:** Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with sample extraction buffer for ~ 10 minutes improves the result.
- d) Samples should be centrifuged for ~10 minutes at 3000g to obtain a sediment-free sample. Floating material such as grain shells can be neglected.
- e) Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

## 7. ASSAY PROCEDURE

### *Principle of the test*

**IDK® TurbiFIT®** is a latex turbidimetric assay for the quantitative detection of haemoglobin (faecal occult blood) in human stool samples.

The **IDK® TurbiFIT®** assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies. Antigens contained in the sample bind with the antibodies' anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample, the mixture will maintain its clear appearance.

### *Test Procedure*

#### **Preparation of the calibration curve**

For calibration use only **IDK® TurbiFIT®** calibrator vials. They are ready-to-use liquid calibrators containing human haemoglobin at different concentrations. The exact concentrations are indicated on the label of each vial.

Calibration of the system at least once a month is highly recommended. The system must be recalibrated when the reagent lot is changed or when the controls are out of the assigned range (see the product specification sheet provided with the kit).

#### **Calibration and control vials**


Fig 1: Prepared calibrators (5 points and blank) and controls for performing the calibration curve.

Label	Concentration	Lid Colour
CAL 0	0 ng/ml	Transparent
CAL 1	50 ng/ml	pink
CAL 2	100 ng/ml	purple
CAL 3	250 ng/ml	blue
CAL 4	500 ng/ml	brown
CAL 5	1 000 ng/ml	white
CTRL 1	Control 1*	green
CTRL 2	Control 2*	red

\* The batch-dependent concentration is given by the product specification sheet.

## 8. QUALITY CONTROL

For quality control use only *IDK® TurbiFIT®* Control 1 (CTRL 1) and *IDK® TurbiFIT®* Control 2 (CTRL 2). These are two ready-to-use liquid controls containing different concentrations of human haemoglobin. The respective concentration for each vial is given by the product specification sheet.

The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running extracted patient samples to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents and the technique must be reviewed.

Since the sample dilution is already considered in the calibration curve, the dilution factor is 1.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## Reference Range

We recommend each laboratory to establish its own reference concentration range.

### Reference range in stool samples

For colorectal cancer screening the following cut-off is estimated:

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If the dilution factor is 100, the conversion factor is:

1 ng/ml = 0.1 µg/g

10 ng/ml = 1 µg/g

## 9. PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

Limit of Detection (LoD): See application note.

### Prozone

Studies have been carried out to a concentration of 1 mg of human haemoglobin/g of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

*Data obtained by the analyser Biolis 24i (Tokyo Boeki).*

### Clinical sensitivity and specificity

An evaluation with 100 stool samples was performed comparing a turbidimetric assay (IDK® TurbiFIT®) and another commercial immunoassay (Eiken FOB Latex®, EIKEN Chemicals, CO.LTD.). The results were as follows:

	Sensitivity	Specificity
IDK® TurbiFIT® vs Eiken FOB Latex®	95.7 %	90 %

In addition, a comparison of **IDK® TurbiFIT®** and **IDK® Hemoglobin ELISA** was performed, with 45 stool samples. The results were as follows:

	Sensitivity	Specificity
<b>IDK® TurbiFIT®</b> vs <b>IDK® Hemoglobin ELISA</b>	<b>100%</b>	<b>94%</b>

The results showed a high sensitivity and specificity to detect human haemoglobin (hHb) using **IDK® TurbiFIT®**.

### *Interferences*

An evaluation was performed to determine the interferences of **IDK® TurbiFIT®**. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

### *Cross reactivity*

An evaluation was performed to determine the cross reactivity against other faecal markers occasionally present in faeces, such as: human calprotectin, bovine and pig haemoglobin, bovine and human transferrin.

Cross reactivities were not found for bovine and pig haemoglobin at concentrations below 0.2 µg/ml, 0.9 µg/ml and 250 µg/ml, respectively. Thus, no cross reactivities were found at physiological concentrations.

## **10. LIMITATIONS**

**IDK® TurbiFIT®** should only be used in human stool samples. **IDK® TurbiFIT®** must be used with properly extracted faecal samples.

The results obtained determine the concentration of human haemoglobin in faecal samples. A positive result should be contrasted with additional diagnostic procedures to determine the cause and source of the occult blood in stool.

If symptoms persist further diagnostic procedures should be carried out. Negative results do not exclude bleeding, as some polyps and types of colorectal cancer may bleed intermittently or not during certain stages of the disease.

This assay may be less sensitive for the detection of upper gastrointestinal bleeding as blood is degraded during the gastrointestinal passage.

Patients samples should not be collected during the menstrual cycle, if they have bleeding haemorrhoids, blood in urine or constipation.

## 11. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- A trained person in the turbidimetric technique and the use of the auto-analyser is required. Manual application of the test kit is not possible.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Reagents should not be used if the primary packaging is damaged or opened.
- Do not interchange different lot numbers of any kit component within the same assay.
- Only **new** empty reagent vessels are strictly recommended, because washed or recycled/reused reagent vessels may influence the test result.

## 12. DISPOSAL

- Liquid test components, pipets tips, tubes etc. are to be treated as ordinary laboratory waste, unless otherwise stated. The solutions should be discarded in a proper container after testing following local regulations.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. Other potentially infectious materials (e.g. sample collection container) must be disposed in accordance with official regulations.

## 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, eye protection and mask. Do not eat, drink or smoke in the working area.

- In the event of warranty claims, notifications of defects or serious incidents, notify Immundiagnostik AG immediately, but at the latest within 14 days after the problem occurring. Please send the product together with a written complaint (including e-mail, address and telephone number) to Immundiagnostik AG.
- If stored properly (2–8°C), Immundiagnostik AG guarantees, after the first opening, no loss of activity of the kit components for the following three months and at least 20 opening-closing cycles.
- Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- **IDK® TurbiFIT®** is a trademark of Immundiagnostik AG.

## 14. REFERENCES

1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. *Cochrane Database Syst Rev.* 2000;(2)CD001216.
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756–759.
3. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. *Ann Intern Med.* 1997; 126: 811–822.
4. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. *Ann Intern Med.* 1997; 126: 808–810.


## 15. ANNEX - COMPATIBLE TURBIDIMETRIC ANALYSERS

<b>Abbott</b>	Alinity c-series
<b>Abbott</b>	Architect c-series
<b>Awareness</b>	Chemwell-T
<b>Beckman Coulter</b>	AU 480 / 680
<b>Biobase</b>	BK Mini

<b>Biosystems</b>	A15
<b>Binding Site</b>	Optilite
<b>DiaSys</b>	Respons 910
<b>Erba</b>	XL-180 / XL-200
<b>Furano</b>	CA800
<b>Mindray</b>	BS200E
<b>Roche</b>	Cobas c111 / c311 / c501 / c502
<b>Siemens</b>	Atellica
<b>Tecom</b>	TC220
<b>Tokyo Boeki</b>	Biolis 24i / 30i / 50i
<b>Ortho Clinical Technologies</b>	Vitros 5600

The application notes for use our turbidimetric products with the various devices can be obtained on request from our technical support.

**Used symbols:**

	Temperature limitation		Catalogue number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		Consult instructions for use

## **Immundiagnostik AG**

Stubenwald-Allee 8a  
64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: +49 6251 70190-363

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)