



Data Sheet

SIALYLATED MUCIN-TYPE O-GLYCANS

ANTIBODY, MONOCLONAL

Catalog no.:	A 1053.1 / A 1053.2
Immunogen:	Purified mucus glycoproteins (mucins) from amniotic fluid
Host:	Mouse Balb/c
Clone no.:	FW17
Isotype:	IgM
Matrix:	Cell culture supernatant, purified, PBS pH 7.4
Specificity:	Sialylated mucin-type O-glycans, oncofetal antigen
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA, immunocytochemistry, immunohistochemistry (paraffin sections) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	05 July 2016

For research use only

Publishing research using A 1053? Please let us know so that we can cite your publication as a reference.

