

Data Sheet

HUMAN WISP3 / CCN6

ANTIBODY, POLYCLONAL

Catalogno.:	A 2101.1 / A 2101.2
Immunogen:	Recombinant human WISP3-flag
Synonyms:	WNT1-inducible-signaling pathway protein 3, CCN family member 6
Swiss-Prot No:	O95389
Gene Information:	Gene Name: WISP3, CCN6 GeneID: 8838
Host:	Rabbit
Matrix:	Serum
Specificity:	Human WISP3
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	11 August 2016

For research use only

Publishing research using A 2101? Please let us know so that we can cite your publication as a reference.

