

Data Sheet

HUMAN CALCITONIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 4111.1 / A 4111.2
Immunogen:	Synthetic human Calcitonin, BSA-conjugated (CGNLSTCMLGTYTQDFNKFHTFPQTAIGVGAP)
Swiss-Prot No:	P01258
Gene Information:	Gene Name: CALCA, CALC1 GeneID: 796
Host:	Sheep
Matrix:	Serum
Specificity:	Human Calcitonin There was no cross reactivity obtained with human Katalcalcin or human Calcitonin Gene related Peptide.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:80 000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	19 August 2016

For research use only

Publishing research using A 4111? Please let us know so that we can cite your publication as a reference.

