

## Data Sheet

# HUMAN CALCITONIN

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 4112.1 / A 4112.2
<b>Immunogen:</b>	Synthetic human Calcitonin, BSA-conjugated (CGNLSTCMLGTYTQDFNKFHTFPQTAIGVGAP)
<b>Swiss-Prot No:</b>	P01258
<b>Gene Information:</b>	Gene Name: CALCA, CALC1 GeneID: 796
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Calcitonin  There was no cross reactivity obtained with human Katalcalcin or human Calcitonin Gene related Peptide.
<b>Contents:</b>	20 µl / 100 µl (lyophilized)  Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA (1:15 000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved)  Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	19 August 2016

### For research use only

**Publishing research using A 4112? Please let us know so that we can cite your publication as a reference.**

