

## Data Sheet

# HUMAN CALCITONIN GENE-RELATED PEPTIDE 2

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 4612.1 / A 4612.2
<b>Immunogen:</b>	Synthetic human Calcitonin gene-related peptide 2, bTG-conjugated (ACNTATCVTHRLAGLLSRSGGMVKSNTFVPTNVGSKAF)
<b>Synonyms:</b>	Beta-type CGRP, Beta-CGRP, Calcitonin gene-related peptide II, CGRP-II
<b>Swiss-Prot No:</b>	P10092
<b>Gene Information:</b>	Gene Name: CALCB, CALC2 GeneID: 797
<b>Host:</b>	Goat
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Calcitonin Gene-related Peptide There was no cross reactivity obtained with human and salmon Katalcalcin or Calcitonin.
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA (1:3000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 July 2016

### For research use only

**Publishing research using A 4612? Please let us know so that we can cite your publication as a reference.**

