



Data Sheet

HUMAN CALCITONIN GENE-RELATED PEPTIDE 2 SYNTHETIC

Catalog No.:	A 4612AG.1 / A 4612AG.2
Name:	Synthetic human Calcitonin gene-related peptide 2
Synonyms:	Beta-type CGRP, Beta-CGRP, Calcitonin gene-related peptide II, CGRP-II
Sequence:	ACNTATCVTHRLAGLLSRSGGMVKSNFVPTNVGSKAF
MW:	3.8 kDa
Swiss-Prot No:	P10092
Gene Information:	Gene Name: CALCB, CALC2 Gene ID: 797
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	-20°C (lyophilized); - 20°C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 June 2019

For research use only

Publishing research using A 4612AG? Please let us know so that we can cite your publication as a reference.



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 (0) 62 51/70 19 00 · Fax: +49 (0) 62 51/84 94 30 · info@immundiagnostik.com · www.immundiagnostik.com