



Data Sheet

HUMAN CALCITONIN GENE-RELATED PEPTIDE 1 SYNTHETIC

Catalog No.:	A 4613AG.1 / A 4613AG.2
Name:	Synthetic human Calcitonin gene-related peptide 1
Synonyms:	Alpha-type CGRP, Calcitonin gene-related peptide I, CGRP-I
Sequence:	ACDTATCVTHRLAGLLSRSGGVKNNFVPTNVGSKAF
MW:	3.79 kDa
Swiss-Prot No:	P06881
Gene Information:	Gene Name: CALCA, CALC1 Gene ID: 796
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	-20°C (lyophilized); - 20°C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 June 2019

For research use only

Publishing research using A 4613AG? Please let us know so that we can cite your publication as a reference.

