



Data Sheet

HUMAN C-PEPTIDE

SYNTHETIC

Catalog No.:	A 5510AG.1 / A 5510AG.2
Name:	Human C-Peptide trifluoroacetate salt (EAEDLQVGQVELGGPGAGSLQPLALEGSLQ)
Synonyms:	Insulin (aa 57-87)
MW:	3020.3 Da
Swiss-Prot No:	P01308
Gene Information:	Gene Name: INS GeneID: 3630
Purity:	92% (HPLC)
Source:	Synthetic
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C (lyophilized); - 20°C (dissolved)
Last updated on:	13 June 2019

For research use only

Publishing research using A 5510AG? Please let us know so that we can cite your publication as a reference.

