

Data Sheet

HUMAN CARCINOEMBRYONIC ANTIGEN [CEA]

ANTIBODY, POLYCLONAL

Catalog no.:	A 7222.1 / A 7222.2
Immunogen:	Purified human carcinoembryonic antigen
Synonyms:	Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5), Meconium antigen 100, CD_antigen: CD66e
Swiss-Prot No:	P06731
Gene Information:	Gene Name: CEACAM5 GenelD: 1048
Host:	Rabbit
Matrix:	Serum
Specificity:	Human carcinoembryonic antigen
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:1000-2000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	15 October 2015

For research use only

Publishing research using A 7222? Please let us know so that we can cite your publication as a reference.

