



Data Sheet

HUMAN CARCINOEMBRYONIC ANTIGEN [CEA] NATIVE

Catalog no.:	A 7222AG.1 / A 7222AG.2
Name:	Human carcinoembryonic antigen
Synonyms:	Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5), Meconium antigen 100, CD_antigen: CD66e
MW:	180 kDa
Swiss-Prot No:	P06731
Gene Information:	Gene Name: CEACAM5 Gene ID: 1048
Source:	Human tissue
Purity:	> 95% (HPLC)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Reference:	
Last updated on:	7 March 2018

For research use only

Publishing research using A 7222AG? Please let us know so that we can cite your publication as a reference.

