



Data Sheet

3-EPI-25-HYDROXYVITAMIN-D₃

SYNTHETIC

Catalog no.: AD1007AG.2
Name: 3-Epi-25-Hydroxyvitamin-D₃
Synonyms: CAS 73809-05-9



Structure: 

Molecular formula: C₂₇H₄₄O₂

MW: 400.65

Source: Synthetic

Purity: >98%

Appearance: Colourless, amorphous solid

Solubility: Soluble in EtOH/MeOH, DMF and DMSO, well soluble in halogenated hydrocarbons, not soluble in water

Contents: 1 mg (lyophilized)
Resuspend in 99% EtOH p.a. for required concentration

Store at: - 20 °C, protect from heat and light

Last updated on: 14 June 2017

For research use only

Publishing research using AD1007AG? Please let us know so that we can cite your publication as a reference.

