



Data Sheet

HUMAN PRO-C-TYPE NATRIURETIC PEPTIDE [proCNP] (aa 51-79)

SYNTHETIC

Catalog No.:	AE1002AG.1 / AE1002AG.2
Name:	Human proCNP (aa 51-79)
Sequence:	DLRVDTKSRAAWARLLGEHPNARKYKGAN
Swiss-Prot No:	P23582
Gene Information:	Gene Name: NPPC, CNP2 GeneID: 4880
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	-20°C (lyophilized); - 20°C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 June 2019

For research use only

Publishing research using AE1002AG? Please let us know so that we can cite your publication as a reference.

