



## Data Sheet

# HUMAN CHEMOKINE CC-1/3 (aa 62-74)

## SYNTHETIC

<b>Catalog no.:</b>	AE1010AG.1 / AE1010AG.2
<b>Name:</b>	Human Chemokine CC-1/3 (aa 62-74)
<b>Synonyms:</b>	C-C motif chemokine 14, HCC-1/HCC-3, HCC-1(1-74), NCC-2, Small-inducible cytokine A14
<b>Sequence:</b>	KWVQDYIKDMKEN
<b>MW:</b>	1.7 kDa
<b>Swiss-Prot No:</b>	Q16627
<b>Gene Information:</b>	Gene Name: CCL14, NCC2, SCYA14 Gene ID: 6358
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1 % acetic acid for required concentration.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	23 April 2018

**For research use only**

**Publishing research using AE1010AG? Please let us know so that we can cite your publication as a reference.**

