

Data Sheet

HUMAN CHEMOKINE CC-1/3 (aa 1-16)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1011.1 / AE1011.2
Immunogen:	Synthetic human Chemokine CC-1/3 (aa 1-16) (TKTESSSRGPYHPSEC)
Synonyms:	C-C motif chemokine 14, Small-inducible cytokine A14
Swiss-Prot No:	Q16627
Gene Information:	Gene Name: CCL14 GeneID: 6358
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Chemokine CC-1/3 (aa 1-16)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:10 000); RIA (1:10 000); Western Blot (1:10 000); Immunohistochemistry (1:10 000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	08 March 2018

For research use only

Publishing research using AE1011? Please let us know so that we can cite your publication as a reference.

