



Data Sheet

HUMAN CHEMOKINE CC-2 (aa 18-31)

SYNTHETIC

Catalog no.:	AE1013AG.1 / AE1013AG.2
Name:	Human Chemokine CC-2 (aa 18-31)
Synonyms:	HCC-2, C-C motif chemokine 15, Leukotactin-1 (LKN-1), MIP-1 delta, Macrophage inflammatory protein 5 (MIP-5), Mrp-2b, NCC-3, Small-inducible cytokine A15
Sequence:	ENPVVLSNFHFAAD
MW:	1.56 kDa
Swiss-Prot No:	Q16663
Gene Information:	Gene Name: CCL15 Gene ID: 6359
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1 % acetic acid for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2018

For research use only

Publishing research using AE1013AG? Please let us know so that we can cite your publication as a reference.

