



## Data Sheet

# RAT GUANYLIN (aa 78-92)

## SYNTHETIC

<b>Catalog no.:</b>	AE1019AG.1 / AE1019AG.2
<b>Name:</b>	Rat Guanylin (aa 78-92) active peptide
<b>Synonyms:</b>	Guanylate cyclase activator 2A
<b>Sequence:</b>	PNTCEICAYAACTGC
<b>Swiss-Prot No:</b>	P28902
<b>Gene Information:</b>	Gene Name: Guca2a, Guca2 GeneID: 25656
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	14 March 2018

### For research use only

**Publishing research using AE1019AG? Please let us know so that we can cite your publication as a reference.**

