



Data Sheet

HUMAN pro-C-TYPE NATRIURETIC PEPTIDE (proCNP) [aa 1-19]

SYNTHETIC

Catalog No.:	AE1028AG.1 / AE1028AG.2
Name:	Human pro-CNP (aa 1-19)
Sequence:	KPGAPPKVPRTPPAEELAE
MW:	1.98 kDa
Swiss-Prot No:	P23582
Gene Information:	Gene Name: NPPC, CNP2 GeneID: 4880
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	-20°C (lyophilized); - 20°C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 June 2019

For research use only

Publishing research using AE1028AG? Please let us know so that we can cite your publication as a reference.

