



Data Sheet

HUMAN BRAIN NATRIURETIC PEPTIDE [BNP] (aa 1-10) SYNTHETIC

Catalog no.:	AE 1034AG.1 /AE 1034AG.2
Name:	Human BNP (aa 1-10)
Synonyms:	Natriuretic peptides B, Gamma-brain natriuretic peptide
Sequence:	SPKMQGSGC(G)
Swiss-Prot No:	P16860
Gene Information:	Gene Name: GeneID: 4879
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg /1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	2 March 2018

For research use only

Publishing research using AE1034AG? Please let us know so that we can cite your publication as a reference.

