



Data Sheet

HUMAN SPERM-ASSOCIATED ANTIGEN 11B [SPAG11B] ISOFORM C (aa 52-88)

SYNTHETIC

Catalog no.:	AE1041AG.1 / AE1041AG.2
Name:	Human SPAG11B isoform C (aa 52-88)
Synonyms:	Human epididymis-specific protein 2 (He2), Protein EP2, Sperm antigen HE2
Sequence:	LKVVDCCRSEGFQCYCNVYMETQVGYCSKSKKDACCLH
Swiss-Prot No:	Q08648-3
Gene Information:	Gene Name: SPAG11B, EP2, HE2 GeneID: 10407
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	15 March 2018

For research use only

Publishing research using AE1041AG? Please let us know so that we can cite your publication as a reference.

