



Data Sheet

HUMAN β -DEFENSIN 4 (aa 3-39)

SYNTHETIC

Catalog no.:	AE1046AG.1 / AE1046AG.2
Name:	Human β -Defensin 4 (aa 3-39)
Sequence:	ELDRICGYGTARCRKKCRSQEYRIGRCPNTYACCLRK
Synonyms:	hBD-4, β -defensin 104, DEFB-4
Swiss-Prot No:	Q8WTQ1
Gene Information:	Gene Name: DEFB104A, DEFB104, DEFB4 and DEFB104B GeneID: 140596 and 503618
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	9 March 2018

For research use only

Publishing research using AE1046AG? Please let us know so that we can cite your publication as a reference.

