



## Data Sheet

# HUMAN $\beta$ -DEFENSIN 1 (aa 18-26)

## SYNTHETIC

<b>Catalog no.:</b>	AE1049AG.1 / AE1049AG.2
<b>Name:</b>	Human $\beta$ -Defensin 1 (aa 18-26)
<b>Sequence:</b>	PIFTKIQGT
<b>Synonyms:</b>	$\beta$ -Defensin 1 (BD-1), Defensin $\beta$ 1
<b>Swiss-Prot No:</b>	P60022
<b>Gene Information:</b>	Gene Name: DEFB1 GeneID: 1672
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	9 March 2018

**For research use only**

**Publishing research using AE1049AG? Please let us know so that we can cite your publication as a reference.**

