



Data Sheet

HUMAN β -DEFENSIN 1 (aa 1-36)

SYNTHETIC

Catalog no.:	AE1109AG.1 / AE1109AG.2
Name:	Human β -Defensin 1 (aa 1-36)
Sequence:	DHYNCVSSGGQCLYSACPIFTKIQGTCYRGKAKCCK
Synonyms:	β -Defensin 1 (BD-1), β 1
Swiss-Prot No:	P60022
Gene Information:	Gene Name: DEFB1 GeneID: 1672
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	9 March 2018

For research use only

Publishing research using AE1109AG? Please let us know so that we can cite your publication as a reference.

